

Oct 24, 2020

Fluorescent labeling *Bacillus mycoides*

DOI

dx.doi.org/10.17504/protocols.io.bkuukwww


Andreea S¹

¹University of Groningen

iGEM Groningen 2020



r.rajadhyaksha



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bkuukwww

Protocol Citation: Andreea S 2020. Fluorescent labeling *Bacillus mycoides*. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bkuukwww>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Other

The protocol is developed based on literature and has not been tested yet.

Created: September 04, 2020

Last Modified: October 24, 2020

Protocol Integer ID: 41588



Abstract

GFP method is performed in bacteria to localize a desired peptide in the bacteria. In this experiment, we are focusing on the Nlp-14a. In the paper written by (Yi,2017) GFP was mutated on a B. mycoides E18 strain to study the plant-interaction studies. In their paper, they constructed a superfolder green fluorescent protein (sfGFP) and red fluorescent protein (mKate2) in the B. mycoides E18 strain. They tracked the GFP- tagged bacterial strain after inoculating in Chinese cabbage plants in a hydroponic system. They noticed that the bacterial strain rapidly attached to the plant during colonization and formed a matrix.

CITATION


Yi, Y., Frenzel, E., Spoelder, J., Elzenga, J. T. M., van Elsas, J. D., & Kuipers, O. P. (2017). Optimized fluorescent proteins for the rhizosphere-associated bacterium *Bacillus mycoides* with endophytic and biocontrol agent potential. *Environmental microbiology reports*, 10(1), 57-74.

LINK

[10.1111/1758-2229.12607](https://doi.org/10.1111/1758-2229.12607)



Strains and Growth conditions of bacteria




- 1 Routinely culture *Bacillus* strains in Luria-Bertani (LB) medium at  30 °C with aeration at 200 rpm.

- 2


Preparation of *B. mycoides* cells for GFP




- 3 Prepare the *B. mycoides* strain aliquotes for electroporation.



- 3.1 Cultivate the bacterial strain overnight in LB broth at  30 °C and at 180 rpm.

- 3.2 Transfer  1 mL of the overnight culture into  100 mL of LB medium (with 2% [wt/wt] glycine) and incubate it at  30 °C and 180rpm until optical density at 600 nm is 0.4 to 0.7.

- 3.3 Centrifuge the cells and wash the pellets with increasing concentrations of ice-cold glycerol (2.5%, 5%, and 10%). Resuspend this pellet in precooled electroporation buffer (10% glycerol) and shock freeze in liquid nitrogen.





- 4 Add the library vector DNA in an amount of  2 µg to the cells, and perform electroporation. The settings for electroporation are 2.0 kV, 25 IF and 200 X in a 2-mm cuvette using a Bio Rad Gen Pulser II electroporation system (Bio-Rad).

- 5 Add  1 mL of LB medium and grow the cells for  02:00:00 at  30 °C and 150 rpm for recovery and then plate on LB-Cm4 agar.


- 6 After  24:00:00 of growth at  30 °C, harvest the colonies from the plates and pool in LB medium

- 7 Store the libraries at  -80 °C as 15% glycerol stocks.

FACS

- 8 Inoculate the *B. mycoides* strain *mKate2mut* library in  50 mL of LB-Cm4 and grow at  7 or  6 to an OD_{600nm} of 0.3-0.6.
- 9 *B. mycoides* has been seen to show extensive cell-chaining and hence a mild sonication step of 4 rounds of 3 X 10 pulses of 1s with an amplitude of 30% can be applied to disassemble the aggregated cells.
- 10 Sort the cells on a flow cytometer at 20 psi using a  70 micromolar (μM) nozzle at a flow rate of 1.0 with the highest sort precision mode (0– 32-0 sort purity mask).
- 11 Using a sequential gating strategy with FCS height versus widths, followed by SCC height versus width, cellular debris, and chained cells can be excluded.
- 12 To separate the brightest variants choose a cutoff of 3% of the brightest event in the first round of cell sorting and 0.3% of the brightest events in the second round of sorting with the light scatter parameters.

Screening of FP variants and flow cytometry measurements

- 13 After FACS sorting, plate the final fluids containing bright cells on LB-Cm4 plates and grown them overnight at  30 °C .
- 14 Observe the colonies using a fluorescence microscope. Keep the filter setting for GFP as excitation at 460/480nm and emission at 495/540 nm with a 485 nm dichromatic mirror; for RFP, the filter setting can be kept as excitation at 545/580 nm and emission at 610 nm with a 600 nm dichromatic mirror.
- 15 Capture the images on a camera and calculate the intensity of single-cell with Image J software.
- 16 Calculate the total cell fluorescence the formula is: Corrected total cell fluorescence (CTCF) = Integrated Density – (Area of selected cell x Mean fluorescence of background readings)

Citations

Yi, Y., Frenzel, E., Spoelder, J., Elzenga, J. T. M., van Elsas, J. D., & Kuipers, O. P.. Optimized fluorescent proteins for the rhizosphere-associated bacterium *Bacillus mycoides* with endophytic and biocontrol agent potential
[10.1111/1758-2229.12607](https://doi.org/10.1111/1758-2229.12607)