

Nov 21, 2024

FoldSeek screen

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqd563vk5/v1>

James Marshall¹, Eleanor Warren¹, John Bergqvist¹

¹MRC Laboratory of Medical Sciences, Imperial College London

Behavioural Genomics



John Bergqvist

MRC Laboratory of Medical Sciences, Imperial College London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqd563vk5/v1>

Protocol Citation: James Marshall, Eleanor Warren, John Bergqvist 2024. FoldSeek screen. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgqd563vk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 12, 2024

Last Modified: November 21, 2024

Protocol Integer ID: 111966

Keywords: proteins with insecticidal potential, worms under bluelight stimulation, od600, recording worm, insecticidal potential, protein, expressing protein

Abstract

This protocol outlines the screening procedure of a library of transformed *E. coli* expressing proteins with insecticidal potential. The two outputs from this procedure is 1) the OD₆₀₀ as a measurement of how much *C. elegans* have consumed the transformed *E. coli*, and 2) recording worms under bluelight stimulation to capture variations in motility.

Materials

100 mm Petri dishes

Round-bottom Nunclon 96-well plates

2 mL deep 96-well plates

Breathable film cover for 96-well plates

Lysogeny Broth medium

Nematode Growth Medium

M9 medium

S-medium

Kanamycin

Nystatin

IPTG

Bleach solution

MilliQ water

C. elegans (N2)

E. coli (OP50)



Day 1 - Preparing worms

30m

1 Egg laying

Pick 20 adult worms onto six 100 mm plates seeded with  1.5 mL of OP50. Leave on plates for  24:00:00 h .

30m



Day 2 - Removing adult worms & making S-medium

2h 20m

2 Remove adult worms from the 100 mm plates if sufficient eggs have been laid

20m

3 Check if S-medium is available.

If not, make S-medium.

2h

Note

1 L Pre-autoclave S-medium (S-basal)

NaCl	 5.85 g	4th floor
K ₂ HPO ₄	 1 g	4th floor
KH ₂ PO ₄	 6 g	4th floor
MilliQ water	 1 L	4th floor

**Note****1 L Post-autoclave S-medium**

**Kanamycin, Nystatin and IPTG do not last long and so should only be added to smaller aliquotes upon use.*

Cholesterol ([M] 5 mg/mL in ethanol)	1 mL	4th floor
1 M potassium citrate 6	10 mL	4th floor
trace metals solution	10 mL	4th floor
1 M CaCl ₂	3 mL	4th floor
1 M MgSO ₄	3 mL	4th floor
1000x Kanamycin*	1 mL	4th & 6th floor
5000x Nystatin*	200 µL	4th floor
1 M IPTG*	100 µL	4th & 6th floor

If needed, make more 1 M IPTG by mixing 2.38g IPTG in 10 mL sterile water, followed by filtering using a 0.22 µm filter.

Day 4 - Beginning of Bacterial Growth

30m

4 Grow FoldSeek Library of plates

Fill the four pre-made FoldSeek transformed E. coli plates with 100 µL of LB broth+Kanamycin. Cover plates with breathable film cover.

30m



Leave overnight at 37 °C and shaking at 200 rpm.

Note

Make LB+Kanamycin by pipetting 500 µL of 1000x Kanamycin into 500 mL LB.

Day 5 - Inducing bacteria & bleaching worms

5h 30m

5 Dilute overnight cultures 1 in 100

Dilute 5 µL of overnight cultures in 450 µL LB broth+Kanamycin in 2 mL deep-well plates (LB+Kan can be filled by ViaFill using programme '16 DEEP WELL').

30m

**6 Grow diluted cultures to OD₆₀₀ 0.4 - 0.8**

3h

Incubate the plates at 37 °C 230 rpm for about 03:00:00 h until they have reached an OD₆₀₀ between 0.4 - 0.8.

7 Induce bacterial cultures with IPTG at OD₆₀₀ 0.4 - 0.8

30m

Make a 10x solution of 1 M IPTG (1000x) in LB+Kanamycin.

When cultures have grown to an OD₆₀₀ of 0.4 - 0.8, use ViaFill (programme '16 DEEP WELL') to dispense 50 µL of LB+Kan+IPTG in each well. Rinse ViaFill tubes with heated MilliQ water.

Leave overnight at 37 °C and 250 rpm.

Note

4x 96-well plate:

50 µL x 384 wells = 19.2 mL

The ViaFill requires about 5 mL for priming. As such, to ensure there is enough dispensed without risking bubbles, *make up a solution of 40 mL*.

10x IPTG in LB+Kanamycin (1/100 dilution):

400 µL of 1 M IPTG (1000x IPTG) in 40 mL of LB+Kanamycin.

8 Bleach worms and place in diapause

1h 30m

Need around 50-100,000 L1s. To be used on Day 6, so leave spinning overnight.

Day 6. Add worms to bacteria & measure

3h 15m

9 Complete S-medium

10m

Use a 500 mL measuring cylinder aliquote 400 mL of S-medium into a new autoclaved glass container. Add the following components to complete the S-medium:

**Note****400 mL complete S-medium:**1000x Kanamycin  400 μL 4th & 6th floor5000x Nystatin  80 μL 4th floor1 M IPTG  40 μL 4th & 6th floor

Kanamycin, Nystatin and IPTG available on 4th floor 'Toxic' -20°C freezer in box labelled 'Spare antibiotics/antifungals' or in yellow boxes in a -20°C freezer on 6th floor.

Make more 1 M IPTG by mixing 2.38g IPTG in 10 mL sterile water, followed by filtering using a 0.22 μm filter.

10 Pellet & re-suspend bacteria in S-medium

30m

Pellet plates for  00:10:00 min at 500x g. Flick out the LB and replace the medium with  500 μL of S-medium using ViaFill.

Note

1630x rpm for centrifuge 5810R to get 500x g.

11 Centrifuge again and refill with S-medium

5m

Fill the wells with  300 μL of S-medium using ViaFill. Rinse the ViaFill tubes with MilliQ water to ensure no clogging occurs.

12 Resuspend & make triplicates plates

1h

Using a multichannel pipette, resuspend the cultures thoroughly.

Pipette  80 μL into three separate Nunclon round-bottom 96-well plates.

13 Add worms to the plates

30m

Add  20 μL of L1 worms to each well (between 50 - 100 worms).

14 Measure OD

30m

Measure OD₆₀₀ of the plates using a plate reader (with 'James_OD_600' protocol).



File naming format: **FS21_1_JB1_Day1**

- **FS**: FoldSeek
- **21**: plate 2.1
- **1**: plate replicate 1 out of 3 of FS plate 2.1
- **JB1**: experiment repeat for *intials*
- **Day1**: day of the measurement out of 14.

15 Record plates using Hydras

30m

Plates recorded without lids using only one hydra for the entire screen.

Script: **python3 run_Short_bluelight.py -r 02 -f ToxAss_2024_MM_DD_JB1_d1_FS21_1**

- **MM** : month for when the first recording was made
- **DD** : day for when the first recording was made
- **JB1**: experiment repeat for *intials*
- **d1**: day of the measurement out of 14
- **FS21**: FoldSeek plate number 2.1
- **1**: plate replicate 1 out of 3 of FS plate 2.1

The script runs for 30s. When the blue light has stopped, the recording is done. Take out the plate and put in the new one. ***Make sure to click 'clear' under 'configurations' before starting to record the new plate.***

16 Store plates

Place plates in a box with a damp paper towel and store at  20 °C

Day 9 to Day 19 - Measuring OD₆₀₀ and motility

1h

17 Measure OD₆₀₀ using the plate reader

30m

Use the protocol 'James_OD_600'.

18 Record plates using Hydras

30m

Run the script 'run_Short_bluelight.py' to record the plates under bluelight.

Note

14 days in total for measurements (Day 6 - Day 19). Measurements done on Day 6 and from Day 9 to Day 19 of Experiment Days. If the first measurement starts on Monday, then the second day of measurement is on Thursday.

	Exp. Day	(Measurement day)	Day of week
1.	Day 6	(1)	Mon
2.	Day 9	(4)	Thurs
3.	Day 10	(5)	Fri
4.	Day 11	(6)	Sat
5.	Day 12	(7)	Sun
6.	Day 13	(8)	Mon
7.	Day 14	(9)	Tues
8.	Day 15	(10)	Wed
9.	Day 16	(11)	Thurs
10.	Day 17	(12)	Fri
11.	Day 18	(13)	Sat
12.	Day 19	(14)	Sun