



Jun 29, 2023

Generation of ATG3 KO Hela cells stably expressing HaloTag-LC3B

DOI

dx.doi.org/10.17504/protocols.io.5qpvo3xq9v4o/v1

Xuefeng Ren¹

¹Laboratory of James H. Hurley, University of California, Berkeley, CA



Xuefeng Ren

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo3xq9v4o/v1>

Protocol Citation: Xuefeng Ren 2023. Generation of ATG3 KO Hela cells stably expressing HaloTag-LC3B. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvo3xq9v4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 12, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 83216

Keywords: ATG3 KO Hela cells, HaloTag-LC3B, ASAPCRN, generation of atg3 ko hela cell, atg3 ko hela cell, hela cell, lc3b this protocol details generation, lc3b, cell, expressing halotag, halotag

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000350

Abstract

This protocol details generation of ATG3 KO Hela cells stably expressing HaloTag-LC3B.

Attachments



[736-1852.pdf](#)

129KB



Materials

Buffers and reagents:

- Growth broth: LB broth
- StbI3 competent *E. coli* (MacroLab, UC Berkeley)
- DMEM medium with GlutaMAX containing 10% FBS and 10% Pen-Strep.

⊗ pVSV-G addgene Catalog ##138479

⊗ pBS-CMV-gagpol addgene Catalog #35614

⊗ pMRX-IP-HaloTag7-LC3 addgene Catalog #184899

⊗ Qiagen Hi-Speed MidiPrep kit Qiagen Catalog #12643

⊗ DMEM, high glucose, GlutaMAX[®] Supplement Thermo Fisher Catalog #10566016

⊗ Gibco[™] Fetal Bovine Serum value heat inactivated (formerly USDA-approved in North America or quali Fisher Scientific Catalog #A5256801

⊗ Penicillin-Streptomycin (10,000 U/mL) Thermo Fisher Scientific Catalog #15140122

⊗ Opti-MEM[®]; I Reduced Serum Medium Thermo Fisher Catalog #31985070

⊗ TransIT[®]-LT1 Transfection Reagent Mirus Bio Catalog #MIR 2300

⊗ Lenti-X[™] Concentrator Takara Bio Inc. Catalog #631231

⊗ 15-Dimethyl-15-diazaundecamethylene polymethobromide Polybrene Merck MilliporeSigma (Sigma-Aldrich) Catalog #H9268

⊗ Puromycin Dihydrochloride Gold Biotechnology Catalog #P-600

Troubleshooting



Procedures

15m

- 1 Transform each plasmid into Stbl3 competent cells for the propagation. Next day, pick up one colony to inoculate  Overnight culture in LB medium with ampicillin. 
- 2 Midi prep the cultures to purify plasmids (Qiagen).
- 3 Plate 5×10^6 HEK 293T cells on a 10 cm plate in DMEM medium.
- 4 HEK293T transfection:
 - 4.1 Add retroviral packaging plasmids (pVSV-G, pBS-CMV-gagpol) and pMRX-IP-HaloTag7-LC3,  5 μ g each in  1.5 mL warm Opti-MEM medium. 
 - 4.2 Add  45 μ L of TransIT-LT1 transfection reagent (Mirus) and swirl. 
 - 4.3 Incubate at  Room temperature for  00:15:00 . 

 - 4.4 Add  1.5 mL dropwise into 10 cm HEK293T plate. 
 - 4.5 At 72 hours post-transfection, collect retroviral supernatant into a falcon tube.
- 5 Concentrate retroviral supernatant to  1 mL using Lenti-X concentrator (Takara) with the manufacturer instruction.
- 6 Retroviral transduction:
 - 6.1 Plate 1×10^5 ATG3 KO HeLa cells into 12-well plate one day before.



- 6.2 Next day, titrate  100 μL ,  200 μL ,  400 μL of concentrated retroviral solution with  8 $\mu\text{g}/\text{mL}$ Polybrene (Sigma) into target HeLa cells.
- 6.3 At 24 hours post-transduction, remove retroviral supernatant and replace with fresh DMEM complete medium with  2 $\mu\text{g}/\text{mL}$ puromycin (GoldBio).
- 6.4 Perform Puromycin selection for two weeks. Confirm Halo-LC3B positive cells by FACS sorting.