



Jun 28, 2023

🌐 Generation of stable ATG3-mCherry wild-type or mutants HeLa cells with HaloTag-LC3B using lentivirus

DOI

dx.doi.org/10.17504/protocols.io.dm6gp3exjvzp/v1

Xuefeng Ren¹

¹Laboratory of James H. Hurley, University of California, Berkeley, CA



Xuefeng Ren

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp3exjvzp/v1>

Protocol Citation: Xuefeng Ren 2023. Generation of stable ATG3-mCherry wild-type or mutants HeLa cells with HaloTag-LC3B using lentivirus. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp3exjvzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: June 12, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 83237

Keywords: generation of stable ATG3-mCherry wild-type, HaloTag-LC3B, lentivirus, ASAPCRN, mutant hela cells with halotag, hela cells with halotag, using lentivirus, mutant hela cell, generation of stable atg3, lentivirus this protocol detail, hela cell, stable atg3, lc3b, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000350

Abstract

This protocol details the generation of stable ATG3-mCherry wild-type or mutant HeLa cells with HaloTag-LC3B using lentivirus.

Attachments



[737-1853.pdf](#)


117KB




Materials

Buffers and reagents:


- XL10-Gold competent E. coli (home made)
- Stbl3 competent E. coli (MacroLab, UC Berkeley)
- Growth broth: LB broth
- DMEM medium with GlutaMAX containing 10% FBS and 10% Pen-Strep.

 pVSV-G **addgene Catalog ##138479**


 pCMVR8.74 **addgene Catalog #22036**

 pCDH1-CMV-HT-LC3-SV40-Hygro **addgene Catalog #182045**

 Q5 High-Fidelity DNA Polymerase - 500 units **New England Biolabs Catalog #M0491L**


 NheI-HF - 5,000 units **New England Biolabs Catalog #R3131L**


 BamHI-HF - 50,000 units **New England Biolabs Catalog #R3136L**


 EZ-10 Spin Column DNA Gel Extraction Kit **Bio Basic Inc. Catalog #BS654.SIZE.250preps**

 EZ-10 Spin Column Plasmid DNA Min-preps Kit **Bio Basic Inc. Catalog #BS614.SIZE.250preps**


 Qiagen Hi-Speed MidiPrep kit **Qiagen Catalog #12643**

 DMEM, high glucose, GlutaMAX[®] Supplement **Thermo Fisher Catalog #10566016**

 Gibco[™] Fetal Bovine Serum value heat inactivated (formerly USDA-approved in North America or quali **Fisher Scientific Catalog #A5256801**


 Penicillin-Streptomycin (10,000 U/mL) **Thermo Fisher Scientific Catalog #15140122**

 Opti-MEM[®] I Reduced Serum Medium **Thermo Fisher Catalog #31985070**

 TransIT[®]-LT1 Transfection Reagent **Mirus Bio Catalog #MIR 2300**











 Lenti-X™ Concentrator **Takara Bio Inc. Catalog #631231**

 15-Dimethyl-15-diazaundecamethylene polymethobromide Polybrene **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H9268**

 Hygromycin B (50 mg/mL) **Thermo Fisher Catalog #10687010**

Troubleshooting

Cloning pCDH1-CMV-ATG3-mCherry-Hygro lentiviral vector (ATG3 wild-type or mutants)




- 1 Amplify the coding sequence for human ATG3 (1-125), mCherry or ATG3 (126-314) fragment using Q5 High Fidelity DNA Polymerase and purify PCR products using Gel extraction kit (Bio Basic). 
- 2 PCR products in step 1 are used as templates to run the second step PCR to amplify ATG3 (1- 125)-mCherry-ATG3(126-314) and gel extracted. Primers include a NheI restriction site at the 5' end and a BamHI site at the 3' end.
- 3 Digest the ATG3-mCherry PCR product and pCDH1-CMV-HT-LC3-SV40-Hygro plasmid using NheI and BamHI restriction enzymes (New England BioLabs). 
- 4 Purify the digested DNAs using Gel extraction kit (Bio Basic).
- 5 Ligate the digested ATG3-mCherry and linearized pCDH1-CMV-hygro backbone using T4 DNA ligase (New England BioLabs) according to manufacturer instructions.
- 6 Transform ligation product into XL10-Gold competent *E. coli* and plate on Ampicillin resistant plates. Incubate the plate at  37 °C for  16:00:00 . 16h
- 7 Inoculate one single colony in  5 mL LB medium with ampicillin at  37 °C  Overnight . 16h 
- 8 Miniprep plasmids (Bio Basic) and send for sequencing.



HEK293T lentiviral transfection

- 9 Each plasmid for lentiviral transfection is transformed into Stbl3 competent cells for the propagation. Next day, pick up one colony to inoculate overnight culture in LB medium with ampicillin.
- 10 Midi prep the cultures to purify plasmids (Qiagen).
- 11 Plate 5×10^6 HEK 293T cells on a 10 cm plate in DMEM medium.



12 HEK293T transfection:


12.1 Add retroviral packaging plasmids (pVSV-G, pCMV R8.74) and pCDH1-CMV-ATG3-mCherryHygro (ATG3 wild-type or mutants),  5 μg each in  1.5 mL warm Opti-MEM medium. 

12.2 Add  45 μL of TransIT-LT1 transfection reagent (Mirus) and swirl. 

12.3 Incubate at  Room temperature for  00:15:00 . 

12.4 Add  1.5 mL dropwise into 10 cm HEK293T plate. 


12.5 At 72 hours post-transfection, collect retroviral supernatant into a falcon tube.

13 Concentrate retroviral supernatant to  1 mL using Lenti-X concentrator (Takara) with the manufacturer instruction.

Transduction of HeLa cells with pCDH-CMV-ATG3-mCherry lentivirus

14 Plate 1×10^5 ATG3 KO/HaloTag-LC3B HeLa cells into 12-well plate one day before.

15 Next day, titrate  100 μL ,  200 μL ,  400 μL of concentrated retroviral solution with  8 $\mu\text{g}/\text{mL}$ Polybrene (Sigma) into target HeLa cells.

16 At 24 hours post-transduction, remove retroviral supernatant and replace with fresh DMEM complete medium with  0.5 mg/mL hygromycin B (Gibco).

17 After one week of hygromycin B selection, sort HeLa cells by FACS to enrich for mCherry positive cells using Wolf G2 cell sorter (Nanocollect).