

Apr 15, 2023

Generation of stable cell lines using retroviral system

DOI

dx.doi.org/10.17504/protocols.io.81wgbyez1vpk/v1

nguyen.tha¹

¹Laboratory of Michael Lazarou, Walter and Eliza Hall Institute of Medical Research, Parkville, Victoria, Australia



nguyen.tha

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbyez1vpk/v1>

Protocol Citation: nguyen.tha 2023. Generation of stable cell lines using retroviral system. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbyez1vpk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 09, 2023



Last Modified: April 15, 2023

Protocol Integer ID: 80233

Keywords: stable cell lines, retroviral system, using retroviral system, retroviral system this protocol details generation, generation of stable cell line, stable cell line, cell line

Abstract

This protocol details generation of stable cell lines using retroviral system.

Attachments



[698-1486.docx](#)

19KB

Guidelines

Attention

- The HEK293T cells detach very easily, be extra gentle when changing the media.

Materials

Buffers and reagents:

- Polybrene (4 mg/mL)

Growth media:

	A	B
	DMEM with 10% FBS	
	Glucose	4.5 g/l
	GlutaMAX™	1x
	MEM NEAA	1x
	HEPES	25 mM

⊗ 45% D-()-Glucose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G8769**

⊗ GlutaMAX™; Supplement **Thermo Fisher Catalog #35050061**

⊗ MEM Non-Essential Amino Acids Solution (100X) **Thermo Fisher Scientific Catalog #11140050**

⊗ HEPES Buffer 1M Solution Cell Culture Grade MP Biomedicals **Fisher Scientific Catalog #ICN1688449**

⊗ Lipofectamine™; LTX Reagent with PLUS™; Reagent **Thermo Fisher Catalog #A12621**

⊗ Gibco™ Opti-MEM™ I Reduced Serum Medium no phenol red **Fisher Scientific Catalog #11-058-021**

⊗ Millex-HV Syringe Filter Unit 0.45 µm PVDF 33 mm gamma-sterilizable sterilized **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SLHVM33RS**

Troubleshooting

Safety warnings

Attention

- All viral waste must be bleached and left under UV light for at least 30' after viral work in TC hoods before disposal.



Day 1

- 1 Seed NIH HEK293T cells into a 6-well plate (900k cells/well if set up in the morning, 950k cells/well if set up in the afternoon).




Note








Set up 1 well for each construct you wish to generate a virus harvest for, can be scaled up according to your need.




Day 2: The following protocol is designed for one well of the 6-well plate

- 2 Transfect cells with viral and helper vectors using lipofectamine LTX. Combine the following in a 1.5 mL tube:

A	B
viral vector construct (pBMN, pBABE or pMX) containing cDNA of interest	1.5 µg
gag-pol vector	1.0 µg (amount for 1 well)
VSV-G vector	0.5 µg (amount for 1 well)
Opti-MEM (RT)	500 µL

- 3 Add  3 µL of Plus reagent and mix well. Incubate at  Room temperature for  00:05:00 .


- 4 Add  9 µL of Lipofectamin LTX (1:3 ratio of Plus:LTX is standard in the lab but can be adjusted for your own protocol) and vortex for  00:00:15 . Incubate at  Room temperature for  00:20:00 .


- 5 Once the 20 min incubation starts, replace the media in each well with  1 mL DMEM/10% FBS media.
- 6 When the 20 min incubation finishes, add the optimum/liposome mix to the well.

5m





20m 15s












**Note**

Do it gently on the side of the well.

Day 3

- 7 In the morning, remove the old media from the HEK293T cells which may contain viruses (at this stage) into a beaker of bleach and add  1 mL of fresh growth media. The next day, viruses can be harvested for infection. 
- 8 Seed the target cells (about 100k-120k cells) into a 6-well plate if intending to do infection with fresh viruses.

Day 4

- 9 In the late afternoon, collect viral supernatant from HEK293Ts, spin down at max speed for  00:05:00 to pellet debris and filter through 0.45µm syringe filters. Viral particles can freshly be used for infection on the cells plated out on day 3 (see below) or can be frozen at  -80 °C for future use. 5m
- 10 For second harvest, add  1.5 mL fresh growth media back to HEK293T cells for 2 days and harvest again (on Day 6). 
- 11 For infection, harvested viruses are topped up with fresh growth media to make up a total of  2 mL .
- 12 Aspirate the media from the target cells.
- 13 Add the  2 mL of virus-containing media (from step 3) to the target cells. Add polybrene to a final concentration of  8 µg/mL to the well and mix well.  

Days 5 and 6

- 14 The viruses can be removed from the cells into a beaker of bleach after 24 h (Day 5) or 48 h (Day 6) and fresh media can be added to the wells.



- 15 All waste must be treated as viral waste for at least 3 media changes over 3 days post-infection.