Aug 09, 2023

# Generation of stable cell lines using viral infection

DOI

dx.doi.org/10.17504/protocols.io.261ged38yv47/v1

Louise Uoselis<sup>1</sup>

<sup>1</sup>Lazarou Lab, WEHI

Louise Uoselis





DOI: dx.doi.org/10.17504/protocols.io.261ged38yv47/v1

Protocol Citation: Louise Uoselis 2023. Generation of stable cell lines using viral infection. protocols.io <u>https://dx.doi.org/10.17504/protocols.io.261ged38yv47/v1</u>

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: August 09, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 86206

Keywords: ASAPCRN

#### Abstract

Protocol for generation and precipitation of retrovirus, and infection of HeLa cells to generate stable cell lines.

#### Day 1

1 Seed HEK293T cells into a 10cm tissue culture plate (6.1 million cells/plate), seeding one plate per construct you are generating virus for.

### Day 2

2 Transfect cells with viral and helper vectors using lipofectamine LTX. In a 15mL falcon tube, combine the following:

A	В	С
	Reagent	Amount
1	Viral vector construct containing cDNA of interest	9 ug
2	Gag-pol vector	6 ug
3	VSV-G vector	3 ug
4	Opti-MEM	3 mL

- 3 Add  $\_$  18 µL of PLUS reagent and vortex for ~  $\bigcirc 00:00:20$ . Incubate at 5m 20s  $\blacksquare$  Room temperature for  $\bigcirc 00:05:00$ .
- Once the 20 minute incubation starts, replace the media on the HEK cell plates with
   7 mL media/plate
- 6 At the end of the incubation, add the lipofectamine/OptiMEM mix to the plate using a sterile transfer pipette, in a drop wise fashion across the plate. Gently rock the plate to mix, and return the plate to the incubator.

### Day 3

7 In the morning, remove the media from the HEK cell plate and dispose of the media in a beaker containing bleach. Add  $\boxed{\_6}$  mL of fresh growth media to the HEK cell plate.

#### Day 4

- 8 In the late afternoon, collect the media from the HEK cell plate into a falcon tube. Add back <u>A</u> 8 mL fresh growth media to the HEK cell plate.
- 9 Centrifuge the media collected at maximum speed for O0:05:00 at
  If Room temperature to pellet any debris. Pass the media through a 0.45 mm syringe filter into a clean falcon tube. Parafilm the falcon tube, and place store the media at
  If 4 °C in the fridge.

# Day 6

- In the morning, collect the media from the HEK cell plate into a falcon tube. Centrifuge
   the collected media at maximum speed for 00:05:00 at Room temperature.
   Pass the media through a 0.45 mm syringe filter into a clean falcon tube.
- 11 Combine the media collected in step 9 and in step 10 into one tube. Add Lentivirus Precipitation Solution (ALSTEM) to the media at a ratio of 1:4 (Lentivirus precipitation solution : viral media).
- 12
   Gently invert the tube 8 10 times to mix the precipitation solution and media together.
   2d

   Parafilm the tube and leave on a rocker gently rocking at 4 °C for 348:00:00
   2d

# Day 7

13 If intending to infect cells with fresh virus, seed the cells to be infected into a 6 well plate (~220K cells/well).

### Day 8

1h

5m

- 15 Remove the supernatant, disposing of it into a beaker containing bleach. Resuspend the precipitated pellet in 3 6 mL of fresh growth media (3 mL if you require high construct expression).
- 16 Replace the media on the cells to be infected with  $\angle 1.5 \text{ mL}$  fresh growth media/well. Add  $\angle 500 \text{ µL}$  of the precipitated virus to each well to be infected.
- 17 Add polybrene at a final concentration of 8 mg/mL to each well, and rock the plate to mix. Place the cells back into the incubator.
- 18 Any leftover virus can be frozen at **3** -80 °C for use at a future date

# Day 9 - 10

19 After 24 – 48h (depending on the construct and expression required), remove the media from the infected cells and place in bleach to dispose. Add fresh growth media back to the cells.

#### Day 15 - 17

20 After 5 - 7 days after the virus has been removed (and after three passages), using FACS, sort the cells populations to the desired fluorescent expression level.

1h