



Sep 23, 2023

Generation of stable cell lines via lentiviral transduction

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr3e5pvmk/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr3e5pvmk/v1>

Protocol Citation: Elias Adriaenssens 2023. Generation of stable cell lines via lentiviral transduction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.6qpvr3e5pvmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 26, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 84010

Keywords: ASAPCRN, lentiviral transduction, lentivirus system, generation of stable cell line, lentivirus, stable cell line, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

Here, we describe a protocol to generate stable cell lines using a lentivirus system. Please note that necessary safety measures are to be taken in working with lentivirus.

Attachments



[751-1919.pdf](#)

120KB



Materials

MATERIALS

Cell lines

- HEK293T for virus packaging and propagation
- HeLa cells

Plasmids

- Gag/Pol plasmid
- VSV-G plasmid
- Lentiviral vectors (pHAGE-FKBP-GFP-GOI)

Note

Note: We purify plasmids using a QIAGEN Plasmid Maxi kit following the manufacturer's protocols and ensure sterile reagents are used and mixtures prepared in tissue culture hood to avoid contamination.

Media and Reagents

- DMEM
- 10% FBS
- 1% Penicillin-Streptomycin
- 1% non-essential amino acids
- 25 mM HEPES

Transfection media (for HEK293T cells)

- Opti-MEM I Reduced Serum Medium (Gibco)
- Lipofectamine 3000 (ThermoFisher) or PEI Max (MW 40000, Polysciences)
- Polybrene (Sigma)

Safety warnings

⚠ Please note that necessary safety measures must be taken to work with lentivirus.



Packaging lentiviral plasmid into a lentiviral particles for in...

2d 0h 20m

1 Grow HEK293T cells to 60-70% confluency in Growth media in a 6-well Petri Dish.

2 Prepare a transfection mix in a sterile 1.5 ml Eppendorf tube, containing:

A	B
Lentiviral vector with your gene-of-interest	1500 ng
Gag/Pol plasmid	1000 ng
VSV-G plasmid	500 ng
P3000 reagent (Lipofectamine 3000 kit)	5 μ L
OptiMem	125 μ L

3 Prepare Lipofectamine 3000 mixture in a sterile 1.5 ml Eppendorf tube, containing:

A	B
Lipofectamine 3000	5 μ L
OptiMem	125 μ L

4 Incubate each mixture (from steps 2 and 3) separately for ~  00:05:00 at  Room temperature .

5m



5 Mix the two suspensions and incubate at  Room temperature for  00:15:00 .

15m



6 Add the mixture drop-wise to the cells from step 1 using a P1000 sterile pipette.



7 Incubate cells at  37 °C for  24:00:00 .

1d





- 8 Collect the supernatant from the cells (that now contains the lentiviruses) and pass it through a  0.45 μm syringe filter. If needed, add fresh growth medium and collect this too  24:00:00 later.

1d

Lentiviral infection of HeLa cells

1d

- 9 Add  2 μL of polybrene (8 mg/ml) to  2 mL of lentivirus infection media from step 8 to HeLa cells seeded in 6-well plate (at 700.000/well) at 60% confluence.



- 10 Incubate at  37 $^{\circ}\text{C}$ for  24:00:00 .

1d



- 11 Change media to fresh medium and incubate until confluency. Split three times before cells can leave the viral S2 lab.



- 12 Cells can now be passaged and plated for experiments or frozen down for long-term storage in liquid nitrogen.

Note

Freezing media: growth media added with 20% FBS and 10% v/v DMSO.