

Aug 05, 2025

🌐 GETCampy Sequencing protocols

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw4jpwlmk/v1>



Ben Pascoe¹, Henry Badji², Madison Goforth¹, Evangelos Mourkas³, Evariste Bako⁴, Isidore Bonkougou⁴, Modeste T Gampene⁴, M Edith M Nikiema⁴, Barthélemy Zoma⁴, Ousman E Bah², Eustacia J Cassell², Abdoulie K Ceesay², Bubcarr E Ceesay², Ousman Ceesay², Bakary Conteh², Lamin Drammeh², Binta Faye², Abdoulie F Jallow², Fatima Jallow², Ousman Jallow², Samba Juma Jallow², Pa Modou Joof², Modou Kandeh², Mehrab Karim², Jarra Manneh², Abdou Ceesay⁵, Ebrima Fofana⁵, Mamud Jallow⁵, Modou B Jarju⁵, Samba Ba⁶, Dodou Sanyang⁶, Demba B Jallow⁷, Mukaila I Alebiosu⁸, Akosua Bonsu Karikari⁸, Atanyiwoen Brusah⁹, Courage KS Saba⁹, Kaisa Haukka¹⁰, Catherine Matilda Collins¹¹, Adrian W Leach¹¹, Polina Levontin¹¹, Shani UP Ali¹², Frances M Colles¹³, Matthew Hitchings¹³, Martin Antonio¹, Abdul Karim Sesay¹, Ozan Gundogdu¹⁴, Jahangir Hossain², Samuel K Sheppard¹, Bubacarr E Ceesay¹⁵, Bubacarr E. Ceesay¹⁶

¹Ineos Oxford Institute for Antimicrobial Research, Department of Biology, University of Oxford, Oxford, United Kingdom;

²Medical Research Council (UK) Unit The Gambia at the London School of Hygiene and Tropical Medicine, Banjul, The Gambia;

³Zoonosis Science Center, Department of Medical Sciences, Uppsala University, Uppsala, Sweden;

⁴Laboratory of Molecular Biology, Epidemiology and Monitoring of Bacteria and Virus Transmitted by Food (LaBESTA), University Joseph KI-ZERBO, 03 BP 7021 Ouagadougou 03, Burkina Faso;

⁵Department of Livestock Services, Banjul, The Gambia;

⁶Regional health Directorate, Upper River Region, Ministry of Health, The Gambia;

⁷National Agricultural Research Institute, The Gambia;

⁸Department of Clinical Microbiology, School of Medicine, University for Development Studies, Tamale, Ghana;

⁹Faculty of Biosciences, Department of Biotechnology. University for Development Studies, Tamale, Ghana;

¹⁰Department of Microbiology, University of Helsinki, Helsinki, Finland;

¹¹Centre for Environmental Policy, Imperial College London, The Weeks Building, 16-18 Princes Gardens, London SW7 1NE, United Kingdom;

¹²Milner Centre for Evolution, Department of Life Sciences, University of Bath, Bath, UK;

¹³Department of Biology, University of Oxford, Oxford, United Kingdom;

¹⁴Department of Infection Biology, Faculty of Infectious & Tropical Diseases, London School of Hygiene and Tropical Medicine;

¹⁵Medical Research Council at the London School of Hygiene and Tropical Medicine;

¹⁶Medical Research Council Unit, The Gambia, at London School of Hygiene and Tropical Medicine

GETcampy



Ben Pascoe

Ineos Oxford Institute for Antimicrobial Research, Departmen...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw4jpwlmk/v1>

External link: <https://sheppardlab.wordpress.com/projects/get-campy-africa/>

Protocol Citation: Ben Pascoe, Henry Badji, Madison Goforth, Evangelos Mourkas, Evariste Bako, Isidore Bonkougou, Modeste T Gampene, M Edith M Nikiema, Barthélemy Zoma, Ousman E Bah, Eustacia J Cassell, Abdoulie K Ceesay, Bubcarr E Ceesay, Ousman Ceesay, Bakary Conteh, Lamin Drammeh, Binta Faye, Abdoulie F Jallow, Fatima Jallow, Ousman Jallow, Samba Juma Jallow, Pa Modou Joof, Modou Kandeh, Mehrab Karim, Jarra Manneh, Abdou Ceesay, Ebrima Fofana, Mamud Jallow, Modou B Jarju, Samba Ba, Dodou Sanyang, Demba B Jallow, Mukaila I Alebiosu, Akosua Bonsu Karikari, Atanyiwoen Brusah, Courage KS Saba, Kaisa Haukka, Catherine Matilda Collins, Adrian W Leach, Polina Levontin, Shani UP Ali, Frances M Colles, Matthew Hitchings, Martin Antonio, Abdul Karim Sesay, Ozan Gundogdu, Jahangir Hossain, Samuel K Sheppard, Bubacarr E Ceesay, Bubacarr E. Ceesay 2025. GETCampy Sequencing protocols. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvw4jpwlmk/v1>

Manuscript citation:

TBC

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: August 05, 2025



Last Modified: May 18, 2026

Protocol Integer ID: 224086

Keywords: Campylobacter, gastroenteritis, genomics, pathogen surveillance, Whole genome sequencing, 16S sequencing, metagenomics sequencing, genomic surveillance of campylobacter, other enteric pathogens across africa, getcampy sequencing protocol, getcampy sequencing protocols as part, enabling hybrid metagenomic approach, analysing enteric pathogen, assembled genome, hybrid metagenomic approach, enteric pathogen, resolution genome assembly, other enteric pathogen, genomic surveillance, sequencing method, campylobacter, dna extraction protocol, taxonomic profiling, total community dna, oxford nanopore, high molecular weight recovery,

Funders Acknowledgements:

Medical Research Council (UK)

Grant ID: MR/V001213/2

Disclaimer

Protocols for research purposes only

Abstract

As part of the GETCampy-Africa project, this protocol provides a modular workflow for extracting, sequencing, and analysing enteric pathogens from stool and environmental samples. It includes short-read (Illumina MiSeq) and long-read (Oxford Nanopore) sequencing methods, enabling hybrid metagenomic approaches for high-resolution genome assembly. DNA extraction protocols support both total community DNA and high molecular weight recovery. Library preparation workflows include Nextera XT (Illumina), rapid barcoding and ligation kits (ONT), and full-length 16S rRNA sequencing for taxonomic profiling. The protocol also outlines hybrid assembly, binning, and quality control steps for metagenome-assembled genomes (MAGs). Designed for field-lab compatibility and resource-limited settings, these methods support genomic surveillance of *Campylobacter* and other enteric pathogens across Africa.

Image Attribution

CCBY c/o Ben Pascoe

DNA extraction

1 QIAamp DNA Extraction Protocol for *Campylobacter*

Kit: QIAGEN QIAamp DNA Mini Kit

Adapted from: QIAamp DNA Mini Handbook, Appendix D (page 55) and Step 3 (page 33)

Organism: *Campylobacter* spp.

1.1 Sample Preparation

1. Grow *Campylobacter* isolates:

- **Medium:** Chocolate, blood, or Brucella agar
- **Conditions:** 42°C for **48 hours**

Using a sterile loop, harvest **at least one loop-full** of bacterial growth.

Resuspend cells in **180 µL Buffer ATL**.

Incubate at **37°C for 30 minutes**.

1.2 Cell Lysis and Protein Digestion

Add **20 µL Proteinase K** (from the kit).

Mix thoroughly by vortexing.

Incubate at **56°C for 30–60 minutes**.

Briefly spin down to remove condensation from the tube lid.

1.3 RNase Treatment and Lysis Enhancement

Add:

- **4 µL RNase** (100 mg/mL)
- **200 µL Buffer AL** (from the kit)

Vortex to mix completely.

Incubate at **70°C for 10 minutes**.

1.4 DNA Binding

Add **200 µL of 95–100% ethanol** to the lysate.

Vortex to mix.

Transfer the full volume (~600 µL) into a **QIAamp spin column** in a **2 mL collection tube**.

Centrifuge at **8000 rpm for 1 minute** at **room temperature (RT)**.

Discard the flow-through and move the column to a clean 2 mL collection tube.

1.5 Wash Steps

Add **500 µL Buffer AW1** to the column.

Centrifuge at **8000 rpm for 1 minute** at RT.

Discard flow-through.

Add **500 µL Buffer AW2**.

Centrifuge at **14,000 rpm for 3 minutes** to dry the membrane.

1.6 DNA Elution

Place the QIAamp column into a **clean 1.5 mL microcentrifuge tube**.

Add **50 µL Buffer AE** (or nuclease-free water) directly to the column membrane.



Let sit at **room temperature for 1 minute**.
Centrifuge at **8000 rpm for 1 minute** to elute DNA.

1.7 **Storage**

- **Short-term:** Store at **4°C**
- **Long-term:** Store at **-20°C**

1.8 **Notes**

- Ensure all reagents from the kit are at working temperature.
- Eluted DNA is suitable for downstream applications including PCR and sequencing.

Library preparation for short read sequencingUntitled section

2 **Illumina Custom Protocol (Nextera XT on MiSeq)**

Instrument: MiSeq

Library Prep: Nextera XT DNA

Reagent Kit: MiSeq v3

Indexing: Dual Indexing

Software: BaseSpace Sequence Hub (FASTQ generation)

2.1 **Generate Sample Sheet with IEM**

Start IEM Software

- Start → All Programs → Illumina → Illumina Experiment Manager or click IEM icon.

Create Sample Plate

Select **Create Sample Plate** → Choose **Nextera XT**.

Set:

- Unique plate name
- Index reads = 2

Enter sample details per well (ID, i7/i5 index, etc.)

Save as *.nexxt28.plt

2.2 **Create Sample Sheet**

Select **Create Sample Sheet** → Instrument: **MiSeq**.

Application: **FASTQ Only** → Click **Next**.

Set:

- Reagent barcode
- Prep kit: Nextera XT
- Index reads: 2
- Experiment name, date, paired/single-end, cycles

Load plate, select samples, and save sample sheet.

2.3 **Nextera XT Library Prep**

Tagment Genomic DNA

- Thermal cycler: **55°C for 5 min**, hold at **10°C**
- Mix in PCR plate:
- 10 µL TD + 5 µL gDNA (0.2 ng/µL)
- Add 5 µL ATM, mix
- Spin 280 × g, 1 min → Run tagmentation
- Add 5 µL NT, mix → Spin → RT 5 min

2.4 **Amplify Libraries**

- Thermal cycler:
- 72°C 3 min → 95°C 30s
- 12 cycles: 95°C 10s → 55°C 30s → 72°C 30s
- Final: 72°C 5 min → hold at 10°C
- Add i7 and i5 indices (5 µL each), then 15 µL NPM
- Spin → Run PCR

Safe stop: 2–8°C for 2 days or overnight on thermal cycler

2.5 **Clean Up Libraries**

Transfer 50 µL PCR product to new plate.

Add 30 µL AMPure XP beads → Shake 1800 rpm 2 min → RT 5 min

Magnetic stand → remove supernatant.

Wash ×2 with 200 µL 80% EtOH → Remove residual EtOH

Air dry 15 min.

Elute in 52.5 µL RSB → Shake, incubate, magnet, transfer 50 µL to new plate

Safe stop: -25°C to -15°C for 7 days

2.6 **Check Libraries (Optional)**

- Run 1 µL on Agilent Bioanalyzer (HS DNA chip)

2.7 **Normalize Libraries**

Transfer 20 µL to new plate.

Mix 4.4 mL LNA1 + 800 µL LNB1 → Add 45 µL per well

Shake 1800 rpm 30 min → Magnet → Remove supernatant

Wash ×2 with 45 µL LNW1

Add 30 µL 0.1 N NaOH → Shake 5 min

Add 30 µL LNS1 to new plate (SGP)

Resuspend beads → Shake 5 min → Magnet → Transfer to SGP plate

Spin at 1000 × g for 1 min

Safe stop: -25°C to -15°C for 7 days

2.8 **Pool Libraries**

Spin SGP plate.

Transfer 5 µL per sample to PCR tube strip → Pool into tube (PAL)

Dilute to loading concentration.

Store PAL and SGP: -25°C to -15°C (7 days)



2.9 Denature and Dilute Libraries

Prepare HT1

- Thaw at RT, store at 2–8°C.

2.10 Dilute Library

	Library Volume	HT1 Volume
	6 µL	594 µL
	24 µL (XT)	576 µL

- Vortex → Spin → Incubate 98°C 2 min → Ice 5 min

Illumina short read sequencing

3 C-SOP-601: Illumina MiSeq Operation for Whole Genome Sequencing (WGS)

Preparation Before Starting

Clean workspace with 80% ethanol.

Wear PPE and follow local GLP guidelines.

Prepare ice bucket for temporary reagent/sample storage.

Thaw HT1 buffer at room temperature, then store at 4°C until use.

Prepare 0.2N NaOH:

- Mix 200 µL 1.0N NaOH + 800 µL nuclease-free water.
- Use within 8 hours.

Reagent Cartridge Thawing (2 Hours Prior to Run)

Remove cartridge from -20°C.

Submerge base only in room-temp water bath (max line indicated).

Thaw time:

- v3 kits: ~1.5 hours
- v2 kits: ~1 hour

Dry base and invert 10x to mix reagents.

Check positions 1, 2, and 4 for full thawing (no precipitates).

Store at 4°C until ready to load.

MiSeq Instrument Initialisation

Power on MiSeq (switch on back).

Log in to BaseSpace, if networked.

Wait for MiSeq Control Software (MCS) to initialise.

Denature and Dilute DNA Library

Mix:

- 5 µL 4nM DNA library
- 5 µL 0.2N NaOH

Vortex 10 sec, spin 1 min.

Incubate 5 min at RT.

Add 990 µL chilled HT1 → final 20pM denatured library.

Dilute to final loading concentration (e.g., 10pM).

Denature and Dilute PhiX (Optional Spike-in)

Dilute 2 µL 10nM PhiX + 3 µL Tris-Cl/Tween → 4nM.

Mix:

- 5 µL 4nM PhiX
- 5 µL 0.2N NaOH

Vortex 10 sec, spin 1 min, incubate 5 min.

Add 990 µL HT1 → final 20pM PhiX.

For v2 kits: mix 375 µL PhiX + 225 µL HT1 → 12.5pM.

Add 1–5% PhiX to final sample mix depending on library complexity.



MiSeq Run Setup

Ensure ≥100 GB free disk space.

Enable local analysis replication:

- MCS → Run Options → check "Replicate analysis locally".

Prepare and load Sample Sheet via Illumina Experiment Manager (IEM).



Prepare and Load Flow Cell

Retrieve Flow Cell + PR2 buffer from 4°C.

Rinse flow cell with lab-grade water → dry with lint-free lens paper.

Clean glass with ethanol-wetted lens tissue.

Install flow cell into clamp → wait for RFID recognition.



Load Incorporation Buffer (PR2) and Empty Waste

Raise sipper handle.

Load mixed PR2 buffer bottle → confirm RFID.

Empty and replace waste bottle.

Lower sipper into PR2 + waste bottle.



Load DNA Library into Cartridge

Tap cartridge to gather contents.

Pierce foil at "Load Sample".

Pipette 600 µL denatured library into well.

Tap gently to collect sample at base.



Load Cartridge into MiSeq

Wait for "Load Reagent Cartridge" prompt.

Remove wash tray, dry chiller base if needed.
Insert cartridge (Illumina label facing out) until fully seated.
Confirm cartridge RFID is read.



Sample Sheet Confirmation

MCS looks for matching barcode sample sheet.
If not found, manually browse and select.
Click "Save and Continue" → "Next".



Start Sequencing Run

Review:

- Experiment Name
- Workflow
- Read Length

Perform Pre-run Check.

When all checks pass, click "Start Run".



Post-Run Wash (Required)

Follow on-screen prompts to load:

- Wash tray
- Wash bottle
- Optional: MiSeq wash tube (0.01% NaOCl for bleach wash)

Start Wash.

Do not remove flow cell/wash tray during wash.

Leave wash solution in tray/bottle after wash completes.



Monitor Run

- Use Sequencing screen (view-only) or BaseSpace/SAV.
- Monitor: cycles, cluster intensities, quality scores.



Export Data

- Navigate to:
`D:\Illumina\MiSeq Output\Run Folder\Data\Intensities\BaseCalls`
- Download `fastq.gz` files from R1 and R2.
- Or use BaseSpace/local network share.



Instrument Shutdown or Reboot

- Weekly reboot recommended
- MCS → Manage Instrument → Reboot
- For shutdown:
- MCS → Manage Instrument → Shutdown
- Power switch off → wait 10 min → power on



Monthly Maintenance

Week | Action

- 1 PTL (post-run bleach)
 - 2 PTL
 - 3 PTL
 - 4 PTL or Maintenance
- See MiSeq System Guide for detailed cleaning steps.



Troubleshooting Resources

- MiSeq System Guide
- Sample Sheet Quick Guide
- MiSeq Reporter / SAV Software Guide
- Re-Queue Tutorials

High molecular weight DNA extraction

4 C-SOP-1001: HMW DNA Extraction from Bacteria Using NEB Monarch HMW DNA Extraction Kit

Materials Required (Not Supplied)

- Microcentrifuge
- Thermal mixer (1.5 mL + 2 mL blocks)
- Vertical rotating mixer (e.g., HulaMixer)
- Ethanol ≥95%
- Cold PBS or TE/Tris
- Isopropanol (550 µL/sample or 275 µL low input)
- 1.5 mL DNase-free low DNA binding tubes
- Wide-bore pipette tips
- Lysozyme (Gram-neg: 25 mg/mL; Gram-pos: 10 mg/mL in STET)
- Additional lysis agents (e.g., lysostaphin if needed)

Preparation

Clean bench with 80% ethanol. Wear PPE.

Preheat thermal mixers to 37°C and 56°C.

Add ethanol to gDNA Wash Buffer as per label.

Chill PBS/TE/Tris buffer.

Label tubes. Store Proteinase K and RNase A at -20°C.

Step 1: Cell Pellet

- Centrifuge bacterial culture at **12,000 × g for 1 min** in a Monarch Pestle Tube.

Step 2: Lysis

A. Gram-negative Bacteria



Resuspend pellet in **300 µL cold PBS** (150 µL for low input).

Add **10 µL lysozyme (25 mg/mL)**, vortex briefly.

Add **300 µL HMW Lysis Buffer**, invert 5–10×.

Incubate at **37°C** with agitation:

- **2000 rpm** (faster lysis, shorter fragments)
- **500 rpm** (slower lysis, longer fragments)
- Until lysate turns clear (~5 min)

B. Gram-positive Bacteria

Resuspend in **STET buffer + lysozyme (10 mg/mL)**.

Incubate **30 min at 37°C (no agitation)**.

Add **300 µL HMW Lysis Buffer**, invert 5–10×.

Step 3: Proteinase K Digestion

Increase thermal mixer to **56°C** if needed.

Add **20 µL Proteinase K** (10 µL for low input), invert 10–20×.

Choose Homogenization:

- **Thermal Mixer:** Incubate **30 min at 56°C**, 2000 rpm.
- **Rotor-stator:** Homogenize 5–15 sec, then incubate **30 min at 56°C**.

Step 4: RNA Removal

Add **10 µL RNase A** (5 µL for low input), invert 5–10×.

Incubate **10 min at 56°C**, same speed as lysis.

Step 5: Protein Separation

Add **300 µL Protein Separation Solution**, invert for 1 min or mix at 20 rpm.

Centrifuge **10 min at 16,000 × g**.

Carefully transfer **upper DNA phase (~800 µL)** to new 2 mL tube.

Step 6: DNA Binding to Beads

Add **2 DNA Capture Beads** to tube.

Add **550 µL isopropanol** (275 µL low input).

Mix gently:

- Rotate at 10 rpm for 5 min OR
- Invert 25–30× slowly (5–6 sec per inversion)
- Watch for visible DNA wrapping around beads

Step 7: Wash Steps

Remove liquid carefully (avoid touching beads).

Wash with **500 µL gDNA Wash Buffer**, invert 2–3×, then remove.

Repeat wash.



Step 8: Bead Transfer

Transfer beads to **Monarch Collection Tube II + bead retainer**.

Pulse spin (≤ 1 sec).

Move beads to a fresh **2 mL tube**.

Insert bead retainer into **1.5 mL DNA low-bind tube**.

Step 9: Elution

Add **100 μ L Elution Buffer II** (≥ 5 min at 56°C, 300 rpm).

- Optional: use **50 μ L** with more frequent gentle shaking

Pour eluate + beads into bead retainer in 1.5 mL tube.

Step 10: Final Handling

Centrifuge **30 sec at 12,000 $\times$ g**.

Discard beads and retainer.

Pipette eluate 5–10 $\times$ with wide-bore tips to disperse aggregates.

Dissolve DNA:

- **1 h at 37°C**, or
- **24 h at RT**, or
- **Overnight at 4°C**

Pipette again before use.

Storage

- Short-term: **4°C**
- Long-term: **-20°C** (in Elution Buffer: 10 mM Tris, pH 9.0, 0.5 mM EDTA)

Additional Resources

- [UHMW DNA Extraction for Ultra-Long Read Nanopore](#)
- [UHMW DNA Cleanup Protocol](#)

Rapid barcoding for long read sequencing

5 **Rapid Barcoding Workflow**

Kit: ONT Rapid Barcoding Kit 96 (SQK-RBK114.96)

Platform: MinION / GridION / PromethION

Purpose: For barcoded, multiplexed sequencing of DNA using ONT's transposase-based tagmentation

5.1 **DNA Input Requirements**

- **Input amount:** 10–50 ng per sample (in ≤ 2.5 μ L nuclease-free water)
- **Volume of RB reagent per sample:** 2.5 μ L
- **Volume of barcode (BCXX) per sample:** 0.5 μ L
- **Multiplexing:** Up to 96 samples



5.2 Barcoding Reaction

Prepare Reaction Mix

In a PCR strip tube or plate (per sample):

Reagent	Volume (μL)
Genomic DNA	≤2.5
RB (Fragmentation Mix)	2.5
Barcode (BCXX)	0.5
Total	5.5 max

5.3 Mix and Incubate

- Gently flick or pipette to mix.
- Incubate at **30°C for 1 minute**, then **80°C for 1 minute**.
- Immediately place on ice.

5.4 Pooling and Cleanup

Pool 2 μL from each barcoded sample into a single tube.

Total volume should not exceed 10 μL.

Proceed to magnetic bead cleanup:

- Add 1.8× volume of AMPure XP beads.
- Incubate 5 minutes at RT.
- Place on magnet for 5 minutes.
- Wash twice with 70% ethanol.
- Elute in 10 μL Elution Buffer or nuclease-free water.

5.5 Adapter Ligation

Reagent	Volume (μL)
Barcoded DNA	10.0
RAP (Rapid Adapter)	1.0
Total	11.0

- Mix by flicking or pipetting.
- Incubate at **room temperature for 5 minutes**.
- Proceed immediately to loading.

5.6 Loading onto Flow Cell

- Load the library onto a primed flow cell following the Flow Cell Priming protocol.



- Use appropriate loading buffer and loading volume as per the platform used (MinION/GridION/PromethION).
- Avoid bubbles when loading.

5.7 Storage and Stability

- **Barcoded samples:** Use immediately or store at -20°C (short-term).
- **Final library:** Load immediately. Do not freeze post-adaptor ligation.

Long read sequencing using ONT

6 Long Read Sequencing Protocol using Oxford Nanopore Technologies

Workflow: Ligation Sequencing Kit (LSK114)

Platform: MinION / GridION / PromethION

Purpose: For high-yield, long-read sequencing using ONT flow cells

6.1 Input DNA Requirements

- **Amount:** ≥ 1 μg high molecular weight DNA
- **Concentration:** ≥ 50 ng/ μL in ≥ 20 μL volume
- **Buffer:** Nuclease-free water or low-salt TE (avoid EDTA > 0.1 mM)
- **Fragmentation:** Avoid shearing to preserve read length

6.2 DNA Repair & End-Prep

Reagent	Volume (μL)
DNA sample	Up to 48.0
NEBNext FFPE DNA Repair Buffer	7.0
NEBNext FFPE DNA Repair Mix	3.0
NEBNext Ultra II End Prep Reaction Buffer	3.5
NEBNext Ultra II End Prep Enzyme Mix	3.0
Total	~64.5

Mix gently, spin briefly.

Incubate:

- **20°C for 5 min**
- **65°C for 5 min**

Place on ice.

6.3 AMPure XP Bead Cleanup

Add 1 $\times$ volume (e.g., 65 μL) of AMPure XP beads.

Incubate 5 min at RT.

Magnet for 5 min → remove supernatant.



Wash twice with 200 μL 70% ethanol.
Air dry beads for 2 min (avoid over-drying).
Elute DNA in 61 μL nuclease-free water.

6.4 Adapter Ligation

Reagent	Volume (μL)
Eluted DNA	60.0
NEBNext Quick T4 DNA Ligase Buffer (5X)	25.0
Adapter Mix (AMX)	10.0
Quick T4 DNA Ligase	5.0
Total	100.0

Mix by gentle flicking.
Incubate **room temperature for 10–15 minutes**.

6.5 Post-Ligation Cleanup

Add 0.4 $\times$ volume (~ 40 μL) of AMPure XP beads.
Incubate 5 min at RT.
Magnet for 5 min $\rightarrow$ discard supernatant.
Wash twice with 250 μL Long Fragment Buffer (LFB).
Air dry beads (≤ 2 min).
Elute final library in **15 μL Elution Buffer (EB)**.

6.6 Flow Cell Priming (MinION/GridION)

Thaw and mix:

- Flush Buffer (FB)
- Flush Tether (FLT) $\rightarrow$ mix with FB to make priming mix

Open SpotON port and priming port.
Remove ~ 20 μL from priming port to release air.
Slowly inject **800 μL priming mix** via priming port.

6.7 Load Library

Prepare loading mix:

- 37.5 μL Sequencing Buffer (SQB)
- 25.5 μL Loading Beads (LB, mix just before use)
- 12 μL Final library (from step 5)

Mix gently by pipetting.
Load **75 μL** into the SpotON port drop-by-drop.
Close ports and lid.

6.8 Start Sequencing Run

Use **MinKNOW** software.

Choose appropriate flow cell and kit settings (e.g., LSK114).

Set run time (e.g., 24–72 hours depending on goal).

Start sequencing and monitor real-time QC metrics.

6.9 After the Run

- Optionally wash and reuse flow cells using **Flow Cell Wash Kit** (up to 2–3 times).
- Store remaining library at **-20°C** for short-term re-use (avoid freeze-thaw).

6.10 Optional Tips

- For ultra-long reads, minimize pipetting and vortexing.
- Use wide-bore or cut tips during extraction and library prep.
- Consider using Qubit and Femto Pulse for DNA QC and sizing.

16S sequencing for sample taxonomy profiling

7 16S rRNA Gene Amplicon Sequencing (ONT Workflow)

Goal: Taxonomic profiling via full-length 16S rRNA gene sequencing

Platform: MinION / GridION / PromethION

Kit: 16S Barcoding Kit (SQK-16S024 or similar)

Read Type: Long-read, full-length 16S (~1.5 kb)

7.1 Sample Preparation

- Extract total DNA from microbial or environmental samples.
- Use a high-fidelity, contamination-free method (e.g., Qiagen PowerSoil for stool).
- Quantify DNA using Qubit (dsDNA HS).
- Normalize DNA to **10–50 ng per PCR reaction**.

7.2 16S rRNA Gene PCR Amplification

Reagent	Volume (μL)
Template DNA (10–50 ng)	~2–5
Barcoded 16S primers (forward + reverse)	1.0
LongAmp Taq 2X Master Mix (NEB)	25.0
Nuclease-free water	to 50.0

Cycling Conditions (for full-length 16S):

- 95°C for 3 min
- 25 cycles of:
 - 95°C for 20 sec
 - 55°C for 30 sec
 - 65°C for 2 min

- Final: 65°C for 5 min → hold at 4°C
- ✓ Tip: Use primers covering full-length 16S gene (e.g., 27F/1492R).

7.3 Clean Up PCR Products

- Use **AMPure XP beads (1× ratio)**:
- Mix, incubate 5 min RT
- Magnet, wash ×2 with 70% EtOH
- Elute in 10 µL nuclease-free water

7.4 Pooling and Quantification

- Pool 10 µL of each barcoded sample.
- Quantify final pool via Qubit.
- Aim for **50–200 fmol total DNA** (up to 1 µg).

7.5 Adapter Ligation (if not pre-tagged)

Reagent	Volume (µL)
Pooled amplicons	~10
Rapid Adapter (RAP) or AMX	1.0
Incubate	5 min at RT

7.6 Load Library onto Flow Cell

Prime the flow cell (see ONT Flow Cell Priming SOP).

Prepare the loading mix:

- ~12 µL DNA library
- 25.5 µL Loading Beads (LB)
- 37.5 µL Sequencing Buffer (SQB)

Mix gently and load **75 µL** dropwise onto SpotON port.

7.7 Start Sequencing

- Use **MinKNOW** and select:
- Kit: SQK-16S024 (or latest)
- Flow Cell: FLO-MIN106D (R10.4.1 recommended)
- Duration: 12–48 hours depending on depth

7.8 Data Analysis

- Basecalling via **Guppy** or **MinKNOW**
- Taxonomic classification:
- Use **EPI2ME Fastq 16S Workflow**, or
- Tools like **QIIME 2**, **Kraken2**, or **NanoCLUST**
- Generate bar plot of relative abundances per sample

7.9 Storage

- Store final amplicon libraries at **-20°C** short-term.
- Flow cells can be washed and reused (1–2 times).

7.10 **Notes**

- Full-length 16S sequencing provides species-level resolution.
- Minimize chimera formation by using high-fidelity polymerases.
- Always run negative controls to detect contamination.

Hybrid sequencing for metagenomes

8 **Hybrid Metagenomic Sequencing Protocol (Illumina + ONT)**

Sample Collection & DNA Extraction

Collect samples (e.g. stool, soil, rumen, etc.) under sterile conditions.

Store samples at -80°C if not processing within 1 hour.

Extract high molecular weight (HMW) DNA using a suitable kit (e.g. QIAGEN MagAttract HMW DNA Kit).

- Minimise mechanical shearing.
- Aim for fragment sizes >20 kb for ONT compatibility.

Quantify DNA with Qubit and assess quality on TapeStation or agarose gel.

8.1 **Library Preparation**

Illumina (Short-Read) Library Prep

- Input: ~ 100 ng DNA.

Recommended kits:

- Illumina Nextera XT, or
- NEBNext Ultra II FS.

Fragmentation: ~ 300 – 500 bp.

Add adapters and indices.

Clean up with AMPure XP beads.

Validate size distribution with Bioanalyzer or similar.

8.2 **ONT (Long-Read) Library Prep**

- Input: ~ 1 – 2 μg of high-quality HMW DNA.

Recommended kits:

- Ligation Sequencing Kit (SQK-LSK114) – highest accuracy.
- OR Rapid Sequencing Kit (SQK-RBK114.96) – faster prep.

Perform end-repair/dA-tailing if needed.

Ligate adapters.

Clean up and elute DNA for sequencing.

8.3 **Sequencing**

Illumina

- Platform: HiSeq, NovaSeq or NextSeq.
- Paired-end 150 bp reads.
- Target depth: ~ 20 – 40 Gbp per sample.

ONT

- Platform: MinION, GridION, or PromethION.
- Flowcell: R10.4.1 or latest.
- Target depth: ~5–10 Gbp per sample.
- Monitor in real-time, reload flow cell if necessary.

8.4 D. Data Processing & Hybrid Assembly**Quality Control**

- **ONT:** Basecall using Guppy (super accuracy), filter reads < Q10.
- **Illumina:** Trim adapters and low-quality bases using fastp or TrimGalore

Assembly**Short-read only:**

- metaSPAdes
- MEGAHIT

Long-read only:

- Flye
- Raven

Hybrid assembly:

- metaFlye
- HybridSPAdes
- OPREA-MS

Polishing

Use

- Racon (long reads)
- Pilon (short-read) for iterative polishing.

8.5 Binning & Genome Reconstruction**Contig Binning**

Tools:

- MetaBAT2
- MaxBin2
- VAMB
- SemiBin

Inputs: k-mer composition and coverage profiles.

Bin Refinement

Merge/refine bins with

- DAS Tool
- MetaWRAP



MAG Quality Assessment

Use to assess completeness and contamination.

- CheckM
- BUSCO

Taxonomic Assignment

Annotate bins using

- GTDB-Tk

8.6 Downstream Analysis (Optional)

Functional annotation:

- PROKKA / BAKTA
- eggNOG-mapper

Phylogenomics:

- PhyloPhlAn

AMR/virulence:

- ABRicate
- CARD
- VfDB
- AMRfinder Plus

8.7 Best Practices & Notes

- Aim for **>40 Gbp** of Illumina data per sample for deep profiling.
- **ONT long reads** help resolve repeats, rRNA operons, plasmids, and MGEs.
- For low-complexity samples, less sequencing depth may suffice.
- Optional: Use **Hi-C** or read cloud technologies for improved binning.