



Oct 06, 2019

Haematoxylin and eosin staining

 [PLOS One](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.7zg3p3w>

Norfitriah Mohamed Sohaimi¹, Mohd Hair Bejo¹, Abdul Rahman Omar¹, Aini Ideris¹, Nurulfiza Mat Isa²

¹Faculty of Veterinary Medicine, Universiti Putra Malaysia, Serdang, Selangor, Malaysia.;

²Faculty of Biotechnology and Biomolecular Sciences, Universiti Putra Malaysia, Serdang, Selangor, Malaysia.



Norfitriah Mohamed Sohaimi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7zg3p3w>

External link: <https://doi.org/10.1371/journal.pone.0225863>

Protocol Citation: Norfitriah Mohamed Sohaimi, Mohd Hair Bejo, Abdul Rahman Omar, Aini Ideris, Nurulfiza Mat Isa 2019. Haematoxylin and eosin staining. [protocols.io](https://dx.doi.org/10.17504/protocols.io.7zg3p3w) <https://dx.doi.org/10.17504/protocols.io.7zg3p3w>

**Manuscript citation:**

Sohaimi NM, Bejo MH, Omar AR, Ideris A, Isa NM (2019) Molecular characterization of fowl adenovirus isolate of Malaysia attenuated in chicken embryo liver cells and its pathogenicity and immunogenicity in chickens. PLoS ONE 14(12): e0225863. doi: [10.1371/journal.pone.0225863](https://doi.org/10.1371/journal.pone.0225863)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 06, 2019

Last Modified: October 06, 2019

Protocol Integer ID: 28424

Keywords: eosin

Materials

MATERIALS

 Dulbecco's Modified Eagle's Medium (DMEM) Merck MilliporeSigma (Sigma-Aldrich) Catalog #D5796

For detection of intranuclear inclusion bodies, primary CEL cells were grown on sterile coverslips [Adair, 1978]. Briefly, the growth medium was aspirated out from confluent monolayer and washed twice with sterile medium free serum prior to inoculation of 0.1mL viral supernatant. Culture was added with maintenance medium, DMEM and 2% FBS and kept for 48 hours incubation, harvested and fixed in 10% buffered formalin for 5 minutes and stained with haematoxylin and eosin (HE).

