

Dec 11, 2018

## High molecular weight DNA extraction for long read sequencing

DOI

[dx.doi.org/10.17504/protocols.io.uppevmn](https://dx.doi.org/10.17504/protocols.io.uppevmn)

Karen Cristine Goncalves dos Santos<sup>1</sup>, Megan McDonald<sup>2</sup>

<sup>1</sup>Universite du Quebec a Trois-Rivieres; <sup>2</sup>University of Birmingham

High molecular weight DNA extraction from all kingdoms  
Tech. support email: [See@each.protocol](mailto:See@each.protocol)



Karen Cristine Goncalves dos Santos

Université du Québec à Trois-Rivières

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.uppevmn>

**Protocol Citation:** Karen Cristine Goncalves dos Santos, Megan McDonald 2018. High molecular weight DNA extraction for long read sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.uppevmn>

**Manuscript citation:**



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 16, 2018

**Last Modified:** December 11, 2018

**Protocol Integer ID:** 16847

**Keywords:** dna extration from bipolaris sorokiniana spore, bipolaris sorokiniana spore, high molecular weight dna extraction for long read, high molecular weight dna extraction, dna extraction, dna extration, dna, sequencing adaptation

## Abstract

Adaptation of the protocol [dx.doi.org/10.17504/protocols.io.k6qczdww](https://dx.doi.org/10.17504/protocols.io.k6qczdww)

Optimized for DNA extration from Bipolaris sorokiniana spores.

## Guidelines

Optimized for DNA extraction from Bipolaris sorokiniana

**Buffers are best when fresh and not older than 3-6 months. Buffered Phenol:Chloroform:Isoamylalcohol (25:24:1) should not be older than 3 months.**

Critical steps to obtain high quality DNA:

- Do NOT heat samples during DNA extractions! Perform all steps at RT or 4°C as indicated.
- Do NOT incubate samples with KAc for prolonged time periods
- Perform two steps of buffered Phenol:Chloroform:Isoamylalcohol purification to reduce co-purifying metabolites.





## Materials

### MATERIALS

 N-lauryl sarcosine Merck MilliporeSigma (Sigma-Aldrich) Catalog #L5125-50G

### Reagents required

**BUFFER A:** 0.35 M sorbitol

0.1 M Tris-HCl

5 mM EDTA, pH 8

autoclave to sterilize

**BUFFER B:** 0.2 M Tris-HCl

50 mM EDTA, pH 8

2 M NaCl

2% CTAB

autoclave to sterilize

### Other solutions:

Potassium Acetate 5M (KAc precipitate polysaccharides) pH 7.5

Polyvinylpyrrolidone (40000 MW) 10 % [w/v] (Sigma PVP40)

Polyvinylpyrrolidone (10000 MW) 10% [w/v] (Sigma PVP10)

Sodium Acetate (NaAc) 3M pH 5.2

5% Sarkosyl N-lauroylsarcosine sodium salt (SIGMA L5125)

Filter-sterilize

Isopropanol 100%

Ethanol 70%

Buffered Phenol:Chloroform:Isoamylalcohol P:C:I (25:24:1, Sigma P2069)

Autoclave acid washed Sand

### Enzymes

RNAse A or T1 (1000 U/ml, Thermo Fisher EN0541)

Proteinase K (800U/ml, NEB P81072)

## Troubleshooting



## Extraction

1

### Note

This protocol is an adaptation from <https://www.protocols.io/view/high-quality-dna-from-fungi-for-long-read-sequenci-k6qczdww>

2.5 volume of Buffer A 5 mL

2.5 volume of Buffer B 5 mL

1.0 volume of Buffer C 2 mL

PVP 40 10% 1mL

PVP 10 10% 1mL

## Extraction

- 2 Grind tissue/spores with liquid nitrogen in a mortar with sand, use 1g of sand per 100mg of starting material. Grind for 2 mins in 4× 15 sec bursts adding liquid nitrogen after each 15 sec grinding burst.

### Note


Using a sharp scalpal, I scrape the mycelia and the spores from two fully covered petri dishes (~90mm diameter).

- 3 Transfer ground tissue with sand into 50mL Falcon containing lysis buffer. Add 10 uL of RNase A and 100 uL of Proteinase K and mix by inversion until no clumps are visible. Incubate in rotator at 30 rpm for 30 min at RT. ⌚ 00:30:00


- 4 Cool in ice for 5 minutes ⌚ 00:05:00

- 5 Add 2.8 mL (0.2 vol) of KAc 5M, mix by inversion, incubate on ice for max 5 mins ⌚ 00:05:00

- 6 Spin at 4°C and 5000g for 12 mins ⌚ 00:12:00




View of Falcon tube with sand after centrifugation

- 7 Transfer supernatant to fresh Falcon tube containing 17 mL (~1vol) (P:C:I) and mix by inversion.
- 8 Spin at 4°C and 5000g for 12 mins  00:12:00

### Note

With this fungus, there is always an interphase with debris. I take most of the supernatant without taking the phenol phase.

Even if there is still debris at the end of the purification, the AMPure beads will not bind the debris, so to increase your final DNA yield it is better to take as much as the supernatant as possible.



Fungal spore debris visible at the interphase between the phenol (bottom) and aqueous (upper) layers.






9

⇒ go to step #7 undefined


#### Note

Debris is mostly at the interphase, but it also attaches to the tube wall, so it is carried over to the next step.



View of fungal debris after centrifugation of second P:C:I wash

- 10 Transfer supernatant to a new 50 mL falcon and add 10  $\mu$ L of RNase A. Incubate in rotator at 30 rpm for 30 min at RT.
- 11 Add 100  $\mu$ L of Proteinase K and incubate in rotator at 30 rpm for 30 min at RT.

 00:30:00


- 12 Add 1 vol (17mL) of P:C:I, mix by inversion and spin at 4°C and 5000g for 12 mins

🕒 00:12:00


Solution will look milky like this.




After centrifugation, there should be little to no debris in the aqueous layer, as in here.






## DNA precipitation

- 13 Transfer supernatant to a new 50 mL falcon and add 0.1 vol (~1.5 mL) of NaAc and mix by inversion.
- 14 Add 1 vol (~ 17 mL) of cold isopropanol and mix well by inversion.
- 15 Spin at 4°C and 8000g for 30 mins.  00:30:00





Depending on the fungus, the pellet will appear brown, yellow.

- 16 Pour off supernatant, keeping around 1 mL. Use it to dislodge the pellet and transfer it to a 1.5 mL tube.
- 17 Spin in a microcentrifuge at 15 000g for 3 min.  00:03:00
- 18 Discard supernatant and add 1 mL of cold 70% ethanol, flick the tube to help wash the pellet and spin in a microcentrifuge at 15 000g for 3 min.  00:30:00  
 [go to step #18](#) undefined
- 19 Discard the supernatant, spin again to get the remaining ethanol out. Air dry the pellet for 1 min.

## Re-dissolve DNA

- 20 Add 200  $\mu$ L of 10mM Tris pH 8.5 and flick the tube to dislodge the pellet. Let it dissolve over night.



Measure concentration with BR Qubit and absorbance with Nanodrop.

#### Expected result

	Sample*	Qubit (ng/ μl)	Nanodrop (ng/μl)	A260/A2 30	A260/A28 0	Total amount (μg)
	WAI2411	72.9	1923	2.022	2.161	14.58
	WAI3382	89.4	2353	2.006	2.129	17.88

\**Bipolaris sorokiniana* isolates

## AMPure bead cleanup to remove small DNA Fragments

- 21 Vortex AMPure beads and add 0.6 vol of beads to the DNA sample. Mix beads and DNA by flicking the tube.

DO NOT VORTEX BEAD/DNA MIXTURE

Incubate in a rotator 30 rpm for 10min  00:10:00

- 22 Spin the beads/DNA mix for 1-2 sec in a microcentrifuge or table top centrifuge. Place the tube on a magnetic rack. Wait for the beads to pellet against the magnet. Slowly pipet off the supernatant.

#### Note

DNA should now be bound to the AMPure beads, however I often save this supernatant, just in case something goes wrong. If your final elution contains little or no DNA, you can come back to this supernatant and "rescue" the DNA from there by adding an equal volume of beads. If all goes well this supernatant is waste that can be thrown away at the end of the extraction

- 23 Wash beads with fresh cold 70% ethanol, using enough volume to fill the tube. Pipette the ethanol in the tube on the wall opposite of the bead pellet.



24 Pipette off ethanol and repeat wash.

➡ go to step #23 undefined

25 Remove from magnet and spin the beads for 1-2 sec in a microcentrifuge or table top centrifuge.  
Place back on the magnet and let the bead pellet stick. (Should be very quick 1-2 secs).  
Pipette off the remaining ethanol. Do NOT let the ethanol dry.

26 Take the tube from the rack and add 49  $\mu$ L of nuclease-free water and flick the tube to resuspend the beads. Incubate at 37°C for 5 min. ⌚ 00:05:00

27 Spin the tube and place it at the magnetic rack. Wait for the beads to pellet and recover supernatant.  
Measure concentration with BR Qubit and absorbance with Nanodrop.

#### Note

At this step, the ratio A260/A230 falls abruptly, but sequencing occurs without any issue. Repeat clean up if you loose more than 30% of your DNA.  
From here we follow with Oxford Nanopore protocol "1D genomic DNA by ligation (SQK-LSK108)", version GDE\_9002\_v108\_revU\_10Oct2016. For the elution, at each step in the library preparation, we incubate the beads with water/elution buffer at 37°C (5 minutes instead of 2 at the DNA repair and the End-prep steps).