



Sep 07, 2022

Integra Magbead DNA and RNA Extraction for isolated colonies

DOI

dx.doi.org/10.17504/protocols.io.dm6gpjeqjgzp/v1



Sophana Chea¹, Sreyngim Lay¹, Mengheng Oum¹, Gechlang Tang¹, Cheata Hou¹, Manu Vanaerschot², Christina Yek³, Cristina Tato², Jessica Manning³, Vida Ahyong²

¹International Center of Excellence in Research, National Institute of Health, Cambodia;

²Chan Zuckerberg Biohub; ³National Institute of Health



Sophana Chea

International Center of Excellence in Research, National Ins...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpjeqjgzp/v1>

Protocol Citation: Sophana Chea, Sreyngim Lay, Mengheng Oum, Gechlang Tang, Cheata Hou, Manu Vanaerschot, Christina Yek, Cristina Tato, Jessica Manning, Vida Ahyong 2022. Integra Magbead DNA and RNA Extraction for isolated colonies .

protocols.io <https://dx.doi.org/10.17504/protocols.io.dm6gpjeqjgzp/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 02, 2022

Last Modified: September 07, 2022

Protocol Integer ID: 69485

Keywords: Integra, DNA, RNA, Colony, isolated, Extraction, integra magbead dna, rna extraction for isolated colony, rna extraction, rna from isolated colony, isolated colony, quality dna, rna, extraction, generation sequencing, ng, dna

Abstract

This protocol is the process to extract DNA and RNA from isolated colonies. The extracted high-quality DNA or RNA are suitable for Next-Generation Sequencing (NGS).

Guidelines

Adapted from the ZymoBIOMICS MagBead DNA/RNA Kit Manual (Zymo Research, Cat#R2135).

Materials

1. RNase away spray for RNase decontaminants.

✕ RNase AWAY™ Surface Decontaminant **Thermo Fisher Scientific Catalog #7002PK**

2. ✕ ZymoBIOMIC MagBead DNA/RNA **Zymo Research Catalog #R2135**

3. ✕ 100% Molecular grade ethanol

4. ✕ Molecular Grade Isopropanol

5. ✕ Proteinase K w/ Storage buffer 20mg set **Zymo Research Catalog #D3001-2-20**

6. ✕ DNase I Set **Zymo Research Catalog #E1010**

7. ✕ Nuclease-free water **Ambion Catalog #AM9932**

8. 1ml deep well sterile plate.

9. 2ml deep well sterile plate.

10. Hard-shell PCR Plates 96 V-well (Bio-Rad, Cat# HSP9601).

11. PCR Plate Seal, foil (Bio-Rad, Cat# MSF1001).

12. 96S Super Magnet. (ALPAQUA, Cat# A001322)

Equipment

VIAFLO

NAME

96 channel pipette

TYPE

Integra

BRAND

VIAFLO 96



SKU

<https://www.integra-biosciences.com/united-kingdom/en/electronic-pipettes/viaflo-96-viaflo-384#tech-info>

LIN
K

Troubleshooting

Safety warnings

 All steps should be performed at  Room temperature .

Perform the extraction in the extraction room separate from the PCR room.

Respect the Laboratory safety guideline for all steps of the protocol.

Wearing PPE is recommended.






Note

** When reusing tips, make sure to include a bit of extra air aspiration to avoid drops at the bottom of tips when aspirating volumes, and also a bit of extra air blows out at the end of dispensing steps in plates.






Buffer Preparation

30m

- 1
 1. Add  20 mL isopropanol to the MagBead DNA/RNA Wash 1 concentrate.
 2. Add  30 mL isopropanol to the MagBead DNA/RNA Wash 2 concentrate.
 3. Reconstitute lyophilized Proteinase K at  20 mg/mL with Proteinase K Storage Buffer and mix by vortexing. Use immediately or store at  -20 °C .
 4. Reconstitute each vial of lyophilized DNase I with  2.25 mL DNase/RNase-Free water in a conical tube.



30m

Note

For each sample to be treated, prepare DNase I Reaction Mix (scale up proportionally):
Add  45 µL DNase I (reconstituted) and  5 µL DNA Digestion Buffer in a nuclease-free tube.
mix by gentle inversion and place  On ice until ready to use.

Make buffer plates prior to starting protocol



1h

- 2
 1. Pre-make Lysis Buffer plate with  520 µL DNA/RNA Lysis buffer in 1ml deep well plate.
 2. Pre-make Beads plate with  35 µL ZymoBIOMIC MagBinding Beads into 96 V-well PCR plate.






1h

Note

For the Beads plate, make it immediately before starting, <1h prior to starting the protocol, to ensure the beads are kept in suspension.









3. Pre-make DNA/RNA Wash 1 plate with  520 µL MagBead DNA/RNA Wash 1 into 1ml deep well plate. Make it two plates.
 4. Pre-make DNA/RNA Wash 2 plate with  520 µL MagBead DNA/RNA Wash 2 into 1ml deep well plate. Make it two plates.



5. Pre-make 100% Ethanol plate with  1100 μL of 100% Ethanol into a 2ml deep well plate. Make it three plates.
6. Pre-make Prep Buffer plate with  520 μL DNA/RNA Prep Buffer into a 1ml deep well plate.
7. Pre-make water plate with  60 μL Nuclease-free water in a 96 V-well PCR plate. Make it two plates.
8. Spin all plates down for  00:01:00 except for the bead plate. Perform a quick pulse spin down of the bead plate, just enough to get all the liquid down. Centrifuge the rest of the plate at 12 000 rpm for  00:01:00 .

Sample preparation and Proteinase K





31m

- 3
 1. Create a plate map so you know which sample you are adding to each well. Add  50 μL of isolated colonies samples to plate 1 (leave column 12 for water control).
 2. Top up the 1x DNA/RNA Shield to get  750 μL .
 3. Manually add  120 μL of Proteinase K into the 0.2ml 8-strip well.
 4. Use multichannel pipet to add  10 μL of Proteinase K into each sample and mix (plate 1).
 5. Load a set of Integra tips (tip set 1) onto the Integra.
 6. **Program: Pipet/Mix 250ul, 15 cycles, speed 4.** Program the Integra to pipet  250 μL of your samples up and down for  00:01:00 (15 cycles), then incubate at  Room temperature for  00:30:00 . Keep tips.

31m


Sample binding and washing


35m

- 4
 7. **Program: Pipet 250ul.** Add  500 μL total of Lysis Buffer to the sample plate (plate 1).
 8. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix samples and buffer for  00:02:00 . Keep tips.
 9. Aliquot  35 μL of MagBinding Beads into 96 V-well PCR plate.
 10. **Program: Pipet/Mix 20ul, 10 cycles, 2 times, speed 4.** Program the Integra to mix the MagBinding Beads plate, so the beads are fully resuspended.
 11. **Program: Pipet 30ul.** Add  30 μL of MagBinding Beads into the sample plate (plate 1).

35m










12. **Program: Pipet/Mix 250ul, 30 cycles, speed 3.** Program the Integra to mix the sample and MagBinding Beads plate, so the beads are fully resuspended. Continue this Integra Program to mix the sample and MagBinding Beads for  00:20:00 .




13. Transfer the plate/tube to the magnetic stand for  00:05:00 until beads (DNA) have pelleted, transfer the cleared supernatant (RNA) into a new 96 V-well plate.

DNA Purification (Beads)

45m

- 5 14. Change new Integra tips.
15. **Program: Pipet 250ul, 2 times, speed 7.** Dispense a total of  500 μ L MagBead DNA/RNA Wash 1 into sample plate and mix well.
16. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix the Wash 1 buffer with the beads. Keep tips.
17. Place the 96-well magnetic stand underneath the sample plate for  00:02:00 until a bead ring forms.
18. **Program: Manual Pipet 250ul, 2 times, speed 3.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.
19. **Program: Pipet 250ul, 2 times, speed 7.** Dispense a total of  500 μ L MagBead DNA/RNA Wash 2 into sample plate and mix well.
20. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix the Wash 2 buffer with the beads. Keep tips.
21. Place the 96-well magnetic stand underneath the sample plate for  00:02:00 until a bead ring forms.
22. **Program: Manual Pipet 250ul, 2 times, speed 3.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.
23. Change new Integra tips.
24. **Program: Pipet 250ul, 2 times, speed 7.** Dispense a total of  500 μ L 100% Ethanol into sample plate and mix well.
25. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix 100% Ethanol with the beads. Keep tips.
26. Place the 96-well magnetic stand underneath the sample plate for  00:02:00 until a bead ring forms.
27. **Program: Manual Pipet 250ul, 2 times, speed 3.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.
28. Repeat step 24.
29. Dry the beads for  00:10:00 on the magnetic stand.
30. Change new Integra tips.











31. **Program: Pipet 30ul, speed 5.** Dispense a total of  30 μL nuclease-free water into the sample plate.
32. **Program: Pipet/Mix 20ul, 30 cycles, speed 7.** Program the Integra to mix nuclease-free water with the beads. Keep tips.
33. **Program: Manual Pipet 30ul, speed 3.** Transfer the plate to the magnetic stand and pellet the beads for  00:05:00 , then aspirate and dispense the eluted DNA to a new 96 V-well plate.
34. Store DNA sample immediately at  -80 $^{\circ}\text{C}$.

RNA Purification (Supernatant)

45m


- 6 35. Change the new Integra tip.

45m

36. **Program: Pipet 230ul, 3 times, speed 7.** Dispense a total of  690 μL 100% Ethanol to the supernatant.
37. **Program: Pipet/Mix 250ul, 30 cycles, speed 7.** Program the Integra to mix 100% Ethanol with the supernatant. Keep tips.
38. Aliquot  35 μL of MagBinding Beads into 96 V-well PCR plate.
39. **Program: Pipet/Mix 20ul, 10 cycles, 2 times, speed 4.** Program the Integra to mix the MagBinding Beads plate, so the beads are fully resuspended.
40. **Program: Pipet 30ul.** Add  30 μL of MagBinding beads into the sample plate.
41. **Program: Pipet/Mix 250ul, 10 cycles, speed 3.** Program the Integra to mix the sample and MagBinding beads plate, so the beads are fully resuspended. Continue this Integra Program to mix the sample and MagBinding Beads for  00:10:00 .
42. Transfer the plate to the magnetic stand for  00:05:00 until beads have pelleted, then discard the cleared supernatant.
43. **Program: Pipet 250ul, 2 times, speed 7.** Dispense a total of  500 μL MagBead DNA/RNA Wash 1 into sample plate.
44. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix the Wash 1 buffer with the beads. Keep tips.
45. Place the 96-well magnetic stand underneath the sample plate for  00:02:00 until a bead ring forms.
46. **Program: Manual Pipet 250ul, 2 times, speed 3.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.
47. **Program: Pipet 250ul, 2 times, speed 7.** Dispense a total of  500 μL MagBead DNA/RNA Wash 2 into sample plate.




48. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix the Wash 2 buffer with the beads. Keep tips.

49. Place the 96-well magnetic stand underneath the sample plate for  00:02:00 until a bead ring forms.

50. **Program: Manual Pipet 250ul, 2 times, speed 3.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.



51. **Program: Pipet 250ul, 2 times, speed 7.** Dispense a total of  500 μ L 100% Ethanol into the sample plate.

52. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix 100% Ethanol with the beads. Keep tips.

53. Place the 96-well magnetic stand underneath the sample plate for  00:02:00 until a bead ring forms.


54. **Program: Manual Pipet 250ul, 2 times, speed 3.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.

55. Repeat step 51.

56. **DNase I** treatment, use multiple channel pipet to transfer  50 μ L of DNase I Reaction Mix and mix gently for  00:10:00 .


57. **Program: Pipet 250ul, 2 times, speed 7.** Dispense a total of  500 μ L DNA/RNA Prep Buffer into sample plate.

58. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix the DNA/RNA Prep Buffer with the beads. Keep tips.


59. Place the 96-well magnetic stand underneath the sample plate for  00:02:00 until a bead ring forms.

60. **Program: Manual Pipet 250ul, 2 times, speed 3.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.

61. Repeat step 57 to 60.

62. **Program: Pipet 30ul, speed 5.** Dispense a total of  30 μ L nuclease-free water into the sample plate.

63. **Program: Pipet/Mix 20ul, 30 cycles, speed 7.** Program the Integra to mix nuclease-free water with the beads. Keep tips.

64. **Program: Manual Pipet 30ul, speed 3.** Transfer the plate to the magnetic stand and pellet the beads for  00:05:00 , then aspirate and dispense the eluted RNA to a new 96 V-well plate.

65. Store RNA sample immediately at  -80 $^{\circ}$ C .