

Nov 28, 2022

iPSC Transduction

DOI

dx.doi.org/10.17504/protocols.io.q26g7y1y1gwz/v1

Ali Albalakhi^{1,2}, Ning Xia^{1,2}

¹Massachusetts General Hospital; ²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Jackson Schumacher

Massachusetts General Hospital

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7y1y1gwz/v1>

Protocol Citation: Ali Albalakhi, Ning Xia 2022. iPSC Transduction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7y1y1gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 28, 2022

Last Modified: December 12, 2024

Protocol Integer ID: 73279

Keywords: ASAPCRN, ipsc transduction protocol, protocol, transduction

Abstract

This is the iPSC transduction protocol.

Attachments



Transduction Protoco...

38KB

Guidelines

Note:
Cells will need to be maintained in StemFlex medium to reform healthy iPSCs colony-like morphology (4-5 days)
Given that we don't know what viral titer will be toxic to the cells during transduction we will have to optimize the transduction as we progress

Materials

Material needed:

1. Packaged virus
2. Polybrene (diluted 100 μ g/mL)
3. Stem flex (+) warmed to room temperature ~ 10 to 15min
4. Vitronectin-coated plates (coating to final concentration of 1.0 μ g/cm²)

Troubleshooting



Begin transduction with cell line 7026

1

Day 1:

2 Split iPS cells into single cells in complete medium with RVC. (500µl per well)

3 Seed cells at 5.0×10^4 cells in glass-bottom 24-well plates

Day 2: Calculations are for 24-well plates. (If using a 6 well plate, then seed at 1×10^5 cells per well, 850µL of medium the next day with 150µL)

4 Virus infection:

5 Remove the culture medium from cells, add (**300µL**) fresh medium without RVC into each well

6 Add **45µl** virus and 1µg/ml polybrene into cells and mix it. Incubate in 37 degree, 5% CO₂.

	Well 1	Well 2	Well 3	Control
	PLVX mCherry HLA-A2	PLVX IRES- Puro HLA-A2	pLoc MIana	

7 Add Polybrene into the infected well to make the final concentration to be 1ug/ml.

7.1 Stock Polybrene (10mg/mL)

7.2 Make a 100x dilution to make a working concentration of 100µg/mL

$$(10\text{mg/mL}) (V_{\text{stock}}) = (100\mu\text{g/mL})(0.5\text{mL})$$

$$V_{\text{stock}} = 5\mu\text{L of stock into } 495\mu\text{L of DPBS } (-)(-)$$

$$(100\mu\text{g/mL}) (V_{\text{stock}}) = (1\mu\text{g/mL})(0.3\text{mL}) \rightarrow V_{\text{stock}} = \mathbf{3\mu\text{L}}$$
 of working concentration per 300µL of medium added

Day 3:

8 Monitor the cells for the next 48 hours

Day 4 (or 48 hours later):

9 Check IRES-GFP or HLA-A2's expression at least **48 hours** after the transduction.

10 For antibiotic-based selection, the cells will be maintained to reform colonies and subsequently grown in medium with **1µg/mL** of Puromycin to initiate selection