

Oct 23, 2023

JMN-MSMP Notexin Muscle scRNAseq

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwd8r9lmk/v1>

ccherry¹

¹JMN-MSMP

Cellular Senescence Net...



ccherry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwd8r9lmk/v1>

Protocol Citation: ccherry 2023. JMN-MSMP Notexin Muscle scRNAseq. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwd8r9lmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: October 23, 2023

Last Modified: October 23, 2023

Protocol Integer ID: 89762

Keywords: SenNet JMN-MSMP, msmtp notexin muscle scrnaseq murine notexin muscle injury, cell rna sequencing, single cell rna, cell rna, rna, jmn

Abstract

Murine notexin muscle injury followed by single cell RNA sequencing

Skeletal muscle injury and collection

- 1 In a given mouse, both tibialis anterior (TA) muscles were injected with 10 μ l of notexin (10 μ g/ml; Latoxan, France).
- 2 At various days post injury, the mouse was sacrificed and the TAs were collected. Each TA was processed independently to generate single cell suspensions.

Generate single-cell suspensions

- 3 Muscles were enzymatically digested with 8 mg/ml Collagenase D (Roche, Basel, Switzerland) and 10 U/ml Dispase II (Roche, Basel, Switzerland) and then manually dissociated to generate cell suspensions.
- 4 Myofiber debris was removed by filtering the cell suspensions through a 100 μ m and then a 40 μ m filter (Corning Cellgro # 431752 and # 431750).
- 5 After filtration, erythrocytes were removed by incubating the cell suspension in erythrocyte lysis buffer (IBI Scientific # 89135-030).
- 6 After digestion, the single-cell suspensions were washed and resuspended in 0.04% BSA in PBS at a concentration of 10^6 cells/ml. A hemocytometer was used to manually count the cells to determine the concentration of the suspension.

Single-cell RNA library preparation and sequencing

- 7 Single-cell RNA-sequencing libraries were prepared using the Chromium Single Cell 3' reagent kit v3 (10x Genomics, Pleasanton, CA) following the manufacturer's protocol.
- 8 Cells were diluted into the Chromium Single Cell A Chip to yield a recovery of 6,000 single-cell transcriptomes with <5% doublet rate.
- 9 Libraries were sequenced on the NextSeq 500 (Illumina, San Diego, CA).