

May 29, 2025

## Lentiviral Production

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr8ydolmk/v1>

Camille Goldman<sup>1</sup>

<sup>1</sup>Icahn School of Medicine at Mount Sinai



Camille Goldman

Icahn School of Medicine at Mount Sinai, University of Oxf...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.6qpvr8ydolmk/v1>

**Protocol Citation:** Camille Goldman 2025. Lentiviral Production. protocols.io  
<https://dx.doi.org/10.17504/protocols.io.6qpvr8ydolmk/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** August 29, 2024

**Last Modified:** May 29, 2025

**Protocol Integer ID:** 106634

**Keywords:** ASAPCRN

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's**

Grant ID: ASAP-024297

**NASA**

Grant ID: 80ARC022CA004

**NIH/NINDS**

Grant ID: R01NS114239

**NIH/NINDS**

Grant ID: UH3NS115064

**NIH/NIA**

Grant ID: T32AG04968

**NIH/NINDS**

Grant ID: F31NS13090

**New York State Department of Health**

Grant ID: NYSTEM-C32561GG

## Abstract

Protocol for producing lentivirus



## Day -1: Preparing HEK Cells

- 1 Grow HEK cells in Dulbecco's Modified Eagle Medium (DMEM) with 10% bovine calf serum at 37°C in a 5% CO<sub>2</sub> incubator. Use 25 mL for a 15cm plate
- 2 Passage  $\sim 10^7$  HEK cells per 15cm dish coated in 0.1% gelatin

## Day 0: Transfecting HEK Cells

- 3 Change media on HEK cells, just prior to the following steps
- 4 Combine 14.7 ug pMDLg/pRRE (Addgene #12251), 5.7 ug pRSV.Rev (Addgene #12253), and 7.9 ug of pMD2.G (Addgene #12259) per 15 cm plate
- 5 Add 22.5 ug of transfer vector per 15 cm plate
- 6 Add 1 mL of 278 uM CaCl<sub>2</sub> and mix thoroughly
- 7 Add 1mL of 2x BBS solution (280 mM NaCl, 50mM BES, 1.5mM Na<sub>2</sub>HPO<sub>4</sub>; pH = 6.95) dropwise. Invert to mix.
- 8 Incubate for 1 minute at room temperature
- 9 Add dropwise to HEK cells, distributing throughout the plate. Shake the plate on the bench in an X-shape to distribute the plasmid mixture thoroughly

## Day 1-3: Collecting Virus

- 10 Day 1: After 24 hours, change the media
- 11 Day 2: Collect conditioned media after another 24 hours



- 12 Day 3: Collect conditioned media after another 24 hours. Combine with the media collected on day 2

### Day 3-4: Concentrating Virus

- 13 Centrifuge media at 2000xg for 10 minutes to pellet debris
- 14 Filter media through a 0.22um PES filter
- 15 Add PEG8000 to a final concentration of 5%
- 16 Add NaCl to a final concentration of 0.15M
- 17 Invert multiple times to mix well and incubate over night at 4°C
- 18 Centrifuge at 3000xg for 20 min to pellet viral particles
- 19 Discard the supernatant
- 20 Resuspend the pellet in DMEM/F12 (or media of your choice) to 1% of the original supernatant volume
- 21 Store 50uL aliquots at -80°C

