



Feb 06, 2025

Version 2

## Lentivirus preparation for neuronal transdifferentiation V.2

 Nature Cell Biology

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx317wg8j/v2>

Ching-Chieh Chou<sup>1</sup>, Judith Frydman<sup>1</sup>

<sup>1</sup>Department of Biology, Stanford University



**Ching-Chieh Chou**

Stanford University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.kxygx317wg8j/v2>

**External link:** <https://doi.org/10.1038/s41556-025-01623-y>

**Protocol Citation:** Ching-Chieh Chou, Judith Frydman 2025. Lentivirus preparation for neuronal transdifferentiation. [protocols.io https://doi.org/10.17504/protocols.io.kxygx317wg8j/v2](https://doi.org/10.17504/protocols.io.kxygx317wg8j/v2) Version created by [Ching-Chieh Chou](#)

**Manuscript citation:**

Chou, CC., Vest, R., Prado, M.A. *et al.* Proteostasis and lysosomal repair deficits in transdifferentiated neurons of Alzheimer's disease. *Nat Cell Biol* **27**, 619–632 (2025). <https://doi.org/10.1038/s41556-025-01623-y>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 06, 2025

**Last Modified:** February 06, 2025

**Protocol Integer ID:** 119733

**Keywords:** ASAPCRN, Lentiviruse, lentivirus preparation for neuronal transdifferentiation, lentivirus preparation, lentivirus, neuronal transdifferentiation, protocol for the preparation, preparation

**Funders Acknowledgements:**

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

## Abstract

This is a protocol for the preparation of Lentiviruses.

## Preparation of culture medium and plasmids

1. HEK293T medium @4C lasts for 1 month. A total of 500 mL:
  - 435 ml DMEM (High glucose, GlutaMAX) (Thermo Fisher Scientific)
  - 50 ml Fetal bovine serum (FBS) (Thermo Fisher Scientific)
  - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
  - 5 ml HEPES (1M; Thermo Fisher Scientific)
  - 5 ml Sodium pyruvate (100 mM; Thermo Fisher Scientific)
  - Sterilized by filtration through a 0.22  $\mu$ m filter
2. Virus medium @4C lasts for 1 month. A total of 500 mL:
  - 470 ml DMEM (High glucose, GlutaMAX) (ThermoFisher Scientific)
  - 10 ml Fetal bovine serum (FBS) (Thermo Fisher Scientific)
  - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
  - 5 ml MEM NEAA (100X; Thermo Fisher Scientific)
  - 5 ml HEPES (1M; Thermo Fisher Scientific)
  - 5 ml Sodium pyruvate (100 mM; Thermo Fisher Scientific)
  - 0.5 ml beta-mercaptoethanol (Thermo Fisher Scientific)
  - Sterilized by filtration through a 0.22  $\mu$ m filter
3. OptiMEM+GlutaMAX (Thermo Fisher Scientific)
4. Lipofectamine 2000 (Thermo Fisher Scientific)
5. Plasmids:
  - Transfer: FUW-M2rtTA (Addgene plasmid# 20342), pTet-O-Ngn2-puro (Addgene plasmid# 52047), Tet-O-FUW-Ascl1 (Addgene plasmid# 27150), Tet-O-FUW-Brn2 (Addgene plasmid# 27151), Tet-O-FUW-Myt1l (Addgene plasmid# 27152)
  - Packaging: psPAX2 (Addgene plasmid# 12260)
  - Envelope: pMD2.G (Addgene plasmid# 12259)
6. Millex syringe 0.22 & 0.45  $\mu$ m filter

## Culturing HEK293T cells

2. Culture HEK293T cells in a 10-cm dish or T-75 flask and place in a 37°C, 5% CO<sub>2</sub> incubator
3. Monitor cell growth daily and passage when the cells reach 80–90% confluence (~every 2–3 days).

- 4 Before passaging, pre-warm the complete DMEM and 0.05% Trypsin-EDTA at 37°C.
- 5 Aspirate the old medium and rinse the cells once with 5 mL PBS.
- 6 Add 3–5 mL 0.05% Trypsin-EDTA and incubate at 37°C for 3 min until cells detach (tap the flask gently if needed). Add 5–10 mL complete DMEM to neutralize trypsin.
- 7 Transfer the cell-containing solution into a conical tube and centrifuge at  $300 \times g$  for 5 min.
- 8 Aspirate the supernatant and resuspend the pellet in fresh DMEM. Plate cells at the desired density.

## Generation of Lentiviruses

- 9 Seed  $6 \times 10^6$  HEK293T cells in a 10-cm dish for each transfer plasmid and culture to reach 80% confluence next day.
- 10 Prepare DNA and Lipofectamine 2000 dilutions.
  - 5 µg Transfer plasmid
  - 4 µg Packaging plasmid
  - 2.5 µg Envelope plasmid
  - DNA diluted into 300 µL OptiMEM
  - 23 µL Lipofectamine 2000, diluted into 300 µL OptiMEM
- 11 Vortex DNA and Lipofectamine dilutions, centrifuge for 5 sec and incubate at room temperature for 5 mins.
- 12 Combine DNA and Lipofectamine dilutions, vortex, centrifuge and incubate the mixture for additional 15 mins.
- 13 Remove old HEK293T medium, rinse the cells with OptiMEM once and then add 5 mL OptiMEM to each dish.
- 14 Add a total of 600 µL of Lipofectamine/DNA mixture to the cells dropwise across the whole dish, gently shake in vertical  $\times$  horizontal direction a few times and return the dishes to the incubator.



- 15 After 6-hr incubation, remove the transfection medium and add 9 mL of fresh Virus medium to each dish.
- 16 For additional 24-hr incubation, collect the viral particle-containing medium from each dish and store at 4°C.
- 17 After medium collection, add fresh 8 mL of pre-warm Virus medium to each dish.
- 18 After 24-hr incubation, collect the medium and pool it with the first collection.
- 19 The cells can now be bleached and disposed of accordingly.
- 20 Seal the pooled medium with parafilm and centrifuge in the swinging bucket at 400xG for 5 mins at 4°C to pellet cell debris.
- 21 Carefully take the cleared supernatant and pass through a 0.45 µm filter.
- 22 Viral supernatant is stored at 4°C until use and good for 2 weeks.