

Dec 31, 2019

Version 2

## Lentivirus Production V.2

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.bawxiffn](https://dx.doi.org/10.17504/protocols.io.bawxiffn)

Addgene The Nonprofit Plasmid Repository<sup>1</sup>

<sup>1</sup>Addgene



Addgene The Nonprofit Plasmid Repository

Addgene



### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bawxiffn>

External link: <https://www.addgene.org/protocols/lentivirus-production/>

**Protocol Citation:** Addgene The Nonprofit Plasmid Repository 2019. Lentivirus Production. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bawxiffn>

**Manuscript citation:**



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 31, 2019

**Last Modified:** December 31, 2019


**Protocol Integer ID:** 31415

**Keywords:** Lentivirus, Production, Cell Culture, addgene, 293T, packaging cells, Transfection, virus, lentiviral vector, polyethylenimine, PEI, lentivirus production this protocol, lentivirus production, lentivirus, addgene protocol page, protocol

## Abstract

This protocol is for Lentivirus production. To see the full abstract and additional resources, visit the [Addgene protocol page](#).

## Sample Data



**Figure 1:** Lenti-X 293T cells were transfected with the GFP-expression plasmid **pRosetta** using  $\mu\text{g}$  total DNA to  $\mu\text{g}$  PEI ratios of 1:1, 1:2, 1:3 and 1:6. The 1:2 and 1:3 total DNA:PEI  $\mu\text{g}$  ratios provided high transfection efficiencies as measured by the highest proportion of GFP positive cells without limiting cell growth. Left panels: bright field images; right panels: GFP channel images.



## Guidelines

### Workflow Timeline

**Day 0:** Seed 293T packaging cells

**Day 1 (pm):** Transfect packaging cells

**Day 2 (am):** 18 hours post transfection. Remove media, replace with fresh media

**Day 3-4 (am):** Harvest virus



## Materials

### Reagents

- DMEM high glucose
- L-alanyl-L-glutamine (or alternative stable glutamine)
- Heat-inactivated FBS
- Low serum medium such as Opti-MEM or Opti-Pro SFM
- Chloroquine diphosphate, 25  $\mu$ M
- PEI, 1 mg/mL
- Microcentrifuge tubes
- 10 cm tissue culture dishes
- Pipettes
- Pipette tips
- Hydrochloric acid
- Sodium hydroxide
- 0.22  $\mu$ m polyethersulfone (PES) filter
- 0.45  $\mu$ m PES filter
- Syringes for filtering

### Equipment

- Biosafety cabinet
- Pipetman
- Pipettors
- Incubator
- pH meter
- Stir plate
- Magnetic Stir Bar

### Reagent Preparation

1. DMEM Complete: 10% v/v FBS and 4 mM L-alanyl-L-glutamine

- To a 500 mL bottle of DMEM high glucose, add 55 mL of heat inactivated FBS and 11 mL of 200 mM L-alanyl-L-glutamine. Store at 4 °C .



## Note

### \*Pro-Tips\*

- Different brands and lots of FBS can promote or inhibit transfection.
- Test a variety of brands and lots of FBS to find one suitable with your protocols. FBS can be purchased already head inactivated or it can be inactivated in the lab by heating to 56 °C for 00:30:00

## 2. 25 mM chloroquine diphosphate

- Dissolve 0.129 g of chloroquine diphosphate salt into 10 mL of sterile water.
- Filter sterilize through a 0.22 um filter.
- Aliquot 50 µL - 100 µL and store at -20 °C .
- Aliquots can be thawed and stored at 4 °C prior to use. Thawed aliquots should be discarded after 1-2 months.

## 3. 1 mg/mL polyethylenimine, linear MW 25,000 Da (PEI)

- Dissolve 100 mg of powder into 100 mL of deionized water.
- While stirring, slowly add hydrochloric acid until the solution clears.
- Check the pH of the solution
- Use hydrochloric acid or sodium hydroxide to adjust the pH to 7.0. Typically the solution will be basic and will need adjustment with hydrochloric acid first.

## Note

### \*Pro-Tip\*

The pH of this solution will drift pretty rapidly upon addition of acid or base. Add only a few drops at a time, allow them to mix and recheck the pH to prevent over or undershooting the desired pH.

- Allow the solution to mix for 00:10:00 and then recheck the pH to ensure that it has not drifted.
- Filter the solution through a 0.22 um membrane.
- Aliquot 500 µL - 1000 µL into sterile tubes.
- Store the tubes at -80 °C .
- After thawing the solution can be stored at 4 °C for up to 2 months. After 2 months, discard the tube and thaw a new working stock.
- The optimal mass DNA:mass PEI ratio will need to be empirically determined for each new batch of 1 mg/mL PEI and for each cell line.

## Troubleshooting






## Safety warnings

! See SDS (Safety Data Sheet) for safety warnings and hazards.



## Before start

### Considerations Before You Start




- The health of the packaging cell line is critical for obtaining high viral titer.
- 293T cells should be split 3 times a week:
  - Monday: Plate  $1 \times 10^6$  cells in a T75 flask in  15 mL DMEM complete.
  - Wednesday: Plate  $1 \times 10^6$  cells in a T75 flask in  15 mL DMEM complete.
  - Friday: Plate  $8 \times 10^5$  cells in a T75 flask in  15 mL DMEM complete.
- Do not add pen-strep to the media.
- Use cells that are below passage 15 for viral production.






## Seeding cells

- 1 Seed 293T packaging cells at  $3.8 \times 10^6$  cells per plate in DMEM complete in 10 cm tissue culture plates.
- 2 Incubate the cells at  37 °C , 5% CO<sub>2</sub> for ~  20:00:00 .

## Transfection

- 3 Gently aspirate media, add  10 mL fresh DMEM complete containing  25 micromolar (μM) chloroquine diphosphate and incubate ~  05:00:00 .

### Note

For  10 mL of DMEM complete, add  10 μL of  25 millimolar (mM) chloroquine diphosphate.

- 4 Prepare a mixture of the 3 transfection plasmids:

| Reagent                      | Amount per 10 cm dish * |
|------------------------------|-------------------------|
| psP AX2                      | 1.3 pmol                |
| pMD 2.G                      | 0.72 pmol               |
| Transfer Plasmid*            | 1.64 pmol               |
| Opti Pro SFM to total volume | 500 μL                  |






\*Plasmid concentrations and ratios should be optimized for each transfer plasmid.

**psPAX2, pMD2.G**

#### Note

##### \*Pro-Tip\*

Endotoxins can inhibit transfection, therefore, plasmid DNA purification should include an endotoxin removal step. For high quality plasmid DNA, the plasmid should also be propagated in an endonuclease negative *E. coli* strain such as NEB stable.

- 5 Dilute the above  500  $\mu\text{L}$  mixture into  500  $\mu\text{L}$  PEI-OptiPro SFM with enough PEI such that the ratio of  $\mu\text{g}$  DNA: $\mu\text{g}$  PEI is 1:3 (  1000  $\mu\text{L}$  total per 10 cm dish).

Using transfer plasmid **pHAGE TRE dCas9-KRAB** (total  $\mu\text{g}$  of plasmid DNA - 27.8  $\mu\text{g}$ ), this would be  83.4  $\mu\text{L}$  of 1 mg/mL PEI in  416.6  $\mu\text{L}$  of OptiPro SFM per 10 cm dish.

#### Note

##### \*Pro-Tip\*




There can be batch to batch variation when making the PEI working stock, therefore the ratio of  $\mu\text{g}$  DNA: $\mu\text{g}$  PEI needs to be empirically determined. Once a batch of PEI is prepared, transfect cells with a fluorescent plasmid using a variety of ratios. Check the cells 1-2 days after transfection to determine what ratio gives the highest percentage of GFP positive cells.

## 5.1



|  | Ratio of DNA:PEI | Amount of DNA ( $\mu\text{g}$ ) | Volume of 1 mg/mL PEI ( $\mu\text{L}$ ) |
|--|------------------|---------------------------------|---|
|  | 1:1              | 18.9                            | 18.9                                    |
|  | 1:2              | 18.9                            | 37.8                                    |
|  | 1:3              | 18.9                            | 56.7                                    |
|  | 1:4              | 18.9                            | 75.6                                    |
|  | 1:5              | 18.9                            | 94.5                                    |
|  | 1:6              | 18.9                            | 113.4                                   |



Refer to this table for a possible range of ratios to test

- 6 Gently add the diluted PEI to the diluted DNA. Add the diluted PEI dropwise while gently flicking the diluted DNA tube. Incubate the mixture for 15-30 min at room temperature.
- 6.1 Gently add the diluted PEI to the diluted DNA. Add the diluted PEI dropwise while gently flicking the diluted DNA tube.
- 6.2 Incubate the mixture  00:15:00 -  00:20:00 at  Room temperature .
- 7 Carefully transfer the transfection mix to the Lenti-X 293T packaging cells.


**Note**

Add the transfection mix dropwise being careful not to dislodge the cells.
- 8 Incubate the cells for  18:00:00 , or until the following morning.
- 9 The following morning, carefully aspirate the media. Replace the media with  15 mL of DMEM complete.
- 10 Incubate the cells.





## Harvest Virus

- 11 Virus can be harvested at 48, 72, and 96 hours post transfection in individual harvests or a combined harvest where all the individual harvests are pooled.

**Note**

If pooling harvests, transfer the harvested media to a polypropylene storage tube and store at  4 °C between harvest.



- 12 Centrifuge the viral supernatant at ~  500 x g for  00:05:00 to pellet any packaging cells that were collected during harvesting.
- 13 Filter supernatant through a 0.45 µm PES filter.
- 14 The viral supernatant can be stored at  4 °C for several hours but should be aliquotted and snap frozen in liquid nitrogen and stored at  -80 °C as soon as possible to avoid loss of titer.