

Sep 29, 2025

LONP1 knockdown in iPSC-derived microglia

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqpb5xvk5/v1>

Maria Jose Perez J.¹

¹Institut Imagine

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqpb5xvk5/v1>

Protocol Citation: Maria Jose Perez J. 2025. LONP1 knockdown in iPSC-derived microglia. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgqpb5xvk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 25, 2025

Last Modified: September 29, 2025

Protocol Integer ID: 228166

Keywords: microglia lonp1 knockdown in ipsc, derived microglia lonp1 knockdown

Funders Acknowledgements:

ASAP

Abstract

LONP1 knockdown in iPSC-derived microglia

- 1 Obtain shLONP1 plasmid (target sequence: CCAGTGTTCGAAGAAGACCAA; TRCN0000046793) from the MISSION[®] shRNA library (Sigma-Aldrich) and clone into the pLKO.1-puro vector backbone.
- 2 Obtain non-targeting control shRNA from VectorBuilder (VB010000-0005mme)
- 3 Produce lentiviral particles by transfecting HEK293T cells with shRNA plasmids together with the packaging plasmid psPAX2 (RRID:Addgene_12260) and the envelope plasmid pMD2.G (RRID:Addgene_12259)
- 4 Use TransIT-X2 (Mirus Bio) for transfection according to the manufacturer's instructions
- 5 Quantify the concentration of p24 capsid protein using Lenti-X GoStix Plus (Takara Bio)
- 6 Normalize viral amounts to equal p24 levels
- 7 Add p24-normalized lentiviral particles directly to iPSC-derived microglial cultures without using transfection reagents
- 8 Analyze cells at 24, 48, and 72 hours after shRNA treatment