

May 30, 2025

Manual Wastewater Concentration and Extraction V1.0

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldndjxv5b/v1>

Allison Steedman¹, Mark Zeller¹

¹Scripps Research



Andersen Lab

Scripps Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldndjxv5b/v1>

Protocol Citation: Allison Steedman, Mark Zeller 2025. Manual Wastewater Concentration and Extraction V1.0. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldndjxv5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 30, 2025



Last Modified: May 30, 2025

Protocol Integer ID: 219235

Keywords: wastewater, extraction, ceres, concentration, magmax, rna, dna, manual, magnetic-bead based, dna from collected wastewater, collected wastewater, manual wastewater concentration, wastewater, throughput version of this protocol, rna, extraction, throughput version, protocol, dna

Funders Acknowledgements:

Bill & Melinda Gates Foundation

Grant ID: INV-057213

Abstract

This protocol is designed to manually extract RNA/DNA from collected wastewater.

For a high-throughput version of this protocol, please refer to this protocol: <https://www.protocols.io/view/high-throughput-wastewater-sars-cov-2-detection-pi-81wgb74xyvpk/v1>

Materials

Wastewater, Ceres Nanotrap Beads: SKU#44202-30, MagMax Lysis Buffer: A42361, MagMAX Viral/Pathogen Kit: A42352



Concentration

- 1 Invert wastewater sample to mix
- 2 Aliquot 10mL of wastewater into a 15mL conical tube
- 3 Add 150ul Ceres Nanotrap Microbiome A Particle beads into 10mLs of wastewater and invert to mix
- 4 Incubate at room temp for 15 minutes on slow shaker (or invert every 5 minutes)
- 5 Place 15mL conical on magnetic stand for 5 minutes or until beads are fully separated
- 6 Pipette off the supernatant without disturbing the bead pellet
- 7 Remove from magnet and add 500ul MagMax Lysis Buffer. Resuspend the pellet and transfer to a 1.5mL tube.
- 8 Incubate at room temp for 10 minutes on a slow shaker
- 9 Place 1.5mL tube on magnetic stand for 2 minutes or until beads are fully separated
- 10 Collect supernatant into new 1.5mL tube and discard bead pellet. Continue right to extraction.

Extraction

- 11 Prepare the following mixture in a 1.5mL conical tube:



	Reagent	Volume (ul)
	Binding Solution	500
	Beads	20
	Proteinase K	10
	Sample/Lysate	400

*Reagents (aside from sample) from MagMAX Viral/Pathogen Kit (A42352)

- 12 Vortex to mix
- 13 Incubate for 10 minutes
- 14 Place 1.5mL conical on magnetic rack for 2 minutes or until beads are fully separated
- 15 Discard supernatant
- 16 Take off magnetic rack and add 1000ul MagMax Wash Buffer and resuspend bead pellet
- 17 Place 1.5mL conical on magnetic rack for 2 minutes or until beads are fully separated
- 18 Discard supernatant
- 19 With the tube still on the magnet, add 1000 μ l of 80% EtOH, incubate for 30 seconds, and remove the EtOH wash.
- 20 Repeat previous step for another 80% EtOH wash step
- 21 Leave the plate on the magnet to air dry for 3 min.



- 22 Remove the tube from the magnet. Add 50ul of MagMax Elution solution and pipette and vortex to resuspend
- 23 Incubate at room temperature for 10 minutes
- 24 Place the tube onto the magnet. Transfer the supernatant to the final tube. The sample is ready for downstream analysis or can be stored at -80° C

Protocol references

<https://www.protocols.io/view/high-throughput-wastewater-sars-cov-2-detection-pi-81wgb74xyvpk/v1>