

Sep 30, 2019

Version 1

Marchantia agrobacterium transformation of sporelings in multi-well plates (plus materials info) V.1



Forked from [Marchantia agrobacterium transformation of sporelings in multi-well plates](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.7tbhnn>

Linda Silvestri¹, Eftychis Frangedakis², Susana Sauret-Gueto¹, Marius Rebmann¹


¹Plant Sciences, University of Cambridge, OpenPlant; ²University of Cambridge

OpenPlant Project



Susana Sauret-Gueto

Plant Sciences, University of Cambridge, OpenPlant



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7tbhnn>

Protocol Citation: Linda Silvestri, Eftychis Frangedakis, Susana Sauret-Gueto, Marius Rebmann 2019. Marchantia agrobacterium transformation of sporelings in multi-well plates (plus materials info). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7tbhnin>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 30, 2019

Last Modified: September 30, 2019

Protocol Integer ID: 28227

Keywords: marchantia agrobacterium transformation, marchantia agrobacterium transformation of sporeling, sporeling transformation protocol, sterilised spore, days with agrobacterium, mediated marchantia, ishizaki et al 2008 agrobacterium, agrobacterium, sporeling, media with the appropriate selective antibiotic, liquid media in multiwell plate, multiwell plate

Abstract

A modification of the Ishizaki et al 2008 *Agrobacterium* mediated Marchantia sporeling transformation protocol is used.

Sterilised spores are grown for 5-7 days in ½ strength Gamborg plates prior to co-cultivation for 2 days with agrobacterium in liquid media in multiwell plates. Sporelings are then spread on media with the appropriate selective antibiotic. In about 5 days, positive transformants start to emerge.

Materials

Samples required:

Fertilised spore heads from Cam2 archegonia stored at -80°C or 4 °C (dried using self-indicating silica gel – Fisher Scientific S/0761/53)

GV2260 Agrobacterium Electro-Competent Cells - included in these cells is the pSoup helper plasmid (Hellens et al, 2000)

Plasmid of interest to insert T-DNA into Marchantia cells

Overview of Antibiotics Required for selection of Agrobacterium transformants and Marchantia transformants

Agrobacterium strains:

GV2260 strain is resistant to Rifampicin 10µg/ml + Carbenicillin 50µg/ml

GV3103 strain is resistant to Rifampicin 10µg/ml + Gentamycin 25µg/ml (currently not using this strain in the lab).

Plasmid of interest: Loop pCsA plasmid that confers Spec resistance to Agrobacterium and Hyg resistance to Marchantia plants

- GV2260 transformed with a Loop pCsA plasmid: will be screened in LB plates with Spec 100µg/ml + Rif 10µg/ml + Carb 50µg/ml

- Marchantia sporelings infected by Agrobacterium with a Loop pCsA plasmid will be screened on Gamborg plates with Hyg 20µg/ml + Cefo 100µg/ml

Overview Chemicals and Media required:

- **Milton Tablets** purchased from Boots (1 tablet in 25ml sterile water = 0.05%) or Sodium dichloroisocyanurate from Sigma

- **Sterile Water**

- **Acetosyringone** (3',5'-Dimethoxy-4'-hydroxyacetophenone – Sigma D134406) - 100mM stock solution - dissolved in DMSO (x1000, added to 1/2 GB media plus supplements)

- Liquid **LB media** (LB pH7)

- **Antibiotics** stocks

- Cefotaxime 100mg/ml stock solution, dissolved in water (1000x, added to sterile water)

- **LB plates** (LB pH7 agar 1.5%) + **antibiotics** (see ex above)

- **Gamborg plates** (½ Gamborg B5 with vitamins pH5.8 agar 1.2%) for plating spores

- **Liquid Gamborg media plus supplements** = ½ Gamborg B5 with vitamins pH5.8 (Duchefa Biochemie Cat. G0210) + 0.1% N-Z amino A (Sigma C7290) + 0.03% L-Glutamine + 2% sucrose

- **Gamborg plates** (½ Gamborg B5 with vitamins pH 5.8 agar 1.2%) + **antibiotics** (see example above)

Consumables required:



40µm cell strainer,
2mm electroporation cuvettes,
falcon tubes,
empty 1L/500ml bottles,
6-well and/or 12-well Cell Culture Plates (Corning Costar 3516 and 3513),
Eppendorfs

1 Day -7**Spore preparation:**

Marchantia spores are grown for 7 days in ½ strength Gamborg agar plates (A in Figure).

2 Day -2**Agrobacterium preparation:**

- Inoculate a single colony of agrobacterium, transformed with the plasmid of interest, into 5mL LB media plus antibiotics, and then incubate at 28°C with shaking at 150 rpm for two days.

Day 1

- After two days, centrifuge the 5mL agrobacterium culture for 10 minutes at 2000g
- Re-suspend in 5ml of ½ strength Gamborg's B5 plus vitamins media with 1% sucrose.
- Add acetosyringone to a final concentration of 100 µM
- Incubate with shaking for 6h at 28°C at 150 rpm.

3 Day 1

Using a sterile scalpel transfer the 5-7 days old Marchantia sporelings in 6 mL of liquid ½ Gamborg media in a 50 mL Falcon tube (B and C in Figure) and mix well.

4 Aliquote 3 mL of liquid ½ Gamborg B5 vitamins plus supplements media in each well of a 6-well plate.

5 Aliquote 1 mL of spores in each well (D in Figure).

6 Add acetosyringone to each well to a final concentration of 100 µM.

7 Place the 6- well plate on a shaker at 120 rpm for 2 days at 21°C under continuous light (E in Figure).

- 8 Using a pipette transfer the sporeling in a 70 or 100 µm cell strainer placed in a 50 mL Falcon tube (F and G in Figure).
- 9 Wash the sporelings with 50 mL of sterile water supplemented with 100 µg/ml cefotaxime to remove excessive agrobacterium.
- 10 Plate on ½ strength Gamborg B5 plus vitamins 1.2% agar plates with cefotaxime and antibiotics (H and I in Figure)

11 **Day 7**

After 5-7 days successful transformants start to be visible on the plate (J in Figure)

12

 agro-trans-v2_2.png