


Sep 24, 2019

Marchantia spores sterilisation

DOI

dx.doi.org/10.17504/protocols.io.4wbgxan



Linda Silvestri¹, Eftychis Frangedakis¹, Marius Rebmann², Susana Sauret-Gueto²

¹University of Cambridge; ²Plant Sciences, University of Cambridge, OpenPlant

OpenPlant Project



Eftychis Frangedakis

University of Cambridge, Plant Sciences

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4wbgxan

Protocol Citation: Linda Silvestri, Eftychis Frangedakis, Marius Rebmann, Susana Sauret-Gueto 2019. Marchantia spores sterilisation. protocols.io <https://dx.doi.org/10.17504/protocols.io.4wbgxan>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 28, 2019

Last Modified: September 24, 2019

Protocol Integer ID: 25251

Keywords: marchantia spore sterilisation milton

Guidelines

Use 2-4 sporangia per transformation planned



- 1 Add one Milton tablet into 25 mL of sterile water to prepare the "sterilization solution".
- 2 Add the sporangia into a 1.5 mL centrifuge tube and add 0.5 mL of sterilization solution into the tube (B in Figure)
- 3 Use sterile metal tweezers to crush the sporangia (C in Figure).
- 4 Place a 40 μ M cell strainer into a 50 mL Falcon tube and pour the 0.5 mL suspension of crushed spore heads onto the filter and then wash the strainer with another 1 mL of sterilization solution (D in Figure).
- 5 Once filtered, aliquot the spore solution in 1.5 mL Eppendorf tubes (1.5 mL per tube) and incubate at room temperature for 10 min (the exposure to the sterilization solution should not exceed 25 min) (E in Figure). 10m
- 6 Centrifuge for 2 min at 15800 xg at room temperature, discard the supernatant and re-suspend the spores in sterilization solution (150 μ l per transformation planned) (F in Figure). 2m
- 7 Spread 150 μ l of re-suspended spores on a plate containing growth media (G in Figure), seal the lid with micropore tape and place under continuous light at 21 °C 150 μ mol/m²/s for 7 days (H in Figure). Incubate with the bottom side up, which will prevent rhizoids from growing into the agar, making the removal of sporelings from the plate easier by simply using a sterile scalpel.

8

