


Oct 31, 2020

MGH COVID-19 Effort Blood Processing Protocol for PBMC, Neutrophils and Plasma Isolation

DOI

<https://dx.doi.org/10.17504/protocols.io.bjhnkj5e>



Moshe Sade-Feldman¹, Miguel Reyes², Alexandra-Chloe Villani¹

¹Harvard University; ²Broad Institute

Human Cell Atlas Metho...

Coronavirus Method De...



Anita Broellochs

protocols.io

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bjhnkj5e>

Protocol Citation: Moshe Sade-Feldman, Miguel Reyes, Alexandra-Chloe Villani 2020. MGH COVID-19 Effort Blood Processing Protocol for PBMC, Neutrophils and Plasma Isolation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjhnkj5e>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 10, 2020

Last Modified: February 04, 2021

Protocol Integer ID: 40206

Keywords: COVID-19, blood processing, PBMC, neutrophils, plasma isolation, methods for blood processing, plasma isolation this protocol, pbmc, neutrophil,

Abstract

This protocol outlines the steps and methods for blood processing for PBMC to assist the MGH COVID-19 effort.

Attachments



[Villani_AC_MGH_COVID..](#)

418KB

Guidelines

1. The protocol is designed to have two people processing samples: one person should focus on PBMC isolation while the other person on neutrophil isolation and plasma aliquots.
2. This protocol will require the use of two separate BSCs to allow required "social distancing" between individuals processing samples.
3. This protocol will require the use of two separate centrifuges: one for blood and one for neutrophil processing. Neutrophil processing is optional.
4. The number of patient samples that can be processed in a single batch is limited by centrifuge access. Currently plan to be able to do up to 28 patient samples in a batch.
5. Approximately 8 mLs of blood per patient will be split into two parts: 0.25-0.5 mL of the 8 mL will be used for neutrophil isolation and lysis while the remainder will be used to isolate and cryo-preserve PBMCs and freeze plasma aliquots.



Materials

Materials and equipment required for PBMC Isolation

PPE


- Disposable lab gown
- Face shield
- Face mask
- Double gloves
- Sleeve covers
- Shoe covers

Required Equipment

- Pipet-aid and serological pipets (10mL, 25mL)
- P200 pipet and tips
- Microscope
- Hemocytometer
- Tally counter (hand clicker)
- Calculator
- Centrifuge capable of accomodating at least two 15mL conical tubes
- Alcohol-resistant markers
- Mr. Frosty freezing container containing isopropanol, refrigerated at 4 degrees Celsius

Reagents and materials required per draw

- One 50mL SepMate tube
- One 50mL conical tube
- Two 15mL conical tubes
- One 5mL transfer pipet
- Serological pipets: aspirating, 10mL
- 1.5mL cryopreservation tubes
- Trypan blue at room temperature
- 15mL Ficoll at room temperature (NOTE: Ficoll will remain at room temperature in the hood in aluminium foil)
- 16mL RPMI at room temperature
- 50mL RPMI refrigerated
- ~2.5mL Cryostart (CS10) refrigerated



PBMC Isolation OVERVIEW

Materials and equipment required for optional Neutrophil Isolation and Lysis

- 4 × 15mL conical tubes
- Serological pipets: aspirating, 10mL
- 2 × 1.5mL cryopreservation tubes
- 250u: 10X RBC Lysis Buffer
- EasyStep Direct human Neutrophil Isolation Kit (includes Isolation Cocktail and RapidSpheres)
- Trypan Blue
- 1X PBS
- UltraPure 0.5M EDTA, pH 8
- 200 uL TCL
- 2uL 2-Mercaptoethanol, bME (store in corrosives chemical cabinet)
- Pipet-aid and serological pipets
- P1000 and tips
- P200 and tips
- P20 and tips
- EasySep Magnet for 5mL round-bottom tube
- Microscope
- Hemocytometer
- Tally counter (hand clicker)
- Calculator
- Centrifuge capable of accommodating at least two 15mL conical tubes
- Alcohol-resistant markers

Reagent Ordering Information

- 15 mL SepMate tubes: Stemcell Technologies, 85415; 100 tubes
- 50 mL SepMate tubes: Stemcell Technologies, 85450; 100 tubes

- 15 mL conical tubes: VWR, 21008-918; 50 tubes
- Ficoll: VWR cat no. 95021-205; 6×100mL bottles.
- Cryostor CS10: Stemcell Technologies cat no. 7956; 100mL bottle
- PBS: ThermoFisher cat no. 10010023; 500mL bottle
- RPMI: ThermoFisher cat no. 11835055; 10×500mL bottle
- “Mr. Frosty” freezing container: VWR cat no. 55710-200
- VWR lab marker, alcohol-resistant: VWR cat no. 52877-310
- 1.5 mL Cryovials: VWR cat no. 66008-710
- 5mL polystyrene round-bottom tube: VWR International Catalog #60819-820
- UltraPure 0.5M EDTA, pH8: Life Technologies (Thermo Fisher) Catalog # 15575020
- 10X RBC lysis buffer: Life Technologies (Thermo Fisher) Catalog # 00-4333-54
- EasySep Direct Human Neutrophil Isolation Kit: Stemcell Technologies Catalog # 19666
- 3% Acetic acid with Methylene Blue: Stemcell Technologies Catalog # 07060
- PBS: Life Technologies (Thermo Fisher) Catalog #10010049 (10 bottles)
- Buffer TCL: Qiagen Catalog # 1031576
- 2-Mercaptoethanol, bME: Thermo Fisher Catalog #21985023

Safety warnings

 For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

1. For sample hand-off: upon arriving at the lab entrance (doors are locked), the CRC bringing the sample(s) should call the lab. There will be a box in front of the locked door. The CRC with the blood will put the blood sample in the box and once they see one of the researchers through the window, they will walk away. This will ensure that samples are never left unattended. Once the CRC is gone, the researcher will open the door with proper PPE, and wipe down the box containing the samples with 10% bleach before entering any lab space.
2. Perform all procedures in a BL2+ lab space inside a sterile hood.
3. All plastics (tubes, tips, pipets) should be disposed of in a large beaker filled with 10% bleach.
4. Use a sealed container (e.g.- covered centrifuged buckets, plastic box, or ziploc bag) when transferring samples in and out of the hood. Spray down and wipe the secondary container with 10% bleach followed by 70% EtOH before putting inside the hood. After opening the bag/box, spray down the tube inside with 10% bleach followed by 70% EtOH.
5. When using covered centrifuged buckets, spray with 70% EtOH before putting inside the hood. Upon closing the buckets, spray again the covered buckets with 70% EtOH before moving the buckets outside the BSC (prior to moving the buckets back to the centrifuge).
6. When moving the centrifuge buckets to the centrifuge with samples in them, make sure to use a cart.
7. After finishing your work and/or before a new person starts working in the hood, make sure to clean surfaces first with 10% bleach, wait approximately 10 minutes, then spray with 70% EtOH.

Before start

First thing EVERY morning:

- RPMI at room temperature: Note that before starting work on blood samples, RPMI (20mL per sample) will need to be brought at room temperature. Every morning, bring out at least 2 bottles of RPMI. It may take up to 30 minutes to bring it to room temperature. Also note that only the 2:1 blood dilution step requires **RT RPMI**, the remainder of the protocol uses **cold** RPMI to maintain cell viability.
- Ficoll at room temperature: Ficoll should be used at room temperature. We will keep a ficoll bottom in the BSC in aluminium foil. Every morning, make sure you have at least half a bottle full in the hood. If not, bring a new bottle wrapped with aluminium foil.
- Prepare 1mL TCL + 1% BME in a 1.5mL eppendorf tube that should be enough for the day. This may need to be adjusted based on number of samples received. Make sure to RNaseZap the surface of the chemical hood and pipettes before transferring.
- Prepare two trash buckets with 10% bleach solution (one for each hood).
- Make sure there are enough labels in the label printer.

Receiving batches of blood and transport to the BSC: Blood samples will be processed every 3 hours. Upon receiving a batch of blood samples, in a team of 2 people, bring the bag(s) containing the vacutainer in the BL2+ cell culture space. Spray down and wipe the bag(s) containing the vacutainers with 10% bleach followed by 70% EtOH before putting inside the hood. After opening the bag/box, spray down the tube inside with 10% bleach followed by 70% EtOH.

Sample organizing and recording: In a team of two, one person will organize the vacutainers in the hood in holder, taking care of grouping the 2 vacutainers collected for each individual together. Person #1 will read the information on the vacutainers to person #2 that will take note in a spreadsheet of:

- sample ID (will be already encoded)
- date
- time of blood collection
- blood volume received or at least whether 1 or 2 vacutainers was received
- time at which blood processing started

Labeling tubes: Upon compiling information for EACH patient, person #1 labels the following sets of tubes with patient code for all collected samples:

- one 50mL falcon tube (for RPMI mixing)
- one 50mL SepMate tube (if receive < 3mL of blood, use 15mL SepMate tube)
- one 15mL falcon tube to keep empty (for neutrophil step)

Organize the tubes in the racks for blood transfer.



50ml Rack 2nd row	50ml Sepmate with 16ml Ficol #1	50ml Sepmate with 16ml Ficol #2	50ml Sepmate with 16ml Ficol #3
50ml Rack 1st row	50ml tube Blood+RPMI #1	50ml tube Blood+ RPMI #2	50ml tubeBlood+ RPMI #3

15mL rack 1st row	15ml empty sample #1	15ml empty sample #2	15ml empty sample #3
Rack that holds vacutainer	2X vacutainers sample #1	2X vacutainers sample #2	2X vacutainers sample #3



PBMC Isolation

1 Follow the steps detailed in the "Before Start" section prior to moving forward in this protocol.



2 Invert vacutainer 5x to mix blood.



Note

Do not shake.

3 For transferring blood to the different tubes, proceed **ONLY** one patient sample set at a time and avoid opening and closing vacutainers.

4 Gently open one vacutainer and move whole blood to a separate 15mL conical tube for neutrophil isolation based on the amount of blood received.



- 4-8ml blood - take 5000ul for neutrophil isolation
- 2-3ml blood - take 250ul for neutrophil isolation
- don't take blood for neutrophil isolation if there is < 2ml of blood

Note

Do not put on ice.

Note

This step could be skipped if neutrophil isolation is not needed.

4.1 After completing blood transfer for neutrophil, disinfect the 15mL conical tubes with 10% bleach solution and handoff the tubes to person #2 who can move the tubes to another BSC to proceed with RBC lysis and neutrophil isolation.

5 Using a 10mL serological pipet, transfer the remaining blood in the EMPTY 50mL tube.



Note

Take note of the total volume of blood transfer (you need the number to calculate RPMI volume).

- 6 After finishing the blood transfer, move serological pipet to the bleach bucket.
- 7 For each of the tubes listed in Step 5, make sure to fill with the correct buffer:
 - Fill 50mL conical with (2x volume of blood) RPMI at  Room temperature
 - Fill one 50mL SepMate tube with  16 mL room temperature Ficoll and add through the hole in the middle of the cone inside the tube.
 - If the blood volume is < 3mL, use a 15mL SepMate tube and fill with  4.5 mL Ficoll
- 8 Repeat Steps 2-7 until all blood samples have been transferred to the empty 15mL falcon tube and 50mL falcon tube containing blood.
- 9 Mix all the 50mL tubes by inverting 5x. 
- 10 Using a different 25mL serological pipet for each patient, transfer the diluted blood into the SepMate tube. 
- 11 Slowly pipet down the side of the tube to avoid blood mixing with Ficoll that may come out of the hole in the insert (< 1 mL/s. After finishing the blood transfer, move serological pipet to bleach bucket). 
- 12 Centrifuge at  1200 rpm, 20°C, 00:20:00 with maximum acceleration and the brake on. 
- 13 During centrifugation, prepare the following:
 - Grab Mr. Frosty from  -80 °C freezer if not already in lab refrigerator
 - An ice bucket
 - Pre-fill round-bottom well 96-well plate with  10 µL Trypan Blue and prepare one well per sample (for counting)
 - Label one 15mL conical for PBMC per sample, leave empty
 - Grab RPMI bottle from  4 °C fridge and keep  On ice



- 14 After centrifugation, carefully remove the SepMate tube and verify the PBMC layer is visible.

Note

If no layer is visible, centrifuge again at  1200 rpm, 20°C, 00:10:00 with the brake on.

- 15 Place the SepMate tube into a rack.

- 16 Change the centrifuge temperature to  4 °C and close the lid.

- 17 Using a 10mL serological pipet, collect  5 mL plasma from the tube and transfer into the empty 15mL conical tube.

- 18 Place  On ice until later.

- 19 Remove remaining plasma by aspirating and leave approximately  3 mL above the layer.

- 20 With a transfer pipet, collect and transfer PBMC layer from the SepMate tube and move to the 15mL conical tube.

- 21 Top up the conical tube to 15mL with **cold** RPMI.

- 22 Centrifuge  300 rpm, 4°C, 00:05:00 using max acceleration and max break.



- 23 Remove the supernatant by aspirating, making sure to avoid disturbing the pellet, and place tube  On ice .

- 24 Resuspend the pellet in  1 mL cold RPMI .



25 Place  On ice .

Note

If pellet is not visible, resuspend in  500 μ L or less cold RPMI.

26 Add  10 μ L cell suspension to the round bottom 96-well plate containing  10 μ L Trypan Blue .

27 Count the cells to determine the yield and number of aliquots to be frozen.

Note

See the "Counting PBMCs" protocol.

28 Top up the tube to 10mL with RPMI.

29 Centrifuge the cell suspension at  300 rpm, 4°C, 00:08:00 with max acceleration and max break.

30 During centrifugation:

- gather cell freezing materials: cold Crystor (CS10), refrigerated Mr. Frosty, cryo vials
- print cryotube labels and stick on tubes: label with "Blood PBMC", subject ID, and timepoint "PXX-DY", aliquot number, number of cells, and date
- for each patient, we will freeze aliquots of  0.5 Mass Percent cells each and freeze the extra leftover as a separate aliquot.

Note

If the leftover cell is less than 200K cells, then distribute the 200K cells across all aliquots.

31 Remove supernatant through aspiration.



32 Resuspend pellet in cold Cryostor (0.5mL per 0.5M cell aliquot) using the total volume to be used to cover all aliquots accounted for per sample.

33 Resuspend pellet by gently pipetting up and down 10x.

Note

The Cryostor solution is viscous so be sure to pipette up and down slowly.

34 Pipet  0.5 mL cell suspension intro cryotubes.

35 Transfer to Mr. Frosty and immediately place in  -80 °C freezer for a maximum of

 12:00:00 .

35.1 Record the time the blood sample was frozen.

36 The next day, cells can be moved from Mr. Frosty to a cardboard box (labeled PBMC-1...so on).

Note

Note down on the spreadsheet which box and position the samples are placed.

Note

Frozen cells can be kept between  24:00:00 and  168:00:00 in  -80 °C before being send to the Broad Institute on dry ice.

Counting PBMCs and Determining # of aliquots of PBMCs

37 Make sure the cell suspension is well mixed by pietting up and down.

38 Transfer 10 μ L cell suspension into cryovial containing 10 μ L Trypan Blue and mix with a pipet slowly, up and down 5x.



39 Load 10 μ L diluted cells into one of the hemocytometer reservoirs.

Note

If using an automated cell counter (Countess- Invitrogen or BioRad), make sure to load 10 μ L of the diluted cells into the right plastic cell counter chamber. Insert into the machine and press the "Count Cells" button. Cells will be counted automatically and concentration per ml will be given at the end, including the viability of cells. When using the BioRad counting machine, make sure the lower gate is on 7 μ m to avoid counting red blood cells.

40 Count the number of cells in one corner 4 $\times$ 4 grid (count all cells in the squares).



- PBMC cells should have a light blue hue (due to nucleus staining)
- count the cells in contact with the borders of only two sides of the squares




Figure 1: Hemacytometer grid




Figure 2: Cell counting guidelines

Plasma Protocol

41 Once samples are in the freezer, centrifuge the plasma at 1000 rpm, 4°C, 00:05:00 with max acceleration and max brake.



42 During centrifugation, print cryotube labels and stick on 3 cryotubes: label with "Plasma", subject ID, timepoint "PXX-DY", aliquot number, and date.



- 43 Without disturbing the pellet, aliquot  1.5 mL supernatant into 3 cryovials (4.5mL total).

Note

The remaining can be discarded.

- 44 Freeze the aliquoted plasma at  -80 °C in a normal cardboard box or plastic box (labeled Plasma-1...so on).

- 44.1 Take note on the spreadsheet which box and position the samples are placed.

- 44.2 Take note of the time of when the samples were frozen.

Note

Plasma samples can be shipped to the Broad Institute with the PBMCs.

Neutrophil Isolation (OPTIONAL)

- 45 Before starting the upcoming sections:
- Set centrifuge to room temperature for RBC lysis
 - Have EasySep magnet on hand
 - If the anticoagulant in the vacutainer is not EDTA, add 1uL of 0.5M EDTA to the blood
 - Make PBS+EDTA buffer: 500mL PBS + 1mL of 0.5M EDTA and keep at room temperature
 - Label three round-bottom tubes and one 15,L conical tube with patient ID



RBS Lysis

- 46

Note

This step should be performed while the other person is preparing for Ficoll isolation.

Organize the 15mL conical tubes containing  500 µL whole blood in the hood.

47 Add  8 mL ACK lysis buffer to blood using a 10mL serological pipette and mix the sample 5 times.



48 When done, move the pipette to the bleach bucket.

49 Spin  300 x g, Room temperature, 00:05:00 with max acceleration and max deceleration.



50 Aspirate supernatant.

Note

Be careful not to disturb the pellet.

51 Resuspend cells in  250 µL PBS+EDTA buffer at  Room temperature by pipetting up and down gently 5x.



52 Continue to EasySep isolation.

Neutrophil EasySep Isolation

53

Note

Aft anytime, **DO NOT** cap the 15mL tube.

Vortex RapidSpheres for  00:00:30 .

54 Mix the isolation cocktail with 200uL pipette up and down 10x.



		RapidSphere volume	Isolation Cocktail



			volume
	500-1mL blood	50uL	50uL
	250uL blood	25uL	25uL

55 Add isolation cocktail to sample.



56 Add RapidSpheres to sample.



57 Mix gently by pipetting up and down 10x.



58 Incubate for 00:05:00 at Room temperature .



59 After 00:05:00 , top up to 4 mL with PBS+EDTA buffer at Room temperature .

60 Close the tube tightly and then mix gently by tilting the tube back and forth 5x.



61 Place tube into magnet.

Note

Make sure each 15mL tube sits well at the bottom of it without lid.

62 Incubate at Room temperature for 00:05:00 on magnet.



Note

While on the magnet, label a second set of 15mL tubes.



- 63 After  00:05:00 , using a 5mL serological pipet, transfer supernatant into new 15mL tube (tube 2). 

Note

Angle the pipet such that the tip is opposite the magnet to avoid touching the beads.

Note

Procedure is negative selection so neutrophils will be in supernatant.

- 64 Add  25 μ L RapidSpheres to the new 15mL tube (tube 2) if started with  500 μ L whole blood or  12.5 μ L if started with  250 μ L whole blood.

- 65 Close the tube tightly and mix gently by tilting the tube back and forth 5x (avoid bubbles). 

- 66 Incubate  00:05:00 at  Room temperature . 

- 67 After  00:05:00 , place tube on magnet (without lid).

- 68 Incubate  00:05:00 at  Room temperature on magnet. 

Note

While on magnet, label a third set of 15mL tubes.

- 69 After  00:05:00 , using a 5mL serological pipet, transfer supernatant into new 15mL tube (tube 3). 



70 Place the new tube (tube 3) with the supernatant back on the magnet.

71 Incubate for  00:05:00 at  Room temperature on magnet.



Note

While on magnet, label a 15mL flask tube (tube 4).

72 After  00:05:00, using a 5mL serological pipet, transfer supernatant (~3.5mL) into 15mL tube (tube 4).



73 Spin  300 rpm, Room temperature, 00:05:00 with full acceleration and full deceleration.



74 Aspirate supernatant, being careful not to disturb the pellet.

75 Resuspend in  1 mL PBS+EDTA and mix well.



76 Count the cells to determine the yield and number of aliquots to be frozen.



Note

See "Cell Counting" in the next section.

Cell Counting

77 Add  10 μ L resuspended cells to  10 μ L Trypan Blue .

78 Mix well by pipetting up and down.



79 Add  10 μ L to a hemocytometer.

Note

If using an automated cell counter(Countess- Invitrogen or BioRad), make sure to load 10ul of the diluted cells into the right plastic cell counter chamber. Insert into the machine and press the "Count Cells" button. Cells will be counted automatically and concentration per ml will be given at the end, including the viability of cells. When using the BioRad counting machine, make sure the lower gate is on 4um to avoid counting red blood cells.

80 Count the number of cells with a hemocytometer.



- There are four section of 4×4 grids (outlined in red below)
- Count cells in two of the 4×4 grids (take average between the two grid counts)
- When counting, include cells on the edges of 2 sides of each grid, and exclude cells on the edges of the other 2 sides of each grid (as shown below on right)
- If there are not many cells, count number of cells in all four grids and take the average of all four (will be more accurate)
- Take average of cells counted in 4×4 grids to use as "# cells" for the calculation below.




Figure 1: Hemacytometer grid




Figure 2: Cell counting guidelines

Cell Lysis

81 Desired outcome:

- 2 cryotubes per patient, each tube with 100uL of neutrophil lysate in each
- final concentration of cells in RNA lysis buffer = 1000 cells per ul
- a total of 200K cells to split between two cryotubes

82 Mix 200K cells to a 1.5mL Eppendorf tube by transferring an appropriate volume of the cell suspension.





82.1 Determine fraction of cells to transfer to new tube: $20/(\text{'\# cells'})$.

83 Top up tube of cells to  1.5 mL with PBS+EDTA.

84 Spin  300 rpm, Room temperature, 00:05:00 .

84.1 Make fresh RNA lysis buffer (TCL + 1% bME) according to number of samples received.

Note

Make sure to RNAZap the surface of the hood where the lysis buffer will be used and pipettes to be used and use Rnase free tips for pipetting.

84.2 Print cryotube labels and stick on tubes: label with "Neu Lysate", subject ID and timepoint "PXX-DY", aliquot number, concentration of cells ("1k/uL"), and date.

85 Pipe out supernatant and remove excess with P200 (be cautious).

86 Resuspend in  200 μL Lysis Buffer .

Note

If lower number of cells are received, resuspend cells in lower volumes to reach 1000 cell/1mL.

86.1 If less than 10,000 cells are recovered, resuspend in  10 μL and note cell concentration on the cryotube. Freeze single cryotube.

87 After resuspending cells in  200 μL , mix well and transfer  100 μL to two cryotubes.

88 Place tubes on dry ice immediately.



- 89 After  00:05:00 , the aliquoted lysates can be moved to the  -80 °C in a normal cardboard or plastic box.

Note

Note down on the spreadsheet which box and position the samples are placed and take note of the time when the samples were frozen.