

Jan 11, 2023

MinION Sequencing protocol for Rabies virus of Arctic lineage



Forked from [MinION Sequencing protocol for Rabies virus of Indian subcontinental lineage](https://dx.doi.org/10.17504/protocols.io.3byl4jzb8lo5/v1)

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4jzb8lo5/v1>

pramadaprasad¹, Harsha P.K.¹, Chitra Pattabiraman¹

¹National Institute of Mental Health and Neuro Sciences

NIMHANS



pramadaprasad

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4jzb8lo5/v1>

Protocol Citation: pramadaprasad , Harsha P.K., Chitra Pattabiraman 2023. MinION Sequencing protocol for Rabies virus of Arctic lineage. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4jzb8lo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 23, 2022

Last Modified: January 11, 2023

Protocol Integer ID: 74417

Keywords: MinION, Sequencing, Rabies, whole genomes of rabies virus, protocol for rabies virus, rabies virus, arctic lineage an ont minion, sequencing protocol, sequencing, minion, arctic lineage, whole genome, ont minion, genome, zika, virus

Abstract

An ONT MinION-based sequencing protocol to retrieve whole genomes of Rabies virus of the Arctic lineage. Primers were selected using PrimalScheme (<https://primalscheme.com/>) and the protocol was modified from the Zika and SARS-CoV-2 sequencing protocols using the same tool.

Protocol materials

✕ QIAamp Viral RNA Mini Kit **Qiagen Catalog #52904**

✕ Q5 High-Fidelity DNA Polymerase - 500 units **New England Biolabs Catalog #M0491L**

✕ AMPure XP beads **Beckman Coulter**

✕ Qubit® dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**

✕ NEBNext Ultra II End Prep Reaction Buffer **New England Biolabs Catalog #E7647**

✕ NEBNext Ultra II End Prep Enzyme Mix **New England Biolabs Catalog #E7646**

✕ Native Barcoding Expansion 1-12 (PCR-free) **Catalog #EXP-NBD104**

✕ Native Barcoding Expansion 13-24 (PCR-free) **Oxford Nanopore Technologies Catalog #EXP-NBD114**

✕ NEB Blunt/TA Ligase Master Mix **Catalog #M0367**

✕ Ligation sequencing kit 1D **Oxford Nanopore Technologies Catalog #SQK-LSK109**

✕ Flow Cell Priming Kit **Oxford Nanopore Technologies Catalog #EXP-FLP002**

✕ ONT MinION Flow Cell R9.4.1 **Oxford Nanopore Technologies Catalog #FLO-MIN106D**

✕ LunaScript RT SuperMix Kit - 100 rxns **New England Biolabs Catalog #E3010L**



RNA Extraction

- 1 Extract and purify RNA from Brain tissue homogenate using

 QIAamp Viral RNA Mini Kit **Qiagen Catalog #52904** .

Note

Removal of genomic DNA by on column DNase treatment (optional)
Genomic DNA can be removed by on column

 DNase I (RNase-free) - 5,000 units **New England Biolabs Catalog #M0303L**
treatment for 30 min at room temperature before the wash steps.

cDNA Preparation

20m

- 2

20m

Note

Extracted RNA samples are mixed thoroughly by brief vortex and spin down to collect all liquid at the bottom.

- 2.1 Prepare the master mix for cDNA Prep as follows in a PCR plate/ tubes for 10 μ L reaction.

A	B
Component	Volume
Lunascrit RT Supermix (5X)	2 μ l
Template	8 μ l
Total	10 μ l

 LunaScript RT SuperMix Kit - 100 rxns **New England Biolabs Catalog #E3010L**

**Note**

Prepare all reagents in a biosafety cabinet designated for clean work.

2.2 Incubate the reaction mix following the conditions mentioned below in a thermo cycler:

🌡️ 25 °C for 00:02:00

🌡️ 55 °C for 00:10:00

🌡️ 95 °C for 00:01:00

Hold at 🌡️ 4 °C

Amplicon PCR

5h

- 3 Primers are designed using Primal Scheme. Resuspend the lyophilized oligos fully in Nuclease free water / 1X TE to a concentration of 100µM, vortex thoroughly, Nuclease-free, and spin down.

5h

 RABV_Arctic_primers.csv

- 4 Separate odd and even primer sets. Add 5µl of each odd primer to a 1.5mL LoBind tube labelled Pool 1. Repeat the same for all even primers for Pool 2. These are the 100µM stocks of each primer pool.
- 5 Dilute 100µM stock pools (1:10) in Nuclease free water, to generate 10µM primer stock. A final concentration of 15 nM per primer is used. Here we use a total of 39 primers in Pool 1 and 39 primers in Pool 2. Here 0.75µl of per primer pool (10µM) is required for a 25µl reaction.
- 6 Set up two PCR reactions per sample (for Pool 1 and 2) as follows in tubes or plates. Mix by gentle pipetting and pulse spin to collect liquid at the bottom.



	A	B	C
Component	Reaction Mix 1	Reaction Mix 2	
5X Q5 Reaction Buffer	5µL	5µL	
10 mM dNTPs	0.5µL	0.5µL	
Q5 Hot Start DNA Polymerase	0.25µL	0.25µL	
Pool 1 (10µM)	0.75µL	0µL	
Pool 2 (10µM)	0µL	0.75µL	
Nuclease-free water	18.5µL	18.5µL	



Q5 High-Fidelity DNA Polymerase - 500 units **New England**
Biolabs Catalog #M0491L

- 7 Add 2.5µl cDNA to each of the PCR reactions, mix by gentle pipetting and spin the tube to collect all liquid at the bottom of the tube.
- 8 Set up the following PCR conditions in the thermal cycler.

Stage	0 °C	00:00:00	Cycles
Heat Activation	98 °C	00:00:30	1
Denaturation	98 °C	00:00:15	25-35
Annealing	65 °C	00:05:00	25-35
Hold	4 °C	Indefinite	1



Citation

Quick J, Grubaugh ND, Pullan ST, Claro IM, Smith AD, Gangavarapu K, Oliveira G, Robles-Sikisaka R, Rogers TF, Beutler NA, Burton DR, Lewis-Ximenez LL, de Jesus JG, Giovanetti M, Hill SC, Black A, Bedford T, Carroll MW, Nunes M, Alcantara LC Jr, Sabino EC, Baylis SA, Faria NR, Loose M, Simpson JT, Pybus OG, Andersen KG, Loman NJ (2017)

. Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples.

Nature protocols.

<https://doi.org/10.1038/nprot.2017.066>

LINK

PCR clean-up

- 9 Label the tubes/plates and pool the PCR reactions for each samples into a single LoBind tube .

	A	B
	PCR pool 1	25µL
	PCR pool 2	25µL
	Nuclease Free water	0µL
	Total	50µL

Note

When higher viral loads are anticipated based on Ct values < 20 then reduce the amount of amplicons to 15-20µL each and make up the volume to 50µl with Nuclease free water.

- 10 Clean-up the amplicons using  AMPure XP beads **Beckman Coulter** .



Thaw the AMPure beads to room temperature and vortex thoroughly to ensure that the beads are well suspended/ until the solution is uniformly brown in colour.

- 10.1 Add 0.4x volume (0.4:1) of AMPure beads to the sample tube and mix gently by either flicking or pipetting. For example add 20,µl AMPure beads to a 50 µl reaction.
- 10.2 Spin the tubes to collect entire liquid at the bottom and incubate the tubes at room temperature for 00:05:00.
- 10.3 Place the tubes on magnetic stand and incubate for 00:02:00 or until the beads have settled and the supernatant is completely colourless.
- 10.4 Carefully remove and discard the supernatant, without disturbing the bead pellet.
- 10.5 Add 200µl of 70% ethanol or enough to cover the pellet and wait for 30s.Carefully remove and discard ethanol delicately without disrupting the bead pellet.
- 10.6  and repeat Spin to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible.
- 10.7 Incubate the tube with open lid for 00:00:30 or until the pellet is dried (if the pellet dries completely it will crack and become difficult to resuspend)
- 10.8 Resuspend pellet in 30µl Elution Buffer (EB), mix gently by either flicking or pipetting and incubate for 00:02:00.Make sure that the pellet is completely resuspended.
- 10.9 Place on magnetic stand and transfer sample to a clean 1.5mL LoBind tube ensuring no beads are transferred into this tube.
- 10.10 Quantify 1µl product using the Qubit Fluorometer using the

 Qubit® dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854** .

 Qubit® dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**

End-Prep Reaction

40m



- 11 In a new PCR plate/ tubes set up the following reaction for each sample. Prepare Master mix and aliquote it into tubes . If the Qubit reading is > 15ng/ul ,

A	B
Component	Volume
PCR dilution from previous step	3.3μL
Ultra II End Prep Reaction Buffer	1.2μL
Ultra II End Prep Enzyme Mix	0.5μL
Nuclease-free water	5μL
Total	10μL

If the Qubit reading is 5ng/ul -15ng/ul

A	B
Component	Volume
PCR dilution from previous step	5μL
Ultra II End Prep Reaction Buffer	1.2μL
Ultra II End Prep Enzyme Mix	0.5μL
Nuclease-free water	3.3μL
Total	10μL

If the Qubit reading is <5ng/ul,



A	B
Component	Volume
PCR dilution from previous step	8.3μL
Ultra II End Prep Reaction Buffer	1.2μL
Ultra II End Prep Enzyme Mix	0.5μL
Nuclease-free water	0μL
Total	10μL

🔗 NEBNext Ultra II End Prep Reaction Buffer **New England Biolabs Catalog #E7647**

🔗 NEBNext Ultra II End Prep Enzyme Mix **New England Biolabs Catalog #E7646**

11.1 Reaction conditions for End Prep:

Incubate

🔥 25 °C for 00:15:00

🔥 65 °C for 00:15:00

🧊 On ice for 00:01:00

Native Barcoding

12 In a new labelled tube/plate prepare the following components.

A	B
Component	Volume
End-preparation reaction mixture	0.75μL
NBXX barcode	1.25μL
Blunt/TA Ligase	5μL

	A	B
	Master Mix	
	Nuclease-free water	3μL
	Total	10μL

 Native Barcoding Expansion 1-12 (PCR-free) **Catalog #EXP-NBD104**

 Native Barcoding Expansion 13-24 (PCR-free) **Oxford Nanopore Technologies Catalog #EXP-NBD114**

 NEB Blunt/TA Ligase Master Mix **Catalog #M0367**

12.1 Incubate at

 25 °C for 00:20:00

 65 °C for 00:10:00

 On ice for 00:01:00

12.2 In a new 1.5mL LoBind tube pool all barcoded reactions together.

12.3 Add 0.4x volume of AMPure beads to the sample tube and mix by either flicking or gentle pipetting. For example add 96μl AMPure beads to 240μl pooled barcoding reactions.

12.4 Briefly vortex and mix the reactions and incubate at room temperature for 00.05.00.

12.5 Place on magnetic stand and incubate for 00:02:00 or until the beads have pelleted and the supernatant is completely clear. Carefully remove and discard the supernatant, being careful not to touch the bead pellet.

12.6 Remove the tube from the magnetic stand and add 200μl SFB and resuspend beads completely using pipette. Spin down to collect liquid at the bottom of the tube and place on the magnetic stand. Remove supernatant and discard.

12.7 Add 200μl of 70% ethanol to wash the pellet. Carefully remove and discard ethanol, being careful not to touch the bead pellet.

12.8 Spin to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible.

- 12.9 With the tube lid open incubate for 00:00:30 or until the pellet is dried (if the pellet dries completely it will crack and become difficult to resuspend).
- 12.10 Resuspend pellet in 30µl Elution Buffer(EB), mix gently by either flicking or pipetting and incubate for 00:02:00.
- 12.11 Place on magnet stand and transfer sample to a clean 1.5mL Eppendorf tube ensuring no beads are transferred into this tube.
- 12.12 Quantify 1µl of the barcoded amplicons using the Qubit Fluorometer.

Equipment

Qubit	NAME
Fluorometer	TYPE
Invitrogen	BRAND
Q33228	SKU
https://www.thermofisher.com/order/catalog/product/Q33228 ^{LINK}	



 Qubit® dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**

Adapter Ligation

25m

- 13 Prepare the following reaction for AMII Adapter ligation in a PCR tube.

	A	B
	Component	Volume
	Barcoded amplicon pool	30µL



	A	B
	NEBNext Quick Ligation Reaction Buffer (2X)	40µL
	Adapter Mix (AMII)	5µL
	Quick T4 DNA Ligase	5µL
	Total	80µL

13.1 Incubate at  25 °C for 00:20:00, possibly in a thermo cycler.

Post adapter ligation Clean-up

14 Add 80µl (1:1) of AMPure beads to the sample tube and mix gently by either flicking or pipetting. Spin to collect all liquid at the bottom of the tube.

14.1 Incubate for 00:05:00 at room temperature.

14.2 Place on magnetic rack and incubate for 00:02:00 or until the beads have pelleted and the supernatant is completely clear. Carefully remove and discard the supernatant, being careful not to touch the bead pellet.

14.3 Add 200µl SFB and resuspend beads completely by pipette mixing. Brief spin to collect all liquid at the bottom of the tube. Remove supernatant and discard.

14.4 Repeat step 14.3 to perform a second SFB wash.

14.5 Brief spin and remove any residual SFB. Add 15µl EB (ONT) and resuspend beads by gentle flicking or pipette mixing.

- 14.6 Incubate at room temperature for 00:02:00 and place on magnetic stand until clear. Transfer final library to a new 1.5mL LoBind tube
- 14.7 Quantify 1µl of the barcoded amplicons using the Qubit Fluorometer.

Equipment

Qubit

NAME

Fluorometer

TYPE

Invitrogen

BRAND

Q33228

SKU

<https://www.thermofisher.com/order/catalog/product/Q33228>^{LINK}



 Qubit® dsDNA HS Assay Kit Thermo Fisher Scientific Catalog #Q32854

Flow cell Priming and Sequencing

- 15 Thaw the following reagents.

Sequencing buffer (SQB)

Loading beads (LB)

Flush buffer (FB)

Flush tether (FLT)

- 15.1 Add 30µl FLT to the FB tube and mix well by vortex and spin down.

- 15.2 Prepare the following library dilution for sequencing in a LoBind tube:

	A	B
	SQB	37.5µL



	A	B
	LB	25.5μL
	Library	12-14μL
	Total	75-80μL

 Ligation sequencing kit 1D **Oxford Nanopore Technologies Catalog #SQK-LSK109**

 Flow Cell Priming Kit **Oxford Nanopore Technologies Catalog #EXP-FLP002**

 ONT MinION Flow Cell R9.4.1 **Oxford Nanopore Technologies Catalog #FLO-MIN106D**

15.3 Flowcell Priming, loading of the library and sequencing are performed based on the following protocol:

[dx.doi.org/10.17504/protocols.io.bbmuik6w](https://doi.org/10.17504/protocols.io.bbmuik6w)

Citations

Step 8

Quick J, Grubaugh ND, Pullan ST, Claro IM, Smith AD, Gangavarapu K, Oliveira G, Robles-Sikisaka R, Rogers TF, Beutler NA, Burton DR, Lewis-Ximenez LL, de Jesus JG, Giovanetti M, Hill SC, Black A, Bedford T, Carroll MW, Nunes M, Alcantara LC Jr, Sabino EC, Baylis SA, Faria NR, Loose M, Simpson JT, Pybus OG, Andersen KG, Loman NJ. Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples.

<https://doi.org/10.1038/nprot.2017.066>