



Jun 27, 2023

Version 2

Mitochondrial Antigen Presentation in RAW macrophages V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmk2x5g3p/v2>



Ahmed Fahmy^{1,2}, Tyler Cannon³

¹Faculté de médecine - Département de pathologie et biologie cellulaire, Université de Montréal, Québec, Canada;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815;

³Department of Microbiology/Immunology, McGill University, Québec, Canada



Lilia Rodriguez

Université de Montréal

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmk2x5g3p/v2>

Protocol Citation: Ahmed Fahmy, Tyler Cannon 2023. Mitochondrial Antigen Presentation in RAW macrophages. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvmk2x5g3p/v2> Version created by [Lilia Rodriguez](#)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 29, 2021

Last Modified: May 31, 2024

Protocol Integer ID: 55491

Keywords: antigen, presentation, Ag, ASAPCRN, mitochondrial antigen presentation in raw macrophage, mitochondrial antigen presentation, day mitochondrial antigen presentation, herpes simplex virus, raw macrophage, antigen presentation, mitochondrial matrix, murine macrophage cell line, hsv1, cell hybridoma, galactosidase assay kit

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP 000525

Disclaimer

Protocol Particulars Video

The video below is a supplement with extra context and tips, as part of the Aligning Science Across Parkinson's (ASAP) Protocol Particulars video interview series, featuring conversations with protocol authors.

<https://www.youtube.com/embed/92Vy76AakF0?si=Jmug8HphQanjC0pH>

Abstract

This protocol details methods for 3-day Mitochondrial Antigen Presentation Assay in a murine macrophage cell line (RAW) that expresses a glycoprotein B (gB) of herpes simplex virus 1 (HSV1) targeted to the mitochondrial matrix (mito-gB). A gB-specific CD8⁺ T cell hybridoma recognizing the gB_{498–505} peptide loaded on MHC class I molecules is also used to monitor antigen presentation through a beta-galactosidase assay kit.



Materials

CELL LINES

- H-2K^bRAW macrophage cell line
- β -galactosidase-inducible HSV gB/K^b-restricted HSV-2.3.2E2 CD8⁺ T cell hybridoma (2E2)

REAGENTS

- DMEM media with 10% [v/v] fetal calf serum, penicillin (FCS) [100 units/ml], and streptomycin [100 μ g/ml])
- RPMI-1640 medium supplemented with 5% (v/v) FCS, glutamine (2 mM), penicillin (100 units/ml), and streptomycin (100 μ g/ml).
- LPS stock (5mg/ml)
- PFA
- PBS
- Wash buffer: DMEM with 0.1M glycine + 10% iFBS
- Stock peptide
- 1x lysis buffer (stock = 5x, with triton, pH = 7.8) in dH₂O
- 1M DTT
- β -galactosidase Assay (CPRG) kit containing lysis buffer, CPRG buffer and CPRG substrate

CONSUMABLES

- 96-well plates
- 50ml conical tubes

EQUIPMENT

- Incubator: 37°C with 5% CO₂
- Centrifuge
- Cell counter
- Plate reader

Safety warnings

 Please refer to Safety Data Sheets (SDS) for health and environmental hazards.



Day 1: cell plating

1d

- 1 Prepare a 24 well plate:
 - Label appropriately (experimental condition – Heath shock (HS), controls, infection controls, etc.)
 - Duplicate labelled wells for the peptide control
- 2 Scrape RAW cells and resuspend (pipette up and down 8x). 
- 3 Count RAW cells and adjust concentration to **1 million cells/ml**.
- 4 Add  1 mL RAW cells to each well. 
- 5 Check if you have good 2E2 hybridoma cells and passage them if necessary.
- 6 Incubate at  37 °C with 5% CO₂ . 

Day 2: Treatment, fixation and hybridoma incubation

1d

- 7 Treatments (Prior to beginning, check the cell's media red/orange = good, yellow = bad).
- 7.1
 1. Heath shock treatment: incubate the plate in the water bath at  42 °C for  00:45:00 7h 45m
 2. Place the plate at  37 °C for  04:00:00
 3. Take the cells and re-plate in a 96 well plate. 1 well of a 24 well plate will be distributed in 6 wells = 3 wells of  250 µL and 3 wells of  50 µL for peptide
 4. Incubate the cells at  37 °C for  03:00:00
- 7.2
 1. LPS treatment:  1.5 µL per ml (dilute  20 µL of LPS stock in  666 µL of media, then add  10 µL per  1 mL of media) 6h
 2. Place the plate at  37 °C for  03:00:00



3. Take the cells and re-plate in a 96 well plate. 1 well of a 24 well plate will be distributed in 6 wells = 3 wells of  250 μL and 3 wells of  50 μL for peptide
4. Incubate the cells at  37 $^{\circ}\text{C}$ for  03:00:00

7.3

1. Bacteria treatment: bacterial infection is MOI 1. Usually  3 μL in  1 mL then add  100 μL per well
2. Place the plate at  37 $^{\circ}\text{C}$ for  03:00:00
3. Take the cells and re-plate in a 96 well plate. 1 well of a 24 well plate will be distributed in 6 wells = 3 wells of  250 μL and 3 wells of  50 μL for peptide
4. Incubate the cells at  37 $^{\circ}\text{C}$ for  03:00:00

6h

8

Prepare  1 % (v/v) PFA in PBS at  Room temperature ( 7.4).

9

Discard supernatant in the 96 well plate and add  50 μL 1% PFA .



9.1

Incubate at  Room temperature for  00:10:00 .

10m



10

Prepare 2E2 cells while RAW cells are in 1% PFA:

10.1

Pour flask into 50 ml conical and centrifuge at  1500 rpm, 00:03:00 .



Note

To keep more 2E2 cells growing, add  50 mL RPMI media back into their flask and incubate at  37 $^{\circ}\text{C}$, 5% CO_2 .

10.2


Count 2E2 cells and adjust concentration to **0.4 million cells/ml**. Make enough for all the wells and split volume in 2 tubes (you will need 1 tube to add to cells in  250 μL of media and the other tube for the peptide control wells with  50 μL of media)



- 11 Add  200 μ L wash buffer to each well in the 96 well plate and discard immediately. ! Make sure to wash the cells gently so as to not lose them. Add buffer against the wall of the well.

**Note**

Wash buffer = DMEM with 0.1M glycine + 10% iFBS.

- 11.1 Repeat 2x more: Add  200 μ L wash buffer to each well in the 96 well plate and discard immediately. **(Wash 1/2)**
- 11.2 Repeat 1x more: Add  200 μ L wash buffer to each well in the 96 well plate and discard immediately. **(Wash 2/2)**
- 12 Add  250 μ L 2E2 cells per well in $\frac{1}{2}$ of the samples (the non-peptide ones). Add an extra row of 2E2 cells alone as a negative control.
- 13 Dilute stock peptide (1 μ g/ml) **5000 fold** in the 2E2 cells (final concentration:  0.2 nanogram per milliliter (ng/mL)).
- 14 Add  250 μ L peptide + 2E2 cell solution to the remaining wells. Add an extra row of 2E2 cells+peptide alone as a negative control.
- 15 Incubate at  37 $^{\circ}$ C , 5% CO₂ for  16:00:00 (NO LONGER) .



16h

**Day 3: Lysis and revealing**

5m

- 16 Prepare 1x lysis buffer (stock = 5x, with triton,  7.8) in dH₂O.
- 16.1 Add  30 μ L 1M stock of DTT in  10 mL lysis buffer .
- 17 Prepare CPRG solution ( 7.8).





- 150 μ L CPRG buffer
- 20.2 μ L water
- 0.046 mg CPRG

18 Centrifuge 96 well plate at 2200 rpm, 00:01:00 .

19 Add 50 μ L lysis buffer to each well.

20 When lysis is done, add 170 μ L CPRG solution / well.

20.1 Incubate either at Room temperature or 37°C for 00:30:00 until the sample turns dark red (37°C speeds up reaction to about ~ 20 minutes for peptide samples).

21 Transfer 150 μ L colored solution to a new plate. Stop the reactions by adding 80 μ L of stop solution to each well

Note

Take care not to transfer debris or make bubbles. Make sure all the wells are stop at the same time.

22 Read in a plate reader at a wavelength of 570–595 nm.

Note

* To stop reaction to leave overnight, incubate at 4 °C .