



Jun 02, 2023

Modified Dixon's Agar (Heitman Lab)

 In 1 collection

DOI


<https://dx.doi.org/10.17504/protocols.io.q26g7rp11vwz/v1>

Heitman Lab¹, Amy Gladfelter¹

¹Duke University



Amy Gladfelter



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7rp11vwz/v1>

Protocol Citation: Heitman Lab, Amy Gladfelter 2023. Modified Dixon's Agar (Heitman Lab). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7rp11vwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 03, 2018

Last Modified: June 02, 2023

Protocol Integer ID: 17455

Keywords: growth of malassezia species, malassezia species, malassezia, agar, heitman lab, recipe

Abstract

Media recipe for growth of Malassezia species.

Attachments



Malassezia media 1.p...

33KB



- 1 Malt extract .  36 g
- 2 Mycological peptone (Oxoid L40)  10 g
- 3 Dessicated ox bile (Oxoid) .  10 g
- 4 Agar  20 g
- 5 Tween 60 (heated before use)  10 mL
- 6 Glycerol, 50%  4 mL
- 7 dH₂O  986 mL
- 8 Sterilize by autoclaving at  121 °C for  00:15:00 .
- 9 Cool to  50 °C .