

Jul 12, 2022

Version 1

Monkeypox virus whole genome sequencing using combination of NextGenPCR and Oxford Nanopore V.1



DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj6155lk5/v1>

Matthijs Welkers¹, M. Jonges¹, antonvandenouden²

¹AmsterdamUMC location AMC, Dept. of Medical Microbiology & Infection Prevention, Amsterdam, The Netherlands;

²Molecular Biology Systems B.V., Goes, The Netherlands



Anton van den Ouden

De_Senex_Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj6155lk5/v1>

Protocol Citation: Matthijs Welkers, M. Jonges, antonvandenouden 2022. Monkeypox virus whole genome sequencing using combination of NextGenPCR and Oxford Nanopore. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj6155lk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: June 30, 2022

Last Modified: July 12, 2022

Protocol Integer ID: 65663

Keywords: Monkeypox virus, amplicon-based tiling PCR, Nanopore sequencing, Oxford Nanopore sequencing, Monkeypox, oxford nanopore rapid genomic surveillance of monkeypox virus, monkeypox virus whole genome, use of direct oxford nanopore sequencing, direct oxford nanopore sequencing, sequencing protocol, genomic surveillance, monkeypox virus, combination of nextgenpcr, oxford nanopore rapid, genome, whole genome, amplicon pool, nextgenpcr, single nanopore flowcell, based amplicon, amplicon size, pcr, consensus sequence

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://www.protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://www.protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Rapid genomic surveillance of monkeypox virus (MPXV) can provide valuable insights in order to guide public health interventions. Current sequencing protocols make use of direct Oxford Nanopore Sequencing. However, the obtained depth is a limiting factor which prevents multiplexing samples on a flowcell making sequencing very costly. Here, we provide the protocol for a PCR-based amplicon tiling approach (inspired by SARS-CoV-2 Midnight Protocol by Nikki Freed et. al. and the ARTIC network) for MPXV consisting of a total of 88 primer sets divided over 2 amplicon pools. The amplicon size is ~2,5kB. Our approach will increase the coverage (depth) significantly and allow for multiplexing up to 20 samples on a single Nanopore flowcell.

In our experience clinical samples can be successfully sequenced with CT-values <25. Homopolymer regions will remain an issue in our approach, requiring manual curation of obtained consensus sequence genomes.

Primer pool preparation

- 1 If required, resuspend lyophilised primers to a concentration of 100µM each

Note

Primers for this protocol were designed by Martin Schou Pedersen (Department of Clinical Microbiology, Rigshospitalet, Copenhagen University Hospital, Copenhagen, Denmark) using **Primal Scheme** to generate overlapping 2500bp amplicons.

Primers are available as two pre-mixed pools from **Monkeypox - MBS - Ultra Fast NEXTGENPCR®**

- 1.1 Primers used to generate 2500 bp amplicons are here:

name	pool	seq	size	%gc	tm
monkeypox-2500_1_LEFT	1	AAAAATGTGTGACC CACGACCG	22	50	62,07
monkeypox-2500_1_RIGHT	1	CCGGGAACCTTACGC TTTCAGAT	22	50	60,86
monkeypox-2500_3_LEFT	1	GTTAACGATGCGCA CAATCTCG	22	50	61,01
monkeypox-2500_3_RIGHT	1	TAGTGAGAGCGAGA GTGACAGT	22	50	60,47
monkeypox-2500_5_LEFT	1	AATAGTCTGTAGACC TTTATCGTCGT	26	38,46	59,9
monkeypox-2500_5_RIGHT	1	ACTGCTAGAATCCG GTTTCAGATG	23	47,83	60,43
monkeypox-2500_7_LEFT	1	CACTGTAAGCATGT CCGTACCA	22	50	60,53
monkeypox-2500_7_RIGHT	1	TGAGAACGAGCTCT TCAAACACT	23	43,48	60,43
monkeypox-2500_9_LEFT	1	AGGCTATGTTTCGC CCATCATC	22	50	60,99
monkeypox-2500_9_RIGHT	1	GTCCTTTACGATGA GCTCAAATGT	24	41,67	59,68



	monkeypox-2500_11_LEFT	1	TGCTACTCCATAAC TACCACGC	22	50	60, 27
	monkeypox-2500_11_RIGHT	1	AGTTTCGTCGATAGT ACTGTGTGT	24	41,67	60,1
	monkeypox-2500_13_LEFT	1	TCCTTATGAAGATGA TGTTTGGCG	24	41,67	59, 74
	monkeypox-2500_13_RIGHT	1	CCCTCCTGGAGAAC GACAGTTA	22	54,55	61,3 3
	monkeypox-2500_15_LEFT	1	TGGAAGCGAATGAT CCGGAAAA	22	45,45	60, 8
	monkeypox-2500_15_RIGHT	1	TCCGTGGTTTCTAG TGGGTGTA	22	50	60, 94
	monkeypox-2500_17_LEFT	1	ACCTTGGCTGTCTC ATTCAATAGG	24	45,83	60, 95
	monkeypox-2500_17_RIGHT	1	TGAATGGCTGTCTG CAAAAGGT	22	45,45	60, 93
	monkeypox-2500_19_LEFT	1	AGGCTTCCAAAAAT TTTTTCATCCGT	25	36	60, 84
	monkeypox-2500_19_RIGHT	1	ACGTCGCTGTAATA GACAAGGC	22	50	60, 91
	monkeypox-2500_21_LEFT	1	CCCTAGGACGAACT ACTGCCAT	22	54,55	61,4 6
	monkeypox-2500_21_RIGHT	1	TTGTGCTGCTCTTAT CGTCTGA	22	45,45	59, 95
	monkeypox-2500_23_LEFT	1	AAAAACCCTAGTAT TCTTCCATCGC	25	40	59, 79
	monkeypox-2500_23_RIGHT	1	AACGGTATGTTACG GTTTGCCA	22	45,45	60, 99
	monkeypox-2500_25_LEFT	1	CTCGCCATTTTCGAC ATCTGGAT	22	50	60, 98
	monkeypox-2500_25_RIGHT	1	CGGGACCAAATGTA GTCAAGCT	22	50	60, 8
	monkeypox-2500_27_LEFT	1	ACGCGTTCACTATC TCCAGAGA	22	50	61,1 2



	monkeypox-2500_27_RIGHT	1	TACGCACGCTTCTC CTACCTTA	22	50	61,1 2
	monkeypox-2500_29_LEFT	1	TTGACTTTTTGGTC CACTTTTCCA	24	37,5	59, 8
	monkeypox-2500_29_RIGHT	1	ATATTCGTGACACTG TGCAACG	22	45,45	59, 76
	monkeypox-2500_31_LEFT	1	TCCGGACATGATGG TAAAGACC	22	50	60, 01
	monkeypox-2500_31_RIGHT	1	AACGAATTCTGCGT CTCGTTCA	22	45,45	61,0 4
	monkeypox-2500_33_LEFT	1	TAGGCTCACCGATG ATCATTGG	22	50	60,1 4
	monkeypox-2500_33_RIGHT	1	AACACAGCATCCAA CTGAGCAT	22	45,45	61
	monkeypox-2500_35_LEFT	1	ACAGGGGCAATGTT TACCACAA	22	45,45	60, 88
	monkeypox-2500_35_RIGHT	1	CTAGACGCCACGGG GTTTAAAA	22	50	61,0 5
	monkeypox-2500_37_LEFT	1	TTGTTTCGTCAACA AGTTGGATGA	24	37,5	59, 92
	monkeypox-2500_37_RIGHT	1	CGGATACCAGAGTG ATAATTTCCGT	25	44	60, 89
	monkeypox-2500_39_LEFT	1	CCGCATTGGTGTTT CGATCTTA	22	50	61,1 8
	monkeypox-2500_39_RIGHT	1	TGAACCTGAGGCAT GGAAAAGG	22	50	61
	monkeypox-2500_41_LEFT	1	CCACAGATTCCAAT TATCAGTTGGC	25	44	60, 83
	monkeypox-2500_41_RIGHT	1	AGACGACTCTCCAA AGATAATTGGT	25	40	60, 2
	monkeypox-2500_43_LEFT	1	TGTACAGGTACCTC CATCATTAGGA	25	44	60, 73

	monkeypox-2500_43_RIGHT	1	TTGGTTGTCTGACTT CCCAGTTG	22	50	61,18
	monkeypox-2500_45_LEFT	1	TCCTGAAAACGATG ATGGCAATC	23	43,48	59,87
	monkeypox-2500_45_RIGHT	1	AACTCTTCGAAGTG AGGATCGAT	23	43,48	59,56
	monkeypox-2500_47_LEFT	1	CTCCCGGATCACGA TTTTGTCT	22	50	60,6
	monkeypox-2500_47_RIGHT	1	GAACATATAGCGAC GCCACCAA	22	50	61,23
	monkeypox-2500_49_LEFT	1	TTGCATCTACATCAT CCGTGGA	22	45,45	59,75
	monkeypox-2500_49_RIGHT	1	AATGGAAGCCGTGG TCAATAGC	22	50	61,45
	monkeypox-2500_51_LEFT	1	TCTCTGTAGTCGAC GCTCTCAA	22	50	60,79
	monkeypox-2500_51_RIGHT	1	ACGGCCGGAAATAG TTAAGAGAC	23	47,83	60,68
	monkeypox-2500_53_LEFT	1	GTTGTATGGCATTG CGCAGAAA	22	45,45	60,85
	monkeypox-2500_53_RIGHT	1	CAAGGATGGTGTTT GTGTTGGC	22	50	60,98
	monkeypox-2500_55_LEFT	1	CTGACAATGTACTG GGCCATGT	22	50	60,8
	monkeypox-2500_55_RIGHT	1	ACATCATCGGAGGA TAATACGCTAA	25	40	59,96
	monkeypox-2500_57_LEFT	1	TTGGGAGAACTTAA GCGGCAAG	22	50	61,31
	monkeypox-2500_57_RIGHT	1	AAACGATAAGAGTG GCCGCTTG	22	50	61,68
	monkeypox-2500_59_LEFT	1	AAGATTGCGGCTAA TTGCTTCG	22	45,45	60,4



	monkeypox-2500_59_RIGHT	1	GAGGGAATTGACTC GCGAAAGA	22	50	60,85
	monkeypox-2500_61_LEFT	1	ACAGAACAATTAGA GCGGCAGG	22	50	61,12
	monkeypox-2500_61_RIGHT	1	ACACGATGCGACAA TGTATAGACT	24	41,67	60,52
	monkeypox-2500_63_LEFT	1	GACGATGATGATTG ATCACTATTACACA	28	35,71	60,19
	monkeypox-2500_63_RIGHT	1	AATCCATCCATTGC CGTCTGAT	22	45,45	60,34
	monkeypox-2500_65_LEFT	1	TCAATCCCAAACCC AAAACCGT	22	45,45	61,14
	monkeypox-2500_65_RIGHT	1	CCCAGTAAGCAACT CCATAGCA	22	50	60,28
	monkeypox-2500_67_LEFT	1	ACTTTCGAGGTTATT GGTTGTGGA	24	41,67	60,77
	monkeypox-2500_67_RIGHT	1	GCATACGCTACTCC AGAGAACG	22	54,55	61,03
	monkeypox-2500_69_LEFT	1	TGATGCACTAACGA GAAAATTAGAAGG	27	37,04	60,42
	monkeypox-2500_69_RIGHT	1	ACTTAAACCACCAT CAAAAATCCATGT	27	33,33	60,7
	monkeypox-2500_71_LEFT	1	GGTGGAGTCGTTAA AGGTGACA	22	50	60,4
	monkeypox-2500_71_RIGHT	1	TGCCTTG CATGTGA TAAGACCT	22	45,45	60,21
	monkeypox-2500_73_LEFT	1	ATTGGATTACGGT GGGTCATG	22	50	61,13
	monkeypox-2500_73_RIGHT	1	TCACAGACAGCATT TGGATCCA	22	45,45	60,14
	monkeypox-2500_75_LEFT	1	ATTCGATCGTCATG GGCATAGT	22	45,45	59,88



	monkeypox-2500_75_RIGHT	1	TGTATCTGAATCCAT GTTAGTAGTAAGCA	29	34,48	60,78
	monkeypox-2500_77_LEFT	1	GTTGGGACTGACAG ATGTGTTCT	23	47,83	60,75
	monkeypox-2500_77_RIGHT	1	TGTATCGCATTCCA CCCTTTCC	22	50	60,86
	monkeypox-2500_79_LEFT	1	GATAGATCAGTGGG TGTCCATGAT	24	45,83	60,04
	monkeypox-2500_79_RIGHT	1	GTGTTGGGTACGAC CGCTTATA	22	50	60,34
	monkeypox-2500_81_LEFT	1	CACCTGATGGTCTG GACATACC	22	54,55	60,34
	monkeypox-2500_81_RIGHT	1	ACTACGTCCTTTTG CCATTGCA	22	45,45	61,26
	monkeypox-2500_83_LEFT	1	CCACATTGGCTAGA GGAATGCC	22	54,55	61,51
	monkeypox-2500_83_RIGHT	1	TGATAAGCGACGCC ATTCATGT	22	45,45	60,92
	monkeypox-2500_85_LEFT	1	ACTAAATCTCCTTC ATGCTCTCTCAC	26	42,31	60,85
	monkeypox-2500_85_RIGHT	1	ACCTGCTCGGTTAC TTCTGTGT	22	50	61,79
	monkeypox-2500_87_LEFT	1	CCAAGCTAAGCGAC TACCATCT	22	50	60,08
	monkeypox-2500_87_RIGHT	1	TGATGCAATTGTCT GACAACCTAGA	25	40	60,9

Primers for Pool 1

	name	pool	seq	size	%gc	tm



	monkeypox-2500_2_LEFT	2	TGTTCTACACCCTG ATGCTCCT	22	50	61,0 1
	monkeypox-2500_2_RIGHT	2	TCCACCCACCTTTC TTGAAATGA	23	43,48	60, 38
	monkeypox-2500_4_LEFT	2	GTAGCAGTAGTTGG TGCATGGT	22	50	60, 8
	monkeypox-2500_4_RIGHT	2	TGTGTCCTCTCCTC TTATAACATCG	25	44	60, 08
	monkeypox-2500_6_LEFT	2	AGCGTTGACTTATG GACTCTGG	22	50	60, 27
	monkeypox-2500_6_RIGHT	2	TACCTATCCAACGA CAGGCACT	22	50	61,0 7
	monkeypox-2500_8_LEFT	2	TTGCGGACATGTTA CACTCCTT	22	45,45	60, 41
	monkeypox-2500_8_RIGHT	2	ACTATGGATCCCCA CCTTGA	22	50	60, 75
	monkeypox-2500_10_LEFT	2	TCGCCGTCATTTCT CCAAAGAA	22	45,45	60, 73
	monkeypox-2500_10_RIGHT	2	TCTGTTGTTTACCA CTCAGCGG	22	50	60, 99
	monkeypox-2500_12_LEFT	2	GGAACCGTTTTCTG ACCGTACT	22	50	60, 78
	monkeypox-2500_12_RIGHT	2	AGTCAGGTCTTGAA GGCTACCA	22	50	60, 95
	monkeypox-2500_14_LEFT	2	TGATCCAAACCCTT GATCTCCTC	23	47,83	60, 06
	monkeypox-2500_14_RIGHT	2	ACGGATTTTCAGATG GCCATTGA	22	45,45	60, 54
	monkeypox-2500_16_LEFT	2	GGCTGCTCCTGTTC TTGTAGTC	22	54,55	61,1 1
	monkeypox-2500_16_RIGHT	2	GATAACGCCAAAAT CGCTGCTC	22	50	61,0 3
	monkeypox-2500_18_LEFT	2	AAATTCGCGCCAC AATTCATC	22	45,45	60, 91
	monkeypox-2500_18_RIGHT	2	TCGCCGTTTCATTT TCAACAGC	22	45,45	61,0 3

	monkeypox-2500_20_LEFT	2	AGAAATGCCAAATC TATAAGAAAAGTCC T	29	31,03	60, 27
	monkeypox-2500_20_RIGH T	2	CCTTTATCAACAAG GAAAGCGTGT	24	41,67	60, 34
	monkeypox-2500_22_LEFT	2	TCGTATTGTGGTTAT ATGGCTACAATT	27	33,33	59, 56
	monkeypox-2500_22_RIGH T	2	TGAATTGTTGCAAC GGTTTCCA	22	40,91	60, 01
	monkeypox-2500_24_LEFT	2	TCAGTCGTTCTAAC TCCTTTGCT	23	43,48	59, 93
	monkeypox-2500_24_RIGH T	2	CACGCTTCTATGTT GCCGTCTA	22	50	60, 91
	monkeypox-2500_26_LEFT	2	AGACAGAATATCGT GAACAGGTGG	24	45,83	60, 7
	monkeypox-2500_26_RIGH T	2	TGTTTCGACTGGAG AATCATCCA	23	43,48	59, 99
	monkeypox-2500_28_LEFT	2	TAACTCCAGGCCGT TTGTTTCC	22	50	61,2 5
	monkeypox-2500_28_RIGH T	2	TTGTGTACCAGAAC TCCACCTAAA	24	41,67	59, 92
	monkeypox-2500_30_LEFT	2	CTGCCACGTTAGAG GATGACAG	22	54,55	60, 92
	monkeypox-2500_30_RIGH T	2	ACTAACGTTTCTTA GCGGAGGC	22	50	60, 85
	monkeypox-2500_32_LEFT	2	CAAGACGTTAGAGA CAAGAGACGT	24	45,83	60, 63
	monkeypox-2500_32_RIGH T	2	CAACGCCACAGATT TCTGGAGA	22	50	61,0 5
	monkeypox-2500_34_LEFT	2	GCTATTTAAATGGGT GCCGCAG	22	50	60, 72
	monkeypox-2500_34_RIGH T	2	GGTGATGATCCTTG ACGGAAGA	22	50	60, 01



	monkeypox-2500_36_LEFT	2	GGCCGCCATCATGA TCCTATTC	22	54,55	61,4 4
	monkeypox-2500_36_RIGHT	2	TTACCGCCTTCTGG ATAACCTG	22	50	60, 01
	monkeypox-2500_38_LEFT	2	AGGTGGTGGAAGCTC CTATTGGA	22	50	60, 68
	monkeypox-2500_38_RIGHT	2	CACCGCTTCGAAAC CATGAAAC	22	50	61,0 9
	monkeypox-2500_40_LEFT	2	TCACGTCAGCGGCA TCTAAATT	22	45,45	60, 86
	monkeypox-2500_40_RIGHT	2	TTCATGTGAACTT TGTCCTTTCT	25	36	59, 73
	monkeypox-2500_42_LEFT	2	AGCCCGTAAATGCA ATCAGTGA	22	45,45	60, 8
	monkeypox-2500_42_RIGHT	2	GCCGTAAACCAAG CGAATACA	22	45,45	60, 02
	monkeypox-2500_44_LEFT	2	ACGTGTACTGTATCG ACCGGAT	22	50	61,1 8
	monkeypox-2500_44_RIGHT	2	ACGGGTTTCAGAAAT ATCGACGT	22	45,45	60, 01
	monkeypox-2500_46_LEFT	2	CCAAGATCAAAAGA CACGCACG	22	50	60, 84
	monkeypox-2500_46_RIGHT	2	TTGATGATGTGGAA GGGTCTGC	22	50	60, 8
	monkeypox-2500_48_LEFT	2	AGATGGGCCCCGTTT TCTGAATA	22	50	61,1 5
	monkeypox-2500_48_RIGHT	2	TGTAGCTGTTGTAGA CATAACGGTA	25	40	59, 97
	monkeypox-2500_50_LEFT	2	GCTACTTCGTCGAT GGAAACCA	22	50	60, 85
	monkeypox-2500_50_RIGHT	2	TCCTTAAATCTGGT GCCGTTGT	22	45,45	60, 41
	monkeypox-2500_52_LEFT	2	AACCAAAAAGTCAC ACGCTCCA	22	45,45	61,3 8



	monkeypox-2500_52_RIGHT	2	TTCTATGCAGGATCT CCCGAAG	22	50	59,55
	monkeypox-2500_54_LEFT	2	GAGAACATAATGCC GCCGTAGT	22	50	60,98
	monkeypox-2500_54_RIGHT	2	TGACGTACATCCAG GAGAACCT	22	50	60,74
	monkeypox-2500_56_LEFT	2	CACACACGGCAGAA AAACCATC	22	50	61,03
	monkeypox-2500_56_RIGHT	2	GTTCCGTTCCCATC ATAGTCGT	22	50	60,34
	monkeypox-2500_58_LEFT	2	GAAACGGAATCGGT AGATCGTCT	23	47,83	60,49
	monkeypox-2500_58_RIGHT	2	CATAGCGTCTCCGG ATTCCAAG	22	54,55	61,04
	monkeypox-2500_60_LEFT	2	ACTCGACGAGCTCA CGTTTAAG	22	50	60,84
	monkeypox-2500_60_RIGHT	2	GTTTCGACGATTAAC GGAGAGCA	22	50	60,91
	monkeypox-2500_62_LEFT	2	GCTTCGCGTTTAGT CTCTGGAT	22	50	60,91
	monkeypox-2500_62_RIGHT	2	TCGATGCCTGTAAA GGGGAAAC	22	50	60,8
	monkeypox-2500_64_LEFT	2	ACCATCATCATAGCA TGCGACT	22	45,45	60,14
	monkeypox-2500_64_RIGHT	2	GTGTTTGTTGCGT TATTGCCA	22	45,45	60,98
	monkeypox-2500_66_LEFT	2	TAATAAGTTCGAGG ATGCCGCC	22	50	60,73
	monkeypox-2500_66_RIGHT	2	TTTTCCATGGACTT GTTCAACGT	23	39,13	59,56
	monkeypox-2500_68_LEFT	2	ATGTCTCGTGGGGC ATTAATCG	22	50	60,99



	monkeypox-2500_68_RIGHT	2	ACCGGATTCATCGT CGTAACAA	22	45,45	60,27
	monkeypox-2500_70_LEFT	2	AGACTAGTGTATGTG GAAATGTCATAGA	28	35,71	60,24
	monkeypox-2500_70_RIGHT	2	TCGGATTATAGCTAA GGACTAGATTCTG	27	40,74	60,21
	monkeypox-2500_72_LEFT	2	GCAAAAATCAATGG GTCGTTGGAC	24	45,83	61,93
	monkeypox-2500_72_RIGHT	2	GTGACACCCATTCA TCTGGAGA	22	50	59,95
	monkeypox-2500_74_LEFT	2	TCCTTTTAGTGCTC GACAGTGT	22	45,45	59,82
	monkeypox-2500_74_RIGHT	2	ACATTGTTTGCCAC GTCTTGAT	22	40,91	59,56
	monkeypox-2500_76_LEFT	2	TCTTCCGATATCTAC AAGGATATTCCA	27	37,04	59,61
	monkeypox-2500_76_RIGHT	2	ACGGATGATCTGCA CAGAACTC	22	50	60,6
	monkeypox-2500_78_LEFT	2	CTCATGTTCTTGTGT AATCGCAGT	24	41,67	59,92
	monkeypox-2500_78_RIGHT	2	TGTTCTGCGTCATC TACATCTGA	23	43,48	59,81
	monkeypox-2500_80_LEFT	2	AGCGAGAGATCTAG CAACTAGAGT	24	45,83	60,77
	monkeypox-2500_80_RIGHT	2	TCGAGTCATTTTAC GCACGGTT	22	45,45	61,04
	monkeypox-2500_82_LEFT	2	GCTCAATCTGCCAG GATCAAGT	22	50	60,87
	monkeypox-2500_82_RIGHT	2	TCAATGGAGCAGGA AAATGGGT	22	45,45	60,41
	monkeypox-2500_84_LEFT	2	GACCTCACAACACA GTGCAAGA	22	50	61,18

	monkeypox-2500_84_RIGHT	2	CCAGCTAACATAAG AGCCAATCTCA	25	44	61,07
	monkeypox-2500_86_LEFT	2	AAAACCATGATGTG ATAAAGCTCTGT	26	34,62	60,01
	monkeypox-2500_86_RIGHT	2	CCATTGGATGGTGC ATGTGGT	21	52,38	61,34
	monkeypox-2500_88_LEFT	2	CCGGGAACCTTACGC TTTCAGAT	22	50	60,86
	monkeypox-2500_88_RIGHT	2	AAAAATGTGTGACC CACGACCG	22	50	62,07

Primers pool 2

- 2 If you have ordered each primer independently and need to generate primer pool stocks: add  5 μL of each primer from Pool 1 to a  1.5 mL Eppendorf labeled "Pool 1 (100 μM)" and each primer from Pool 2 to a  1.5 mL Eppendorf labelled "Pool 2 (100 μM)". These are your 100 μM stocks of each primer pool.

Note

Primers should be diluted and pooled in the **mastermix cabinet** which should be cleaned with decontamination wipes and UV sterilised before and after use.

Multiplex PCR

42m

- 3 In the mastermix hood set up the multiplex PCR mastermix as follows in an  1.5 mL Eppendorf PCR tube:

Components

Nuclease Free Water

volume (pool 1)

 3.94 μL

Volume (pool 2)

 3.94 μL

Primer Pool 1 (100 μM)

 1.06 μL

Primer Pool 2 (100 μM)

 1.06 μL



NextGenPCR™ Arctic Fox HF Chemistry-2x

10 µL

10 µL

Total

15 µL

15 µL

- 4 Add 5 µL of each DNA sample to the NextGenPCR microplate containing 15 µL

Pool 1. Mix well by pipetting

- Add 5 µL of each DNA sample to the NextGenPCR microplate containing 15 µL

Pool 2. Mix well by pipetting**Note**

The **extraction and sample addition cabinet** should be cleaned with decontamination wipes and UV sterilized before and after use.

- 5 Set-up the following program on the NextGenPCR™ Instrument #10001, (Molecular Biology Systems B.V., The Netherlands):

42m

Step	Temperature	Time	Cycles
Heat Activation	98 °C	00:01:00	1
Denaturation	98 °C	00:00:10	35
Annealing and Extension	65 °C	00:01:00	35

Total PCR time is **42** minutes**Pooling**

- 6 Label a 0.2 mL PCR tube for each sample and combine the two pools from the individual PCR reaction as follows:

Component	Volume
-----------	--------



Pool 1 PCR reaction	 10 μ L
Pool 2 PCR reaction	 10 μ L
Total	 20 μ L

individual sample Bead cleaning (Optional)

- 7 Ampure XP Bead Cleanup. Add a total of  30 μ L of beads to  20 μ L of pooled samples.
- 7.1 Vortex or resuspend beads thoroughly to ensure they are well resuspended, the solution should be a homogenous brown colour.
- 7.2 Incubate for 5 minutes at room temperature
- 7.3 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear. 2m
- 7.4 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 7.5 Add  150 μ L of freshly prepared room-temperature [M] 80 % volume volume ethanol to the pellet.
- 7.6 Keeping the magnetic rack on the benchtop, rotate the bead-containing tube by 180°. Wait for the beads to migrate towards the magnet and re-form a pellet. Remove the ethanol using a pipette and discard.
- 7.7 Repeat step 7.5 and 7.6
- 7.8 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette
- 7.9 With the tube lid open incubate for  00:01:00 or until the pellet loses it's shine (if the pellet dries completely it will crack and become difficult to resuspend) 1m



7.10 Remove the tube from the magnetic rack. Resuspend Pellet in  20 μL Nuclease free water mix gently by flicking and incubate at room temperature for  00:02:00

2m

7.11 Place on magnet for beads to collect on the magnetic side of the holder and the supernatant to become clear and then transfer  18 μL sample to a clean 96 wells Microplate ensuring no beads are transferred into this tube.

Quantification

8 Prepare a mastermix of Qubit™ working solution for the required number of samples and standards. The Qubit dsDNA kit requires 2 standards for calibration.

8.1

Per sample:

Qubit® dsDNA HS Reagent  1 μL μL

Qubit® dsDNA HS Buffer  199 μL μL

8.2 Aliquot Qubit™ working solution to each tube:

- standard tubes requires 190 μL of Qubit™ working solution
- sample tubes require anywhere from  180 μL to  199 μL (depending how much sample you wish to add)

The final volume in each tube must be 200 μL once sample/standard has been added.

8.3 Add 10 μL of standard to the appropriate tube.

8.4 Add 1–20 μL of each user sample to the appropriate tube.

8.5 Mix each tube vigorously by vortexing for 3–5 seconds.

- 8.6 Allow all tubes to incubate at room temperature for 2 minutes, then proceed to “Read standards and samples”.
- 8.7 On the Home screen of the Qubit™ 3 Fluorometer, press DNA, then select 1X dsDNA HS as the assay type. The Read standards screen is displayed. Press Read Standards to proceed.
- 8.8 Insert the tube containing Standard #1 into the sample chamber, close the lid, then press Read standard. When the reading is complete (~3 seconds), remove Standard #1.
- 8.9 Insert the tube containing Standard #2 into the sample chamber, close the lid, then press Read standard. When the reading is complete, remove Standard #2.
- 8.10 The instrument displays the results on the Read standard screen. For information on interpreting the calibration results, refer to the Qubit™ Fluorometer User Guide, available for download at thermofisher.com/qubit.
- 8.11 Press Run samples.
- 8.12 On the assay screen, select the sample volume and units:
 - Press the + or – buttons on the wheel, or anywhere on the wheel itself, to select the sample volume added to the assay tube (from 1–20µL).
 - From the unit dropdown menu, select the units for the output sample concentration (in this case choose ng/µL).
- 8.13 Insert a sample tube into the sample chamber, close the lid, then press Read tube. When the reading is complete (~3 seconds), remove the sample tube.
- 8.14 The top value (in large font) is the concentration of the original sample and the bottom value is the dilution concentration. For information on interpreting the sample results, refer to the Qubit™ Fluorometer User Guide.
- 8.15 Carefully record all results and store run file from the Qubit on a memory stick.



- 8.16 All negative controls should ideally be 'too low' to read on the Qubit machine, but MUST be < 1ng per ul. If your negative controls >1ng per ul, considerable contamination has occurred and you must redo previous steps.

Normalisation

- 9 Adjust the amount of DNA in the tube to be 100 ng total per sample in 7.5 μ L molecular grade water. For example if your PCR reaction is at 100ng/ul, add 1ul of the PCR reaction to 6.5ul of molecular grade water. Input to the Rapid Barcoding kit will vary depending on the amplicon length but we have determined 50-200 ng works for efficient barcoding of this amplicon length. If there is under 100ng or you do not know the concentration, simply use all 7.5 μ L of the pooled PCR reaction. Use the full 7.5 μ L of the negative control, even if there is no detectable DNA in the PCR reaction.

Rapid barcoding

9m 30s

- 10 Add  7.5 μ L of each diluted PCR reaction to the labeled PCR tube.
Set up the following reaction from each sample:

Component	Volume
DNA amplicons	 7.5 μ L
Fragmentation Mix RB01-12	 2.5 μ L
Total	 10 μ L

- 10.1 Mix gently by flicking the tube, and spin down.

- 10.2 Incubate the reaction in a PCR machine:

2m 30s

 30 °C for  00:01:00 80 °C for  00:01:00 4 °C for  00:00:30

- 10.3 Pool all barcoded samples, noting the total volume

- 11 Bead Cleanup. Use a 1:1 ratio of sample to beads.



- 11.1 Vortex SPRI beads thoroughly to ensure they are well resuspended, the solution should be a homogenous brown colour.
- 11.2 Add an equal volume (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting. For example add  50 μL room temperature SPRI beads to a  50 μL reaction.
- 11.3 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 11.4 Incubate for  00:05:00 at room temperature. 5m
- 11.5 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear. 2m
- 11.6 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 11.7 Add  200 μL of freshly prepared room-temperature [M] 80 % volume volume ethanol to the pellet.
- 11.8 Keeping the magnetic rack on the benchtop, rotate the bead-containing tube by 180°. Wait for the beads to migrate towards the magnet and re-form a pellet. Remove the ethanol using a pipette and discard.
- 11.9 And repeat ethanol wash (steps 11.7-11.8)
- 11.10 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette
- 11.11 With the tube lid open incubate for  00:01:00 or until the pellet loses it's shine (if the pellet dries completely it will crack and become difficult to resuspend) 1m
- 11.12 Remove the tube from the magnetic rack. Resuspend Pellet in  10 μL 10 mM Tris-HCl pH 8.0 with 50 mM NaCl, mix gently by flicking and incubate at room temperature for  00:02:00 2m



11.13 Place on magnet and transfer sample to a clean  1.5 mL Eppendorf tube ensuring no beads are transferred into this tube.

11.14 Add  1 μ L of RAP (from the SQK-RBK110.96 kit) to  10 μ L cleaned, barcoded DNA, mix gently by flicking the tube, and spin down.

11.15 Incubate the reaction for  00:05:00 at room temperature.

5m

11.16 The prepared library is used for loading into the MinION flow cell according to Oxford Nanopore Rapid Barcoding (SQK-RBK110.96) protocol. Please refer to the Oxford Nanopore Rapid Barcoding SQK-RBK110.96 protocol at this stage. Store the library on ice until ready to load.

MinION sequencing

12 Start the sequencing run using MinKNOW.