

Feb 16, 2026

Multi-viral target Direct Detection and Nanopore Sequencing (DDNS)

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl48r8rvo5/v1>

Alison Tedcastle¹, Erika Bujaki¹, Thomas Wilton¹, Javier Martin¹

¹MHRA



Alison Tedcastle

MHRA

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl48r8rvo5/v1>

Protocol Citation: Alison Tedcastle, Erika Bujaki, Thomas Wilton, Javier Martin 2026. Multi-viral target Direct Detection and Nanopore Sequencing (DDNS). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl48r8rvo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 09, 2026



Last Modified: February 16, 2026

Protocol Integer ID: 242917

Keywords: Norovirus, Rotavirus, Hepatitis A, Hepatitis E, SARS-CoV-2, PMMoV, RT-PCR, Nanopore sequencing, Amplicon sequencing, Multi-pathogen, minION, PCR, Environmental surveillance, Wastewater surveillance, DDNS, Enterovirus D68, EVD68, Enterovirus 71, EV71, multiple viral pathogens from wastewater, magmax viral rna isolation large volume extraction protocol, viral rna extraction, magmax viral rna isolation, sensitive identification of multiple different viral pathogen, viral rna, nanopore sequencing, extracted rna, wastewater sample, multiple viral pathogen, multiple different viral pathogen, wastewater, enterovirus d68, enterovirus, oxford nanopore, nanopore, sensitive nested pcr approach, norovirus gi, extraction, large volume extraction protocol, target specific amplification for direct detection

Abstract

This protocol describes the target specific amplification for direct detection of multiple viral pathogens from wastewater extracted RNA. This enables sensitive identification of multiple different viral pathogens of interest for both surveillance and research purposes.

Targets include Enterovirus D68, Enterovirus 71, Rotavirus, Hepatitis A and E, Norovirus GI and GII and SARS-CoV-2.

Wastewater samples should be concentrated using an appropriate method and viral RNA extracted following the method; MagMAX Viral RNA Isolation Large volume extraction protocol using 1.2 mL of concentrate for each extraction (see Viral RNA extraction below). This method then utilises highly-sensitive nested PCR approaches where amplicons are subsequently pooled and purified ready for library preparation for Oxford Nanopore or Illumina sequencing. Barcoding the amplicons prior to nanopore sequencing allows multiplexing of many samples, reducing turn around time and cost per sample.

Materials

- Nuclease-Free Water (Promega, Cat. No. P1193 or Merck, Cat. No. 3098)
- SuperScriptTM III One-Step RT-PCR System with PlatinumTM Taq High Fidelity DNA Polymerase (Invitrogen, Cat. No. 12574035)
- DreamTaq Hot Start PCR Master Mix (Thermo Scientific, Cat. No. K9011)
- Ethanol (96-100%)
- 0.2ml PCR tube strips
- 1.5 mL LoBind tubes (Eppendorf, Cat. No. 022431021)
- Gel electrophoresis agarose gel kit and buffers or TapeStation DNA kit
- Molecular weight marker (such as 1kb DNA Ladder, Promega, Cat. No. G571A)
- AMPure magnetic beads (Beckman Coulter, Cat. No. A63880)
- Qubit Broad Range dsDNA kit (Thermo Scientific, Cat. No. Q33853)



Before start

Viral RNA extraction

We recommend using the MagMAX Viral RNA Extraction method. Please follow the appropriate protocol depending on your starting material and manual or automated workflow:

Purification of viral RNA from environmental samples using MagMAX Viral RNA Isolation Kit - <https://www.protocols.io/view/large-volume-viral-rna-extraction-using-magmax-vir-81wgbzdr3gpk/v3>

Using freshly extracted viral RNA is recommended whenever possible to facilitate successful amplification of long genome fragments.

Repeated freeze-thaw of samples or purified RNA can affect viral nucleic acid integrity, reducing amplification efficiency.



Multi-viral RT-PCR

- 1 Wipe down lab bench surfaces and pipettes with 70% ethanol and an RNase cleaning agent, such as RNase Zap.
- 2 Prepare separate individual master mix preparations for each target using Superscript III One-Step RT-PCR System with Platinum Taq High Fidelity for each target (See Appendix 1). The reaction volumes are detailed in the tables below and should be prepared for the number of samples together with positive, negative controls plus 10% (for pipetting loss), where appropriate.

Note

Volume of nuclease free water can be reduced to allow larger volumes of RNA to be processed. **Each sample should be tested in duplicate.**

Note

PPMoV is used as an internal control for the monitoring of sample quality.

3

	A	B
	Target	Volume (μL)
	EV D68, EV71 or Hep E	
	2x Reaction Mix	12.5
	Anti-sense primer	1
	Enzyme Mix	1
	Water	4.5
	Aliquot	19
	RNA	5
	Total volume	24



4

	A	B
	Target	Volume (µL)
	Rotavirus, Hep A, Norovirus or PMMoV	
	2x Reaction Mix	12.5
	Anti-sense primer	1
	Sense primer	1
	Enzyme Mix	1
	Water	4.5
	Aliquot	20
	RNA	5
	Total volume	25

5

	A	B
	Target	Volume (µL)
	SARS-CoV-2	
	2x Reaction Mix	12.5
	Anti-sense primer	0.5
	Sense primer	0.5
	Enzyme Mix	1
	Water	5.5
	Aliquot	20
	RNA	5
	Total volume	25



- 6 Vortex the master mix briefly and spin down to gather contents at the bottom of the tube. Aliquot the required volume of master mix as shown in the individual tables above (depending on reaction) into each pre-labelled PCR tube/strip.
- 7 Close the lids of the PCR tubes/strips and transfer the rack to the template addition room.
- 8 Add 5 μ L of positive control or sample RNA to the appropriate reaction tubes. Add 5 μ L of nuclease free water for the negative control.
- 9 Place PCR sample tubes in the Thermocycler and run the below conditions relevant to the target.

10

A	B
Profile: EV D68, EV71, Hep E	
50°C for 30 min	
Remove tubes, spin down sample and add 1 μ L of forward primer into each tube.	
94°C for 2min	
Cycles 42	Step 1: 94°C for 15 secs
	Step 2: 55°C for 30 secs
	Step 3: 68°C for 1 min 30 secs
68°C for 5 min	
10°C hold	

11

A	B
Profile: Rotavirus, Hep A, Norovirus, PMMoV, SARS-CoV-2	
50°C for 30 min	

A	B
94°C for 2min	
Cycles 42	Step 1: 94°C for 15 secs
	Step 2: 55°C for 30 secs
	Step 3: 68°C for 1 min 30 secs
68°C for 5 min	
10°C hold	

Note

At this point you can store the PCR products at -20°C. These can be stored indefinitely before continuing to step 12, however it is advised to continue the PCR protocol as soon as possible.

Second round PCR

- 12 Wipe down lab bench surfaces and pipettes with 70% ethanol and an DNase cleaning agent, such as DNase Zap.
- 13 Prepare a separate master mix for each target using DreamTaq Hot Start PCR Master Mix (See Appendix 1). The reaction volumes are detailed in the tables below and should be prepared for the number of samples together with positive, negative controls plus 10% (for pipetting loss), where appropriate.

14

A	B
All targets	Volume (µL)
2x Reaction Mix	12.5
Forward primer	1
Reverse primer	1
Water	8.5
First PCR product	2
Total volume	25



- 15 Vortex the master mix briefly and spin down to gather contents at the bottom of the tube. Aliquot 23 μ L of the master mix into each PCR tube/strip.
- 16 Close the lids of the PCR tubes/strips and transfer the rack to the template addition room.
- 17 Add 2 μ L of the first round PCR product to the appropriate reaction tubes. Add 2 μ L of nuclease free water for the negative control.
- 18 Place PCR sample tubes in the Thermocycler and run the below conditions.
- 19 Run the following cycling conditions.

A	B
Profile: All targets	
94°C for 2min	
Cycles 40	Step 1: 94°C for 15 secs
	Step 2: 50°C for 30 secs
	Step 3: 68°C for 1 min
68°C for 5 min	
10°C hold	

- 20 Check all samples including positive and negative controls from the second round PCR on a 1% agarose gel. The expected band size for each target is in Appendix 1. Proceed to the next step with the confirmed positive reactions for purification and quantification of the DNA.

Amplicon Pooling

- 21 Pool 10 μ L from each duplicate amplicon or 20 μ L from the only positive reaction (depending on successful amplification as noted in the gel) for each virus sample.

Note

We present a simplified approach here directed by electrophoresis results with purifying and quantifying the pooled reactions. Alternatively, each amplicon can be purified and quantified individually for accurate equimolar pooling.

Preparing input DNA

- 22 Allow AMPure beads to warm to room temperature then add an equal volume of beads to the individually pooled PCR samples and pipette gently to mix.
For example, for a 20µL pooled sample add 20µL of resuspended beads
- 23 Incubate at room temperature for 5 minutes.
- 24 Spin down and place on a magnetic rack until clear and colourless, then pipette off the supernatant.
- 25 Whilst on the magnet, wash the beads with 100µL of 80% Ethanol, remove and discard.
- 26 Repeat Step 25.
- 27 Spin down the tubes and place back on the magnet to remove any residual Ethanol.
Allow to air dry for about 30 seconds.
- 28 Remove the tubes from the magnetic rack and resuspend the pellet in 20µL nuclease free water and incubate at room temperature for 2 minutes.
- 29 Pellet the beads on the magnet and remove 20µL (or as much as possible) of each sample into a clean tube.
- 30 Quantify the purified DNA with a Qubit Broad Range dsDNA kit.
- 31 In brief, create a master mix of 200µL (199µL buffer, 1µL Qubit reagent) for each sample + 2 standards + 10 %.



- 32 For the two standards, add 190 μ L master mix to 10 μ L of standard, and for samples add 198 μ L master mix to 2 μ L sample.
- 33 Vortex each sample or standard once they are added to the master mix and incubate for 2 minutes before quantifying.

Note

Important - be sure to select the broad range kit on the Qubit

- 34 Record the concentration of DNA for each sample.
- 35 Transfer 200 fmol of DNA per sample into a clean 96-well plate and adjust volume to 12.5 μ L with nuclease free water. Mix gently by pipetting.
For example, 200 fmol of a 900bp fragment = ~110ng*

Note

Each amplicon is different in size; therefore, each will require different amounts of DNA to equal 200 fmol.

- 36 Once aliquoted, spin down briefly.
- 37 Continue to step 5 in the native barcoding protocol for direct detection nanopore sequencing. <https://dx.doi.org/10.17504/protocols.io.e6nvwjko9lmk/v1>

Appendix 1: Primers (All primers should be at a working concentration of 10 μ M)

38



A	B	C	D	E
Target	Primer name		Sequence	Reference
EVD-68 (~900bp)	EVD-68 F0(1)		YGCTAARTCAATYAATGC YAATGTYGG	
	EVD-68 R0(1)		CTTRCATCTTGCTATYTG RTGTYTTCC	
	EVD-68 F1(2)		CATGRRGCAGAGGCAGC YTACC	
	EVD-68 R1(2)		CAAARACYCCTCCGAAG CCTG	
EV71 (~900bp)	EV71 F0(1)		TCHGTYACCCTTGTRATA CCRTGG	
	EV71 R0(1)		GARTTRCARTARTACACY CCTGTYTG	
	EV71 F1(2)		GGRGAYAGRGTGGCMGA TG	
	EV71 R1(2)		GTCYTCCCADACDAGRT TYGC	
Rotaviruses (~900bp)	Beg9(1)		GGCTTTAAAAGAGAGAA TTTCCGTCTGG	1, 2
	VP7R(1,2)		AACTTGCCACCATTTTTT CC	
	VP7F(2)		ATGTATGGTATTGAATATA CCAC	
Hepatitis A (~550bp)	2820F(1)		YACAAGRAGAACAGGRA AYATTCAG	
	3788R(1)		BACAGARTGATGRAAHC CAAACATC	
	2904F(2)		GGAGAYAARACRGATTCH ACWTTTGG	
	3469R(2)		CATCYTTCATTTTCYGTCC ATTTYTCATC	
Hepatitis E (450bp)	HE040_R(1)		CCCTTRTCCTGCTGAGC RTTCTC	3
	HE044_F(1)		CAAGGHTGGCGYTCKGT TGAGAC	



	A	B	C	D	E
		HE110.2_F(2)	HE110-2.1*	GYTCKGTTGAGACCTCY GGGGT	
			HE110-2.2*	GYTCKGTTGAGACCACG GGYGT	
			HE110-2.3*	GYTCKGTTGAGACCTCT GGTGT	
		HE041_R(2)		TTMACWGTCRGCTCGC CATTGGC	
	Norovirus GI (~300bp)	COG1F(1)		CGYTGGATGCGNTTYCA TGA	4, 5
		G1-SKR(1,2)		CCAACCCARCCATTRTAC A	
		G1SKF(2)		CTGCCCCGAATTYGTAAT GA	
	Norovirus GII (~300bp)	COG2F(1)		CARGARBCNATGTTYAGR TGGATGAG	4, 5
		G2-SKR(1,2)		CCRCCNGCATRHCCRTT RTACAT	
		G2SKF(2)		CNTGGGAGGGCGATCGC AA	
	SARS-CoV-2 (~600bp)	Covid-501minus2 F(1)		GATCTCTGCTTTACTAAT GTCTATGCAG	
		Covid-501R0alt2(1)		CCTGGTGTTATAAACTG ACACCA	
		Covid-501minus1F B(2)		GAGGTGATGAAGTCAGA CAAATYGC	
		Covid-501R4alt(2)		TCAAGAATSTCAAGTGTC TGTGGATC	
	PMMoV (~500bp)	CP/sF(1)		ATGGCATACACAGTTACC AGT	6
		CP/aR(1)		TTAAGGAGTTGTAGCCCA CGTA	

(1) Primers used for Primary RT-PCR

(2) Primers used for Second round PCR

* Prepare primer mixes to a working concentration of 10µM:

1. HE110-2.1, HE110-2.2 and HE110-2.3



Protocol references

1. Gouvea, V. et al. (1990). Polymerase chain reaction amplification and typing of rotavirus gene 9 (VP7). *Journal of Clinical Microbiology*, 28(2), 276–282.
2. M. Iturriza-Gómara, D. Cubitt, U. Desselberger, J. Gray. Amino acid substitution within the VP7 protein of G2 rotavirus strains associated with failure to serotype.
3. Hitoshi Mizuo et al. Polyphyletic Strains of Hepatitis E Virus Are Responsible for Sporadic Cases of Acute Hepatitis in Japan. *J Clin Microbiol*. 2002 Sep; 40(9): 3209–3218.
4. Kojima, S., T. Kageyama, S. Fukushi, F. B. Hoshino, M. Shinohara, K. Uchida, K. Natori, N. Takeda, and K. Katayama. 2002. Genogroup-specific PCR primers for detection of Norwalk-like viruses. *J. Virol. Methods* 100:107–114.
5. Santiso-Bellón, C., Randazzo, W., Pérez-Cataluña, A., Vila-Vicent, S., Gozalbo-Rovira, R., Muñoz, C., Buesa, J., Sanchez, G., & Rodríguez Díaz, J. (2020). Epidemiological Surveillance of Norovirus and Rotavirus in Sewage (2016–2017) in Valencia (Spain). *Microorganisms*, 8(3), 458.
6. Çaglar BK, Fidan H, Elbeaino T. 2013. Detection and molecular characterization of Pepper mild mottle virus from Turkey. *J. Phytopathol*. 161: 434 –438.
7. MagMAX Viral RNA Isolation Large volume extraction protocol
[dx.doi.org/10.17504/protocols.io.81wgbzdr3gpk/v3](https://doi.org/10.17504/protocols.io.81wgbzdr3gpk/v3)
8. Native Barcoding of amplicons for DDNS
[dx.doi.org/10.17504/protocols.io.e6nvwjko9lmk/v1](https://doi.org/10.17504/protocols.io.e6nvwjko9lmk/v1)