

Nov 08, 2023

Neuronal transdifferentiation from human primary adult fibroblasts

 [Nature Cell Biology](#)

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlko6d5v5r/v1>

Ching-Chieh Chou¹, Judith Frydman¹

¹Department of Biology, Stanford University



Ching-Chieh Chou

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlko6d5v5r/v1>

External link: <https://doi.org/10.1038/s41556-025-01623-y>

Protocol Citation: Ching-Chieh Chou, Judith Frydman 2023. Neuronal transdifferentiation from human primary adult fibroblasts. protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlko6d5v5r/v1>

**Manuscript citation:**

Chou, CC., Vest, R., Prado, M.A. *et al.* Proteostasis and lysosomal repair deficits in transdifferentiated neurons of Alzheimer's disease. *Nat Cell Biol* **27**, 619–632 (2025). <https://doi.org/10.1038/s41556-025-01623-y>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 08, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 90649

Keywords: ASAPCRN, Neuron, Transdifferentiation, neuronal transdifferentiation from human primary adult fibroblast, human primary fibroblasts into neuron, converted neuron, neuronal transdifferentiation, human primary adult fibroblast, human primary fibroblast, cortical glutamatergic neuron, characteristics of cortical glutamatergic neuron, neuron, direct conversion, protocol for the direct conversion

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This is a protocol for the direct conversion of human primary fibroblasts into neurons using a combination of transcription factor- and small molecule-based approach. The majority of converted neurons showing characteristics of cortical glutamatergic neurons.

Preparation of culture medium and reagents

1. Fibroblast medium @4°C lasts for 1 month. A total of 500 mL:
 - 435 ml DMEM (High glucose, GlutaMAX) (Thermo Fisher Scientific)
 - 50 ml Fetal bovine serum (FBS) (Thermo Fisher Scientific)
 - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
 - 5 ml MEM NEAA (100X; Thermo Fisher Scientific)
 - 5 ml Sodium pyruvate (100 mM; Thermo Fisher Scientific)
 - 0.5 ml beta-mercaptoethanol (55mM; Thermo Fisher Scientific)
 - Sterilized by filtration through a 0.22 µm filter
2. Neuronal reprogramming medium @4°C lasts for 1 month. A total of 500 mL:
 - 240 mL DMEM/F-12 (Thermo Fisher Scientific)
 - 240 mL Neurobasal Medium (Thermo Fisher Scientific)
 - 10 mL B-27 Supplement (50X) (Thermo Fisher Scientific)
 - 5 mL N-2 Supplement (100X) (Thermo Fisher Scientific)
 - 1.25 mL GlutaMAX Supplement (Thermo Fisher Scientific)
 - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
 - Sterilized by filtration through a 0.22 µm filter
3. Neuronal maturation medium @4°C lasts for 1 month. A total of 500 mL:
 - 480 mL BrainPhys Neuronal Medium (STEMCELL Technologies)
 - 10 mL B-27 Supplement (50X) (Thermo Fisher Scientific)
 - 5 mL N-2 Supplement (100X) (Thermo Fisher Scientific)
 - 1.25 mL GlutaMAX Supplement (Thermo Fisher Scientific)
 - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
 - Sterilized by filtration through a 0.22 µm filter
4. Hexadimethrine Bromide (Polybrene) (Sigma-Aldrich)
5. Forskolin (Sigma-Aldrich)
6. Dorsomorphin (Tocris)
7. SB 431542 (Tocris)
8. XAV939 (Stemgent)
9. Doxycycline (Cayman)
10. Puromycin (Thermo Fisher Scientific)

11. BDNF Recombinant Protein (Peprotech)
12. NT-3 Recombinant Protein (Peprotech)
13. Poly-L-ornithine hydrobromide (Sigma-Aldrich)
14. Human rhLaminin-521 (Thermo Fisher Scientific)
15. Human Vitronectin (Thermo Fisher Scientific)

Transduction/induction/selection of human fibroblasts

- 2 **Day -1:** Plate human fibroblasts at 200K cells per well of a poly-L-ornithine coated 6-well plate.
 - 2.1 Dilute sterile-filtered 0.01% poly-L-ornithine in PBS at a 1-to-10 ratio.
 - 2.2 Coat the plate and incubate in a 37°C incubator for 2 hr.
 - 2.3 Wash the plate with H₂O for 5 min on a shaker and repeat for three times. Airdry the plate before plating the cells.
 - 2.4 Rinse the cells with sterile 1x PBS to remove all complete medium.
 - 2.5 Add 0.05% Trypsin-EDTA to cover the bottom of the dish. Incubate at 37°C with 5% CO₂ for ~5 minutes. Cells should round up and become dislodged.
 - 2.6 Neutralize trypsin-EDTA activity by adding the complete media.
 - 2.7 Gently pipette to resuspend the cells and transfer the cell-containing medium to a 15 mL conical tube.
 - 2.8 Centrifuge at 200xG for 5 min to pellet the cells.

- 2.9 Remove the supernatant. Cell numbers are determined for each sample by a hemocytometer.
- 2.10 Add fresh culture medium to adjust the cell concentration to 200K cells per well and each well contain 2 mL of culture medium.
- 3 **Day 0:** Infect cells with 2 mL of viral supernatants per well of a 6-well plate.
 - 3.1 Make virus dilutions by adding each viral supernatant to cold Fibroblast medium. 100 to 400 μ L of each viral supernatant (M2rtTA, Ngn2-puro, Ascl1, Brn2, Myt1l) depending on viral titers.
 - 3.2 Combine these viral supernatants and add cold Fibroblast medium to reach a total of 12 mL for one 6-well plate.
 - 3.3 Add 4 μ g/mL Polybrene to the virus dilutions.
 - 3.4 Remove old Fibroblast medium and incubate cells with virus dilutions for 24 hr.
- 4 **Day 1:** Remove virus-containing medium and add fresh Fibroblast medium plus 1 μ g/mL doxycycline.
- 5 **Day 2:** Change to fresh Fibroblast medium plus doxycycline and puromycin.
- 6 **Day 4:** Coat glass coverslips or plates with Vitronectin/rLaminin-521 and replat the transduced cells.

Note

The conversion efficiency and cell survival are much better after replating as compared with no replating.

- 6.1 Coat with vitronectin (1:100 dilution of stock concentration 0.5 mg/mL in PBS) for 1 hr at room temperature. No need to rinse for vitronectin-coating coverslips or plates.

- 6.2 Coat with rhLaminin-521 (1:100 dilution of stock concentration 100 µg/mL in DPBS with Ca^{2+} and Mg^{2+}) for 2 hr at 37°C.
- 6.3 In parallel, trypsinize cells by 0.05% Trypsin-EDTA for 5 min at 37°C and perform PSA-NCAM selection using magnetic-based cell sorting.
- 6.4 Replate the PSA-NCAM+ cells in Fibroblast medium plus doxycycline. Generally ~50K cells per cm^2 .
- 7 **Day 5:** Switch to Neuronal reprogramming medium plus small molecules.
 - 5 µM Forskolin
 - 2 µM Dorsomorphin
 - 10 µM SB 431542
 - 2 µM XAV939
 - 1 µg/mL Doxycycline
- 8 Change half of the culture medium every 2 to 3 days.
- 9 **Day 12,** Switch to Neuronal reprogramming medium supplemented with small molecules and neurotrophic factors.
 - 5 µM Forskolin
 - 2 µM Dorsomorphin
 - 10 µM SB 431542
 - 2 µM XAV939
 - 1 µg/mL Doxycycline
 - 10 ng/mL BDNF
 - 10 ng/mL NT-3
- 10 Change half of the culture medium every 2 to 3 days.
- 11 **Day 20,** switch to Neuronal maturation medium supplemented with small molecules and neurotrophic factors.
 - 5 µM Forskolin
 - 2 µM Dorsomorphin
 - 10 µM SB 431542
 - 1 µg/mL Doxycycline
 - 10 ng/mL BDNF
 - 10 ng/mL NT-3



Note

Suggest half medium change to gradually replace the neuronal reprogramming medium with maturation medium. This reduces the exposure of neurons to the air.

- 12 Change half of the culture medium every 3 to 4 days.
- 13 Neurons will mature on and after 5 weeks in culture and be ready for experiments
- 14 Neurons can be cultured for about 9 weeks.