

Dec 08, 2023

## Perturb-seq characterizing regulators of T cell function

DOI

[dx.doi.org/10.17504/protocols.io.dm6gp3qjdvzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gp3qjdvzp/v1)

Sean R. McCutcheon<sup>1</sup>, Adam M. Swartz<sup>1</sup>, Michael C. Brown<sup>1</sup>, Alejandro Barrera<sup>1</sup>, Christian McRoberts Amador<sup>1</sup>, Keith Siklenka<sup>1</sup>, Lucas Humayun<sup>1</sup>, Maria A. ter Weele<sup>1</sup>, James M. Isaacs<sup>1</sup>, Andrea R Daniel<sup>1</sup>, Timothy E. Reddy<sup>1</sup>, Andrew Allen<sup>1</sup>, Smita K. Nair<sup>1</sup>, Scott J. Antonia<sup>1</sup>, Charles A. Gersbach<sup>1</sup>

<sup>1</sup>Duke University

Andrea R Daniel: This protocol was adapted from the work of Sean McCutcheon and colleagues in the Gersbach lab at Duke University.

Gersbach Lab



Andrea R Daniel

Duke University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp3qjdvzp/v1>



**Protocol Citation:** Sean R. McCutcheon, Adam M. Swartz, Michael C. Brown, Alejandro Barrera, Christian McRoberts Amador, Keith Siklenka, Lucas Humayun, Maria A. ter Weele, James M. Isaacs, Andrea R Daniel, Timothy E. Reddy, Andrew Allen, Smita K. Nair, Scott J. Antonia, Charles A. Gersbach 2023. Perturb-seq characterizing regulators of T cell function. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp3qjdvzp/v1>

**Manuscript citation:**

McCutcheon, S.R., Swartz, A.M., Brown, M.C. *et al.* Transcriptional and epigenetic regulators of human CD8<sup>+</sup> T cell function identified through orthogonal CRISPR screens. *Nat Genet* (2023). <https://doi.org/10.1038/s41588-023-01554-0>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 02, 2023

**Last Modified:** December 08, 2023

**Protocol Integer ID:** 91713

**Keywords:** Perturb-seq, CRISPRi, CRISPRa, T cells, scRNA-seq, seq assay, characterizing transcriptional regulator, seq, perturb, cell

**Funders Acknowledgements:**

NIH

Grant ID: HG012053

## Abstract

This protocol describes methods for a Perturb-seq assay characterizing transcriptional regulators of T cell function.

## Materials

Addgene: pLV hU6-gRNA hUbC-dSaCas9-KRAB-T2A-Thy1.1 (Addgene 194278) and pLV hU6-gRNA hUbC-VP64-dSaCas9-VP64-T2A-Thy1.1 (Addgene 194279).

## Troubleshooting

## gRNA library cloning

- 1 Oligonucleotide pools containing 40-gRNA sequences and constant regions for polymerase chain reaction (PCR) amplification were synthesized by Twist Bioscience.
- 2 gRNA amplicons were gel extracted, PCR purified and input into 20 µl Gibson reactions (5:1 molar ratio of insert to backbone) with 200 ng of Esp3I digested and 1 × solid-phase reversible immobilization (SPRI)-selected (Beckman Coulter) plasmid backbone.
- 2.1 Addgene: pLV hU6-gRNA hUbC-dSaCas9-KRAB-T2A-Thy1.1 (Addgene 194278) and pLV hU6-gRNA hUbC-VP64-dSaCas9-VP64-T2A-Thy1.1 (Addgene 194279).
- 3 Gibson reactions were purified using ethanol precipitation and transformed into Lucigen's Endura ElectroCompetent Cells.
- 4 Transformed cells were cultured overnight and plasmids were isolated using Qiagen Midi Kits.

## Transfections for high-titer lentiviral production

- 5 Plate  $1.2 \times 10^6$  or  $7 \times 10^6$  HEK293T cells in a 6 well plate or 10 cm dish in the afternoon with 2 mL or 12 mL of complete opti-MEM (Opti-MEM<sub>000</sub> I Reduced Serum Medium supplemented with 1x Glutamax, 5% FBS, 1 mM Sodium Pyruvate, and 1x MEM Non-Essential Amino Acids).
- 6 The next morning, transfect HEK293T cells with 0.5 µg pMD2.G, 1.5 µg psPAX2, and 0.5 µg transgene for 6 well plates or 3.25 µg pMD2.G, 9.75 µg psPAX2, and 4.3 µg transgene for 10 cm dishes using Lipofectamine 3000.
- 7 Exchanged media 6 hours after transfection and collect and pool lentiviral supernatant at 24 hours and 48 hours after transfection.

## Transduction of primary human T cells

- 8 Centrifuged lentiviral supernatant at 600g for 10 min to remove cellular debris.
- 9 Concentrate lentivirus to 50–100× the initial concentration using Lenti-X Concentrator (Takara Bio).

- 10 Transduce T cells at 5–10% v/v of concentrated lentivirus at 24 h post-activation. For dual transduction experiments, T cells were serially transduced at 24 h and 48 h.

## scRNAseq

- 11 CD8<sup>+</sup>CCR7<sup>+</sup> T cells from three donors were transduced with CRISPRi and CRISPRa mini-TF gRNA libraries.
- 12 T cells were expanded for 10 days and then stained and sorted for Thy1.1<sup>+</sup> cells (a marker to identify transduced cells).
- 12.1 Antibody: Thy1.1-PE, clone OX-7, Mouse IgG1,  $\kappa$ , 1:300, Flow cytometry, StemCell Technologies, 60024PE
- 13 An SH800 FACS Cell Sorter (Sony Biotechnology) was used for cell sorting and analysis.
- 14 For antibody staining of Thy1.1 cells were collected, spun down at 300g for 5 min, resuspended in flow buffer (1× phosphate-buffered saline (PBS), 2 mM ethylenediaminetetraacetic acid and 0.5% bovine serum albumin) with the appropriate antibody dilutions and incubated for 30 min at 4 °C on a rocker.
- 15 Cells were then washed with flow buffer, spun down at 300g for 5 min and resuspended in flow buffer for cell sorting.
- 16 Sorted cells were loaded into the Chromium X for a targeted recovery of  $2 \times 10^4$  cells per donor and treatment according to the Single Cell 5'-High-Throughput Reagent Kit v2 protocol (10x Genomics).
- 17 SaCas9 gRNA sequences were captured by spiking in 2  $\mu$ M of a custom primer into the reverse transcription master mix, as previously done for SpCas9 gRNA capture<sup>36</sup>. The custom primer was designed to bind to the constant region of SaCas9's gRNA scaffold.
- 18 5'-Gene Expression (GEX) and gRNA libraries were separated using double-sided SPRI bead selection in the initial cDNA clean-up step. 5'-GEX libraries were constructed according to manufacturer's protocol.
- 19 gRNA libraries were constructed using two sequential PCRs (PCR 1: 10 cycles, PCR 2: 25 cycles). The PCR 1 product was purified using double-sided SPRI bead selection at 0.6 × and 2 ×. Twenty percent of the purified PCR 1 product was input into PCR 2. The PCR2 product was purified using double-sided SPRI bead selection at 0.6 × and 1 ×.

- 20 All libraries were run on a High Sensitivity D1000 tape to measure the average amplicon size and quantified using Qubit's dsDNA High Sensitivity assay.
- 21 Libraries were individually diluted to 20 nM, pooled together at desired ratios and sequenced on an Illumina NovaSeq S4 Full Flow Cell (200 cycles) with the following read allocation: Read 1, 26; i7 index, 10; Read 2, 90.
- 22 All oligos are available in Supplementary Table 5 of McCutcheon et al. Nature Genetics, 2023. <https://doi.org/10.1038/s41588-023-01554-0>

## Processing and analyzing scRNA-seq

- 23 CellRanger v6.0.1 was used to process, demultiplex and generate UMI counts for each transcript and gRNA per cell barcode.
- 24 UMI counts tables were extracted and used for subsequent analyses in R using the Seurat<sup>69</sup> v4.1.0 package.
- 25 Low-quality cells with <200 detected genes, >20% mitochondrial reads or <5% ribosomal reads were discarded.
- 26 DoubletFinder<sup>70</sup> was used to identify and remove predicted doublets. All remaining high-quality cells across donors for each treatment (CRISPRi or CRISPRa) were aggregated for further analyses.
- 27 gRNAs were assigned to cells if they met the threshold (gRNA UMI >4). Cells were then grouped on the basis of gRNA identity.
- 28 For differential gene expression analysis, we compared the transcriptomic profiles of cells sharing a gRNA to cells with only NT gRNAs using Seurat's FindMarkers function to test for DEGs with the hurdle model implemented in model-based analysis of single-cell transcriptomics (MAST).
- 29 All significant gRNA-gene links are available in Supplementary Table 3 of McCutcheon et al. Nature Genetics, 2023. <https://doi.org/10.1038/s41588-023-01554-0>
- 30 Upregulated DEGs were input into EnrichR's GO Biological Process 2021 database<sup>71</sup> for functional annotation.



## Protocol references

36. Mimitou, E. P. et al. Multiplexed detection of proteins, transcriptomes, clonotypes and CRISPR perturbations in single cells. *Nat. Methods* **16**, 409–412 (2019).
69. Butler, A., Hoffman, P., Smibert, P., Papalexi, E. & Satija, R. Integrating single-cell transcriptomic data across different conditions, technologies, and species. *Nat. Biotechnol.* **36**, 411–420 (2018).
70. McGinnis, C. S., Murrow, L. M. & Gartner, Z. J. DoubletFinder—doublet detection in single-cell RNA sequencing data using artificial nearest neighbors. *Cell Syst.* **8**, 329–337 (2019).
71. Kuleshov, M. V. et al. Enrichr—a comprehensive gene set enrichment analysis web server 2016 update. *Nucleic Acids Res.* **44**, 90–97 (2016).