

May 02, 2018

Version 2

Phenol/chloroform extraction of DNA from cyanobacteria V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.ptndnme>



Robbie Martin¹, Steven W. Wilhelm¹

¹University of Tennessee, Knoxville

The Aquatic Microbial E...



Robbie Martin

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ptndnme>

Protocol Citation: Robbie Martin, Steven W. Wilhelm 2018. Phenol/chloroform extraction of DNA from cyanobacteria.
protocols.io <https://dx.doi.org/10.17504/protocols.io.ptndnme>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 27, 2018

Last Modified: May 02, 2018

Protocol Integer ID: 11854

Keywords: DNA extraction, cyanobacteria, dna from microcystis aeruginosa, cyanobacteria phenol, chloroform extraction of dna, dna extraction from cell, other cyanobacteria, microcystis aeruginosa, extracting dna, dna extraction, based dna extraction, chloroform extraction, dna, chloroform, phenol, rnase, extraction

Abstract

Phenol/chloroform-based DNA extraction from cells pre-treated with RNase A, lysozyme, proteinase K, and SDS. The protocol was optimized for extracting DNA from *Microcystis aeruginosa*, but works well for other cyanobacteria.

Materials

MATERIALS

☒ Sodium acetate **P212121**

☒ RNase A **Qiagen Catalog #19101**

☒ Chloroform/IAA **Amresco Catalog #X205**

☒ TE buffer

☒ Phenol

☒ 70% Ethanol

☒ Proteinase K **Thermo Fisher Scientific Catalog #EO0491**

☒ Sodium Dodecyl Sulfate, 100gm **Promega Catalog #H5113**

☒ Lysozyme **Bio Basic Inc. Catalog #LDB0308.SIZE.1g**

☒ 100% Ethanol

Safety warnings

! See [MSDS](#) for safety and warnings.

Pre-treatment of cells with RNase A, lysozyme, proteinase K, and SDS

- 1 Concentrate ~4 mL of cellular culture in an Eppendorf tube by centrifugation at 10,000 xg for ~5 min.

 4 µL cellular culture

 00:05:00 centrifugation at 10,000 xg

- 2 Discard supernatant.

- 3 Resuspend pellet in 425 µL of standard TE buffer.

 425 µL TE buffer

Note

Standard TE buffer

-10 mM Tris
-1 mM EDTA
-pH 8

- 4 Add 50 µL of 100 µg/mL RNase A TER buffer (RNase A plus TE buffer) for a final concentration of 10 µg/mL RNase A.

 50 µL RNase A TER buffer

- 5 Add 25 µL of 100 mg/mL lysozyme for a final concentration of 5 mg/mL lysozyme.

 25 µL of 100 mg/mL lysozyme

- 6 Incubate at 37 °C for 20 minutes.

 37 °C incubation

 00:20:00 incubation

- 7 Add 50 µL of 1 mg/mL proteinase K for a final concentration of 100 µg/mL proteinase K.

- 8 Add 50 µL of 10% sodium dodecyl sulfate.

 50 µL 10% sodium dodecyl sulfate

- 9 Incubate at 50 °C for 2 hr.

 50 °C incubation



02:00:00 incubation

Phenol/chloroform extraction

10 Add 250 mL phenol (pH 8.0). Mix gently and completely.

250 µL phenol (pH 8.0)

11 Add 250 mL chloroform/isoamyl alcohol (24:1), mix gently and completely.

250 µL chloroform/isoamyl alcohol (24:1)

12 Spin at maximum speed in benchtop centrifuge for at least 2 minutes.

Note

Longer spin times, (up to 10 min.) may help in separating the phases.

13 Transfer aqueous (top) layer to clean Eppendorf tube without removing any of the organic layer.

14 Remove the last of the aqueous layer along with some of the organic layer (to ensure all aqueous volume is collected) and add it to a new, clean Eppendorf tube for later re-extraction.

Note

This re-extraction will significantly increase DNA recovery.

- 15
- Repeat steps 10, 11 and 12 twice (a total of 3 phenol extractions) or until no visible protein layer is seen.
 - Each time, remove the last of the aqueous layer along with some of the organic layer and add to the 're-extraction' Eppendorf tube (refer to step 14).

[go to step #10](#) At least 3 phenol extractions, or until no visible protein layer is seen.

Note

Protein layer will be white scum at interface.

[go to step #14](#) Add to 're-extraction' Eppendorf tube.

16 Add 500 µL of chloroform and mix gently.(1/2)



🧪 500 μ L chloroform

- 17 Spin at maximum speed in benchtop centrifuge for at least 2 minutes. (1/2)

Note

This will help remove traces of phenol.

🕒 00:02:00 centrifugation

- 18 Transfer aqueous (top) layer to clean Eppendorf tube without removing any of the organic layer. (1/2)

- 19 Repeat steps 16-18 to remove all phenol. (2/2)

➡ go to step #16 Repeat phenol removal

- 20 Extract the mixed aqueous/organic contents of the 're-extraction' Eppendorf tube and follow extraction steps 10-19. Combine the final aqueous layer with that collected above.

➡ go to step #10 extraction steps for aqueous/organic contents

Ethanol precipitation

- 21 Add 0.1 volume of 3 M sodium acetate to the collected aqueous phase and mix gently.

- 22 Add 2 volumes of ice cold 100% ethanol and gently mix well.

Note

Strings of precipitating DNA should become visible.

- 23 Place tube in -80 °C freezer until ethanol mixture is partially frozen (~1 hr.)

🧊 -80 °C freezer

🕒 01:00:00 approximate time for partial ethanol freeze

- 24 Spin at maximum speed in benchtop centrifuge at 4 °C for 30 min.

🕒 00:30:00 centrifugation

🧊 4 °C centrifugation

- 25 Discard supernatant and very carefully aspirate the remaining droplets of liquid from the tube without disrupting the DNA pellet.



- 26 Wash the DNA pellet by adding 500 μ L of ice-cold 70% ethanol to the tube and pipetting it gently, several times.

 500 μ L ice-cold 70% ethanol

- 27 Spin at maximum speed in a benchtop centrifuge at 4 °C for 15 min.

 00:15:00 centrifugation

 4 °C centrifugation

- 28 Discard supernatant and carefully aspirate the remaining droplets of liquid from the tube without disrupting the DNA pellet.

- 29 Place in a heater block at ~37 °C for less than 5 minutes to evaporate all ethanol.

Note

Be careful! DO NOT dry to completion! Alternatively, you can air dry to remove the ethanol.

 00:04:59 ethanol drying

 37 °C heater block

- 30 Re-suspend the DNA pellet in 50-100 μ L of TE buffer or water (depending on downstream use and needed concentration).

 50 μ L TE buffer or water

- 31 Freeze at -20 or -80 °C.