

Jul 02, 2020

Phenazine Oxidizer Enrichment and Isolation

DOI

<https://dx.doi.org/10.17504/protocols.io.bh4tj8wn>

Lev Tsypin¹, Yinon Bar-On², Scott Saunders¹, Jared R Leadbetter¹, Dianne K Newman¹

¹California Institute of Technology; ²Weizmann Institute of Science



Lev Tsypin

Stanford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bh4tj8wn>

Protocol Citation: Lev Tsypin, Yinon Bar-On, Scott Saunders, Jared R Leadbetter, Dianne K Newman 2020. Phenazine Oxidizer Enrichment and Isolation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bh4tj8wn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

This protocol was used during the 2017 Microbial Diversity course at the Marine Biological Laboratory in Woods Hole, MA.

Created: July 01, 2020



Last Modified: July 02, 2020

Protocol Integer ID: 38771

Keywords: phenazine oxidizer enrichment, oxidizing microbe, soil sample, protocol for the enrichment

Abstract

This is a protocol for the enrichment of phenazine-1-carboxylic acid oxidizing microbes from a soil sample.

Medium composition

1 Stock solutions

100x Freshwater salts:

- 1.71 M sodium chloride
- 197 mM magnesium chloride
- 68 mM calcium chloride
- 671 mM potassium chloride

1000x Trace Elements:

- 20 mM HCl
- 7.5 mM FeCl₃ 6H₂O
- 480 uM H₃BO₃
- 500 uM MnCl₂ 4H₂O
- 6.8 mM CoCl₂ 6H₂O
- 1 mM NiCl₂ 6H₂O
- 12 uM CuCl₂ 2H₂O
- 500 uM ZnCl₂
- 23 uM Na₂SeO₃
- 150 uM Na₂MoO₄

1000× 13-vitamin solution

- 1 mM MOPS pH 7.2
- 100 µg/mL riboflavin
- 30 µg/mL biotin
- 100 µg/mL thiamine HCl
- 100 µg/mL L-ascorbic acid
- 100 µg/mL d-Ca-pantothenate
- 100 µg/mL folic acid
- 100 µg/mL nicotinic acid
- 100 µg/mL 4-aminobenzoic acid
- 100 µg/mL pyridoxine HCl
- 100 µg/mL lipoic acid
- 100 µg/mL Nicotinamide adenine dinucleotide
- 100 µg/mL thiamine pyrophosphate
- 10 µg/mL cyanocobalamin

2 mM reduced phenazine-1-carboxylic acid (PCA)

- In a sealed Balch tube or serum vial, sparge 2 mM oxidized PCA with 80:20 (vol:vol) H₂/CO₂ over Pd (II) until yellow-green in color
- Stir with a magnetic stirbar or otherwise agitate during sparging and overnight afterwards.

Enrichment medium (anoxic)

- 45 mM sodium bicarbonate buffer (buffer)
- 10 mM ammonium chloride (nitrogen source)
- 1 mM potassium phosphate (phosphorus source)
- 1 mM reduced phenazine-1-carboxylic acid (PCA, electron donor)
- 100 µM sodium acetate (carbon source)
- 50 µM sodium sulfide (sulfur source)
- 1x trace elements solution
- 1x freshwater salts solution
- 1× 13-vitamin solution

Isolation medium (oxic or anoxic, agar plates)

- 1 mM sodium phosphate buffer pH 7
- 10 mM ammonium chloride
- 100 µM sodium sulfate
- 1 mM sodium acetate
- 10 mM sodium nitrate
- 1x freshwater salts solution
- 1x trace elements solution
- 1× 13-vitamin solution
- 1.5 % agar

Sampling

- 2 Collect topsoil by scooping directly with 50 mL conical tube.
- 3 Grind with mortar and pestle
- 4 Pass through 2 mm sieve to remove large rocks and organic matter.
- 5 Transfer to anaerobic chamber



Inoculation of enrichment cultures

- 6 Prepare 96-well plate with 190 μ L of anoxic enrichment medium per well.
- 7 Add 20 mg soil into each well.
- 8 Prepare Balch tubes with 9 mL of anoxic enrichment medium in each.
- 9 Add 1 g soil into each tube and seal.
- 10 For both the 96-well plates and Balch tubes, watch for a color change over the next days. The yellow-green reduced PCA will turn clear if oxidized.
- 11 Passage candidate PCA-oxidizing enrichments via 1:10 dilutions into fresh medium, maintaining anoxic conditions.

Isolation

- 12 Plate enrichment cultures onto LB/Agar and the solid isolation medium.
- 13 Grow both oxically and anoxically in an anaerobic chamber.
- 14 Pick colonies of different morphologies and assay PCA-oxidizing capability in the enrichment medium as a pure culture.
- 15 Select isolates that oxidize PCA in pure culture for further study.