

Feb 09, 2024

# 🌐 Plaque Assay for *Microcystis aeruginosa* NIES-298 and phage Ma-LMM01

DOI


<https://dx.doi.org/10.17504/protocols.io.e6nvwd1j7lmk/v1>

Laura Smith<sup>1</sup>, Steven W Wilhelm<sup>1</sup>

<sup>1</sup>The University of Tennessee, Knoxville



Laura E Smith



## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.e6nvwd1j7lmk/v1>

**Protocol Citation:** Laura Smith, Steven W Wilhelm 2024. Plaque Assay for *Microcystis aeruginosa* NIES-298 and phage Ma-LMM01. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwd1j7lmk/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working



**Created:** August 29, 2023

**Last Modified:** February 09, 2024

**Protocol Integer ID:** 87114

**Keywords:** plaque assay, cyanophage, *Microcystis aeruginosa*, harmful algal blooms, reproducibly quantifies infection by the microcystis, infecting phage ma, microcystis, plaque assay for the enumeration, phage ma, phage, reproducibly quantifies infection, assay

**Funders Acknowledgements:**

The Simons Foundation

Grant ID: 735077

## Abstract

This protocol describes a plaque assay for the enumeration of phage infecting *Microcystis aeruginosa* NIES-298. It is a single agarose layer assay and reproducibly quantifies infection by the *Microcystis*-infecting phage Ma-LMM01. This protocol has been adapted from work by Moore et al. (2007) as per Lindell (2014).

## Guidelines

The protocol has been designed to work with cultures / lysates in a laboratory setting. If phage that will infect this host are present in natural samples, there is no reason the protocol should not detect them.

## Materials

*Microcystis aeruginosa* NIES 298 is available from the National Institute for Environmental Studies (Japan) Ma-LMM01 was kindly provided by Professor Professor Takashi Yoshida (Kyoto University).



## Prepare CT Agarose

- 1 Make 0.56% agarose in **CT medium** base by adding 0.56 g low melting point agarose / 100 mL CT medium prior to autoclaving. You will need 12 mL of CT agarose per plate. Note: the final agarose concentration in plates will be ~ 0.45% once culture and lysate are added (see below). CT medium is modified from Watanabe and Hiroki, 1997 by adding phosphorus of equal concentrations using sodium phosphate in place of  $\text{Na}_2\text{*beta-glycerophosphate}$ .
- 2 After autoclaving, place in water bath at  $34\text{ }^{\circ}\text{C}$  . Add vitamins once agarose has cooled but not solidified. Make 12 mL aliquots and keep each in a water bath at  $34\text{ }^{\circ}\text{C}$  until ready to plate. Alternatively, CT agarose (without vitamins) can be prepared in advance and microwaved prior to the water bath step.

## Concentrate *Microcystis aeruginosa* cells

- 3 Concentrate log-phase *M. aeruginosa* NIES-298 cells (culture should be no more than 3-4 days old and  $\sim 1\text{-}2 \times 10^6$  cells/mL) by centrifugation at 8000 xg for 10 min (x2). The first centrifugation step will collapse the gas vacuoles; the second will form the pellet. Remove excess medium.
- 4 Resuspend cells in fresh CT medium. Cells should be concentrated to  $\sim 2\text{-}3 \times 10^7$  cells/mL. You will need 2 mL of concentrated cells per plate.

## Dilute Viral Lysate

- 5 Make a serial dilution of virus containing solution to be tested using CT media. Fresh Ma-LMM01 lysates generally contains  $\sim 2\text{-}6 \times 10^7$  pfu/mL, so  $10^{-5}$  - $10^{-7}$  dilutions will get plates in the countable range of  $\sim 20\text{-}200$  pfu (plaque forming units) per plate. PFUs outside this range are difficult to quantify accurately.

## Plate

- 6 Add 1 mL of diluted sample to 2 mL of concentrated cells. For negative controls, add 1 mL of CT medium to 2 mL concentrated cells.
- 7 Pour the 12 mL aliquot of CT low melt agarose into the plate. Add the 3 mL combination of lysate and cells next to the agarose. Swirl gently to mix together to form one homogenous layer.



## Incubate

- 8 Plates can be carefully moved to an incubator adjusted for *Microcystis* growth immediately after pouring, or once they are solidified. Note: plates will take 1-2 hours to solidify.  
  
Important: unlike standard bacterial plaque assays **do not flip plates**. Plates will not be firm enough for flipping for at least one day.
- 9 Incubate at  26 °C and 30-40  $\mu\text{mol photons m}^{-2} \text{s}^{-1}$  continuous light. Plates can be flipped on day 2. Plaques should start to become visible ~day 3, and can be counted through day 5.

## Protocol references

Lindell, D. (2014). The Genus *Prochlorococcus*, Phylum Cyanobacteria. In: Rosenberg, E., DeLong, E.F., Lory, S., Stackebrandt, E., Thompson, F. (eds) The Prokaryotes. Springer, Berlin, Heidelberg. [https://doi.org/10.1007/978-3-642-38954-2\\_125](https://doi.org/10.1007/978-3-642-38954-2_125)

Moore, L.R., Coe, A., Zinser, E.R., Saito, M.A., Sullivan, M.B., Lindell, D., Frois-Moniz, K., Waterbury, J. and Chisholm, S.W., 2007. Culturing the marine cyanobacterium *Prochlorococcus*. *Limnology and Oceanography: Methods*, 5(10), pp.353-362.

Watanabe, M. M. and Hiroki, M. 1997. NIES-Collection List of Strains, Microalgae and Protozoa, 5th ed. Microbial Culture Collection, National Institute for Environmental Studies, Tsukuba, Japan. Available at <https://www.nies.go.jp/kanko/gyomu/pdf/f097-1997.pdf> (last access Feb 7, 2024)