


Jun 17, 2017

Plaque Assay For Screening Viral Concentrates for Bacteriophage

DOI

dx.doi.org/10.17504/protocols.io.igxcbxn



Dr. Steven Wilhelm

The Aquatic Microbial E...



Steven W Wilhelm

The University of Tennessee, Knoxville

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.igxcbxn

Protocol Citation: Dr. Steven Wilhelm 2017. Plaque Assay For Screening Viral Concentrates for Bacteriophage. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.igxcbxn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 15, 2017

Last Modified: March 24, 2018

Protocol Integer ID: 6391

Abstract


Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Culture Preparation


- 1 Grow a 4 mL liquid culture of the bacterial isolate of interest in an appropriate growth medium overnight

 4 mL

- 2 Dilute the culture 10-20% and let it grow for ~4 hrs prior to plaque assay


 04:00:00

- 3 Melt 0.6% top-agar in the microwave and let cool ~10-15 min

 00:15:00

- 4 Pipet 1.5 mL liquid culture into 1.5 mL Eppendorf tubes (or enough liquid culture for 500 μ L for each plate sample that you want to make).

- 5 Centrifuge 2 min at 14,000 rpm


 00:02:00

Note

Some bacterial isolates do not require centrifugation. This is something that needs to be empirically determined for each putative bacterial host.

- 6 Discard the supernatant

- 7 Resuspend the remaining cells in 1.5 mL BG-11 media and mix well with pipet

 1.5 mL


- 8 Pipet 4 mL of molton top agar into sterile culture tubes

 4 mL

- 9 Keep top agar in a 45°C head block to keep in a liquid state


- 10 In 500 μ L of resuspended cells, add 100 μ L of 0.22 μ m filtered viral concentrate

- 11 Incubate for 15 min at an appropriate growth temperature for the bacterium in use

 00:15:00



Plaque Assay

- 12 Pipet infected cells into molton top agar
- 13 Vortex briefly avoiding bubbles
- 14 Pour top agar onto the plate and spread
- 15 Let agar solidify at room temperature ~30-60 min
 01:00:00
- 16 Flip plates over and place in an incubator at an appropriate growth temperature and light conditions for the bacterium in use
- 17 Check plates daily until a confluent lawn appears on plates. Plaques will be obvious clearings.