


Mar 20, 2023

Plate Count Agar

DOI

dx.doi.org/10.17504/protocols.io.eq2ly74delx9/v1



Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly74delx9/v1>

Protocol Citation: Andreas Sagen 2023. Plate Count Agar. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2ly74delx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 09, 2023

Last Modified: March 20, 2023

Protocol Integer ID: 78371

Keywords: Plate count agar, agar, plate counting, counting, CFU, plate count agar plate count agar, microbiological growth medium, pca, bacteria, substrate for bacteria, using pca, aerobic bacteria, living aerobic bacteria, viable bacterial growth, growth of the organism, units per gram, solid sample, liquid sample, fermentable carbohydrate, yeast extract, organism, sample, glucose, gram, substrate, source for vitamin, plate, nitrogen

Abstract

Plate Count Agar (PCA), is a microbiological growth medium to assess viable bacterial growth of a sample. The total number of living aerobic bacteria can be determined using PCA which is a substrate for bacteria to grow on. The medium contains **casein** which provides nitrogen, carbon, amino acids, vitamins and minerals to aid in the growth of the organism. **Yeast extract** is the source for vitamins and **glucose** is the fermentable carbohydrate. **Agar** is the solidifying agent. This is a non-selective medium and the bacteria is counted as colony forming units per gram (CFU/g) in solid samples and (CFU/ml) in liquid samples.

Guidelines

Prepare enough for the necessary number of experiments to perform in order to prevent batch-to-batch variation between experiments. Input chemicals should be of a certain quality in order to prevent contaminants that can impact microbial growth and selection.

Materials

Magnetic stirrer, autoclave, scale, tubes, flask

Protocol materials

⊗ Tryptone Merck Millipore (EMD Millipore) Catalog #T9410

⊗ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

⊗ Glucose Merck MilliporeSigma (Sigma-Aldrich) Catalog #G7021

⊗ Agar Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296

Troubleshooting



Safety warnings

- ⚠ Be sure to wear appropriate PPE when working with antibiotics and chemicals. Take care when working with hot flasks and tubes.



500 mL Plate Count Agar

2h

1 In a bottle, add approx.  400 mL deionized water

5m

2 Measure and add:

15m


 5 g Tryptone

 2.5 g Yeast extract


 1 g Glucose

 15 g Agar


Materials:

 Tryptone Merck Millipore (EMD Millipore) Catalog #T9410


 Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

 Glucose Merck MilliporeSigma (Sigma-Aldrich) Catalog #G7021

 Agar Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296

3 Adjust pH to  7.0 using sodium hydroxide

5m

4 Fill bottle to  500 mL total with deionized water

5m

5 Autoclave liquid at  121 °C for  00:45:00

45m

6 Aliquot liquid in  15 mL and  50 mL units and store refrigerated at  4 °C

15m