

Mar 25, 2019

Version 2

Polymerase Chain Reaction - Comparative genomics of Staphylococcus aureus associated with subclinical and clinical bovine mastitis (Rocha et al., 2019) V.2

DOI

dx.doi.org/10.17504/protocols.io.zhpf35n

Lis LSR Rocha<sup>1</sup>

<sup>1</sup>Universidade Federal de Viçosa



Lis LSR Rocha

Universidade Federal de Viçosa

#### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.zhpf35n

**Protocol Citation:** Lis LSR Rocha 2019. Polymerase Chain Reaction - Comparative genomics of Staphylococcus aureus associated with subclinical and clinical bovine mastitis (Rocha et al., 2019). **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.zhpf35n

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 25, 2019

Last Modified: March 25, 2019

Protocol Integer ID: 21775

Keywords: comparative genomics of staphylococcus aureus, pathogenesis of bovine mastitis, bovine mastitis, clinical bovine mastitis, staphylococcus aureus, causing clinical mastitis, clinical mastitis, mastitis outcome, genes present in the subclinical strain, causing subclinical infection, subclinical infection, clinical isolates from the subclinical one, virulence factor, comparative genomic, clinical isolate, subclinical strain, absence of virulence factor, bacteria, subclinical mastit, strain, virulence, comparative genomic analysis, several nonsynonymous single nucleotide polymorphism, associated gene, characterized strain, pathogenesi, nucleotide sequence, polymerase chain reaction

#### Abstract

Comparative genomics of Staphylococcus aureus associated with subclinical and clinical bovine mastitis (Rocha et al., 2019)

Many efforts have been made to understand the pathogenesis of bovine mastitis to reduce losses and promote animal welfare. Staphylococcus aureus may cause bovine clinical mastitis, but it is mainly associated with subclinical infection, which is usually persistent and can easily reoccur. Here, we conducted a comparative genomic analysis between four strains of S. aureus causing subclinical infection (Sau170, 302, 1269, 1364), previously sequenced by our group, and two well-characterized strains causing clinical mastitis (N305 and RF122) to find differences that could be linked to mastitis outcome. A total of 146 virulence-associated genes were compared and no appreciable differences were found between the bacteria. However, several nonsynonymous single nucleotide polymorphisms (SNPs) were identified in genes present in the subclinical strains when compared to RF122, especially in genes encoding host immune evasion and surface proteins. The comparison of orthologous genes using OrthoMCL identified a membrane transporter in the genomes of the bacteria belonging to the subclinical group, but this finding was not confirmed by polymerase chain reaction (PCR) on a collection of field isolates of S. aureus associated with clinical or subclinical mastits. The secreted and surface proteins predicted by different in silico tools were compared through multidimensional scaling analysis, revealing a high degree of similarity among the six strains. However, differences were seen in the nucleotide sequences of a gene that codes for a hypothetical protein (cl3309) and a lipoprotein (cl3700). These findings were also analyzed by PCR on DNA extracted from field isolates of S. aureus. The lipoprotein, but not the hypothetical protein, was able to separate the clinical isolates from the subclinical ones. These results show that sequence variation among bovine S. aureus, and not only the presence/absence of virulence factors, is an important aspect to consider when comparing isolates causing different mastitis outcomes



#### **Materials**

50 ng of total DNA, 1U of Taq DNA polymerase Cellco Biotec, 0.2 μM of each primer, 0.2 mM deoxynucleotide triphosphate mixture, 1X reaction buffer containing 2.0 mM MgCl2, extra 1.0 mM MgCl2, and Milli-Q water to increase the reaction volume to a final volume of 25  $\mu$ L.

The extra 1 mM MgCl<sub>2</sub> was excluded from the PCR reactions that contained the primers LipoP-F-CS/LipoP-R-C.

Table 1 - Primer Sequences for primers used in this Protocol

cl3309s ubF	TGTTGTAGGAGGAACAAT CC
cl3309s ubR	TTCTAATGTCAGCAACATG C
cl3309c liF	GCTATTCCTAGATGCACT
cl3309c liR	TTTTAAGTATGACATGAAT G
cl3316F	ACGCAAAACCCTTTACTA GT
cl3316 R	GCAACAACTAGTAGGAGT GA
LipoP- F-CS	GYTTTGCGAAAACGTTAG AYATGTA
LipoP- R-C	TGCCTTCATCATTAATTGG ACCAATC
LipoP- F-CS	GYTTTGCGAAAACGTTAG AYATGTA
LipoP- R-CS	GGTAAAYTCAATGTYCTTA TRTCC



# Troubleshooting



# primers cl3309sub F/R

- Initial denaturation: 95.0 °C for 5 min;
- 2 35 cycles of denaturation at 95.0 °C for 45 s,
- 3 Anealing: 55 °C for 45s
- 4 Extension: 72 °C for 45 s
- 5 final extension at 72.0 °C for 10 min

# primers cl3316F/R

- 6 initial denaturation: 95.0 °C for 5 min;
- 7 35 cycles of denaturation at 95.0 °C for 45 s,
- 8 Annealing: 55 °C for 45 s
- 9 Extension: 72 °C for 45 s
- 10 final extension at 72.0 °C for 10 min.

# primers cl3700 - LipoP FCS/RC

11 initial denaturation: 95.0 °C for 5 min;



- 12 35 cycles of denaturation at 95.0 °C for 45 s,
- 13 Annealing: 54 °C for 45 s
- 14 Extension: 72 °C for 45 s
- 15 final extension at 72.0 °C for 10 min.

#### cl33009cli F/R

- 16 initial denaturation: 95.0 °C for 5 min;
- 17 35 cycles of denaturation at 95.0 °C for 45 s,
- 18 Annealing: 45 °C for 45 s
- 19 Extension: 72 °C for 30 s
- 20 final extension at 72.0 °C for 10 min.

### primers cl3700 - LipoP FCS/RCS

- 21 initial denaturation: 95.0 °C for 5 min;
- 22 35 cycles of denaturation at 95.0 °C for 45 s,



- 23 Annealing: 50 °C for 45 s
- 24 Extension: 72 °C for 1min
- 25 final extension at 72.0 °C for 10 min.

# Analyzing the amplified fragments

Analyze the amplicons by electrophoresis in 1X Tris-acetate-EDTA on a 1.0% agarose gel 26 and visualize imagen under UV light after staining with 2 mg.ml-1 ethidium bromide.