


May 20, 2022

🌐 Preparation and use of 12-well plates for the rapid detection of terbinafine-resistant dermatophytes

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6116zvqe/v1>



Khalid El Moussaoui¹

¹AIDS Reference Laboratory, Department of Clinical Microbiology, University Hospital of Liege, 4000 Liege, Belgium



Khalid El Moussaoui

AIDS Reference Laboratory, Department of Clinical Microbiolo...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6116zvqe/v1>

Protocol Citation: Khalid El Moussaoui 2022. Preparation and use of 12-well plates for the rapid detection of terbinafine-resistant dermatophytes . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6116zvqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: February 11, 2022

Last Modified: May 20, 2022

Protocol Integer ID: 58072

Keywords: terbinafine resistant dermatophyte, resistant dermatophyte, well plates for the rapid detection, screening plate, plates for the rapid detection, terbinafine, rapid detection, well plate

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol describes the steps necessary for the preparation and use of screening plates for the rapid detection of terbinafine resistant dermatophytes.

Guidelines

- 1) When you take the medium out of the autoclave, be as quick as possible to add the terbinafine and transfer the medium to the plate wells. Since the medium contains agar, it will solidify fairly quickly.
- 2) Work in sterile conditions as much as possible. The medium is very easy to contaminate and can, if so, falsify the results.



Protocol materials

✕ Sabouraud dextrose broth **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S3306**

✕ MilliQ water

✕ Difco granulated agar **Fisher Scientific Catalog #BD 214530**

✕ Terbinafine hydrochloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8826**

✕ Dimethyl sulfoxide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4540**

Terbinafine preparation

13m

1 Weigh  16.33 mg of

13m

 Terbinafine hydrochloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8826**

and add to an erlenmeyer flask containing  5 mL of

 Dimethyl sulfoxide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4540** . Mix

without heating for  00:10:00 to dissolve. The final concentration of this terbinafine preparation is  3200 µg/mL . Do not forget to add a magnetic bar. Cover each flask with glass wool and aluminium foil.

Equipment

Heated magnetic stirrer

NAME

AG602

TYPE

LMR

BRAND

AG602

SKU



2 Autoclave this flask at  121 °C for  00:30:00

30m



Equipment

Autoclave

NAME

Tuttnauer

BRAND

3850EL

SKU

<https://tuttnauer.com/laboratory-autoclaves/laboratory-research/benchtop-autoclaves/el-d-line>

LINK



- 3 Under hood, transfer  1 mL of this preparation into 2 mL eppendorf tubes and store at  -80 °C for later use.

5m

Equipment

Biological safety cabinet

NAME

AC2-4S8-TU

TYPE

Esco

BRAND

303110

SKU

<https://analys.be/fr-be/sample-handling-in-ivf-labs/ac2-4s8-tu-airstream-plus-class-ii-biological-303110>

LINK



Medium preparation

4h 25m



4 Dissolve  30 g of 8m

 Sabouraud dextrose broth **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S3306**

in  1 L of  MilliQ water and let mix on the heated magnetic stirrer for

 00:05:00 (temperature and mixing speed knob at mid-step)

Equipment

Heated magnetic stirrer

NAME

AG602

TYPE

LMR

BRAND

AG602

SKU



5 Add  15 g of  Difco granulated agar **Fisher Scientific Catalog #BD 214530** and 12m
let dissolve for  00:10:00 while mixing and heating

6 Pour  250 mL of this medium into 4 erlenmeyer flasks and label each with final 5m
terbinafine concentrations (0.2 µg/mL for the first flask, 0.1 µg/mL for the second, 0.05 µg/mL for the third and 0 µg/mL for the last flask). Given the uncertainty of measurement, it is advisable to work from the highest concentration to the lowest. Do not forget to add a magnetic bar. Cover each flask with glass wool and aluminium foil.

7 Autoclave these 4 flasks at  121 °C for  00:30:00 4h



Equipment

Autoclave

NAME

Tuttnauer

BRAND

3850EL

SKU

<https://tuttnauer.com/laboratory-autoclaves/laboratory-research/benchtop-autoclaves/el-d-line>

LINK



Plates preparation

1h 11m

- 8 For the first flask (0.2 $\mu\text{g/mL}$) : add  15.6 μL of a 3200 $\mu\text{g/mL}$ terbinafine preparation (see section number 1). Let mix while heating (temperature and mixing speed knob at mid-step) for  00:02:00 . Caution : do not leave for more than 2 minutes otherwise the medium will solidify. 2m
- 9 For the second flask (0.1 $\mu\text{g/mL}$) : add  7.8 μL of a 3200 $\mu\text{g/mL}$ terbinafine preparation (see section number 1). Let mix while heating (temperature and mixing speed knob at mid-step) for  00:02:00 . Caution : do not leave for more than 2 minutes otherwise the medium will solidify. 2m
- 10 For the third flask (0.05 $\mu\text{g/mL}$) : add  3.9 μL of a 3200 $\mu\text{g/mL}$ terbinafine preparation (see section number 1). Let mix while heating (temperature and mixing speed knob at mid-step) for  00:02:00 . Caution : do not leave for more than 2 minutes otherwise the medium will solidify. 2m
- 11 For the last flask (0 $\mu\text{g/mL}$) : does not contain terbinafine.
- 12 In a 12-well plate (3 rows and 4 columns), the first well of each row does not contain terbinafine. It serves as a growth control. Using an automated pipettor, add  3 mL of 15m

the medium that does not contains terbinafine (flask number 4) to all the first wells of the plate.

Equipment

Motorized Pipette Filler

NAME

Levo Plus

TYPE

DR. LAB

BRAND

7033100100

SKU

http://www.drlab.co/index.php?route=product/product&product_id=177^{LINK}



Equipment

Microwell plate

NAME

Thermo Scientific

BRAND

130185

SKU

<https://www.fishersci.com/shop/products/biolite-multidishes-microwell-plates/12556005>

LINK



- 13 The second wells in each row contain a concentration of 0.05 µg/mL of terbinafine. Using an automated pipettor, add  3 mL of the medium that contains  0.05 µg/mL of terbinafine (flask number 3) to all the second wells of the plate.

15m

Equipment

Motorized Pipette Filler

NAME

Levo Plus

TYPE

DR. LAB

BRAND

7033100100

SKU

http://www.drlab.co/index.php?route=product/product&product_id=177^{LINK}



Equipment

Microwell plate

NAME

Thermo Scientific

BRAND

130185

SKU

<https://www.fishersci.com/shop/products/biolite-multidishes-microwell-plates/12556005>

LINK



- 14 The third wells in each row contain a concentration of 0.1 µg/mL of terbinafine. Using an automated pipettor, add  3 mL of the medium that contains  0.1 µg/mL of terbinafine (flask number 2) to all the third wells of the plate.

15m

Equipment

Microwell plate

NAME

Thermo Scientific

BRAND

130185

SKU

<https://www.fishersci.com/shop/products/biolite-multidishes-microwell-plates/12556005>

LINK



Equipment

Motorized Pipette Filler

NAME

Levo Plus

TYPE

DR. LAB

BRAND

7033100100

SKU

http://www.drlab.co/index.php?route=product/product&product_id=177^{LINK}



- 15 The last wells in each row contain a concentration of 0.2 µg/mL of terbinafine. Using an automated pipettor, add  3 mL of the medium that contains  0.2 µg/mL of terbinafine (flask number 1) to all the last wells of the plate.

15m



Equipment

Motorized Pipette Filler

NAME

Levo Plus

TYPE

DR. LAB

BRAND

7033100100

SKU

http://www.drlab.co/index.php?route=product/product&product_id=177^{LINK}

Equipment

Microwell plate

NAME

Thermo Scientific

BRAND

130185

SKU

<https://www.fishersci.com/shop/products/biolite-multidishes-microwell-plates/12556005>

LINK



- 16 Cover the plates with their lids and let the medium solidify at room temperature. Label the plates with the corresponding terbinafine concentrations for each well and store at

 4 °C for future use.

5m

Use of the plates



- 17 Using a sterile swab, prepare a suspension (sterile water) 0.5 McF of the dermatophyte to be studied.
- 18 Transfer  25 μ L of this suspension to each of the 4 wells of the plate.
- 19 Annotate the plate with the number of the studied strains and incubate at 30°C for 96h before reading the result.