

Jul 22, 2024

Preparation of viral sequencing library for Illumina using NEBNext ultra II

 Forked from [Preparation of viral sequencing library for Illumina using WTA2 and QIAseq FX](#)

 [Japanese Journal of Infectious Diseases](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3xrw1g25/v1>

Kenichi Komabayashi¹

¹Yamagata prefectural institute of public health



Kenichi Komabayashi

Yamagata prefectural institute of public health

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3xrw1g25/v1>

External link: <https://doi.org/10.7883/yoken.jjid.2025.055>

Protocol Citation: Kenichi Komabayashi 2024. Preparation of viral sequencing library for Illumina using NEBNext ultra II. protocols.io <https://dx.doi.org/10.17504/protocols.io.kqdg3xrw1g25/v1>

**Manuscript citation:**

Kenichi Komabayashi, Shuji Chikaoka, Hiroki Awano, Akiko Abe, Dai Sendo, Masayuki Furuyama, Masatoshi Kitamura, Saori Tanabe, Manabu Wada, Naomi Ogawa, Mika Sasaki, Yoko Aoki, Yohei Matoba, Junji Seto, Tatsuya Ikeda, Katsumi Mizuta (2025) Parechovirus A3 infections, including myalgia, in Yamagata, Japan between 2022 and 2023 and their next-generation sequencing analyses. Japanese Journal of Infectious Diseases doi: [10.7883/yoken.jjid.2025.055](https://doi.org/10.7883/yoken.jjid.2025.055)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 12, 2024

Last Modified: July 22, 2024

Protocol Integer ID: 96568

Keywords: Genome Sequencing, illumina, metagenome, nuclease, iSeq100, virus, genome sequence of rna virus, analysis on illumina sequencer, rna virus, illumina sequencer, abundant sequences of viral origin, nucleic acids outside the viral particle, library for illumina, template rna, sequencing library, viral particle, metagenomic approach, rna, viral origin, using nuclease, genome sequence, illumina, nucleic acid, genome, sequences of host, sequencing, stranded dna, iseq100

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This method uses a metagenomic approach to analyze the genome sequence of RNA viruses. Nucleic acids outside the viral particles are reduced using nucleases and extracted to obtain template RNA. Templates are converted to double-stranded DNA, and library preparation is performed for analysis on Illumina sequencers. Analysis data with reduced sequences of host and bacterial origin and abundant sequences of viral origin are obtained, allowing multiple samples to be analyzed even with the throughput of the iSeq100.

This protocol was forked from "Preparation of viral sequencing library for Illumina using WTA2 and QIAseq FX".



Guidelines

There are three advantages to using this method.

- (1) No need to design virus-specific primers
- (2) RNA viral genomes can be analyzed while reducing nucleotides of host and bacterial origin
- (3) 10 or more samples can be analyzed at a time on the iSeq100 (For viral genome less than 10,000 bases)

The method consists of two parts: pretreatment and library preparation.

The pretreatment is intended to increase the content of virus-derived nucleic acids in the sample and facilitate genome analysis. The main point of this method is to reduce host genome, ribosomal RNA, and nucleic acids derived from bacteria in advance, taking advantage of the fact that genomes in viral particles are not easily digested by Nuclease.

NEBNext Ultra II RNA Library Prep Kit for Illumina (E7770) is used for library preparation.

The following points are different from the method described in the attached manual.

- (1) Use half the scale of the protocol in the manual.
- (2) The number of cycles of PCR amplification is increased because the initial RNA input is small.

The amplified libraries are pooled into a tube for multiplex analysis. The size and molar concentration of the pooled libraries are adjusted to obtain a final library that can be applied to the flow cell.

Materials

<Pretreatment >

Equipment

New Steradisc	NAME
0.45µm filter 50pcs	TYPE
Kurabo	BRAND
S-1304	SKU

⊗ Micrococcal Nuclease - 320,000 gel units **New England Biolabs Catalog #M0247S**

⊗ Benzonase® Nuclease 2.5ku **Catalog #70746-4CN**

⊗ High Pure Viral RNA Kit **Roche Catalog #11858882001**

Recipe for 100mL of homemade buffer (1M Tris, 100mM CaCl₂, 30mM MgCl₂, pH8)

1. Dissolve 15.06g of Trizma preset crystal pH7.5(M.W. 150.6) into 70mL of nuclease-free distilled water
2. Adjust to pH 8.0 by adding 4.92mL of NaOH (5N)
pH is measured after the temperature drops to room temperature
3. Add 1.47 g of CaCl₂-2H₂O (M.W.* 147.01) and 0.813 g of MgCl₂-6H₂O (M.W. 203.30)
4. Dissolve, and meth up to 100mL
5. Filtrate through a 0.22-µm filter, dispense into tubes, and store.

*: molecular weight

(optional)

⊗ Qubit RNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32852**

<Library preparation>

⊗ NEBNext Ultra II RNA Library Prep Kit for Illumina - 24 rxns **New England Biolabs Catalog #E7770S**

⊗ Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

⊗ NEBNext Multiplex Oligos for Illumina (Index Primers Set 1) - 24 rxns **New England Biolabs Catalog #E7335S**

⊗ NEBNext Multiplex Oligos for Illumina (Index Primers Set 2) - 24 rxns **New England Biolabs Catalog #E7500S**



⊗ Qubit 1X dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q33230**

⊗ NextSeq PhiX Control Kit **Illumina, Inc. Catalog #FC-110-3002**

Protocol materials

⊗ Micrococcal Nuclease - 320,000 gel units **New England Biolabs Catalog #M0247S**

⊗ Benzonase® Nuclease 2.5ku **Catalog #70746-4CN**

⊗ High Pure Viral RNA Kit **Roche Catalog #11858882001**

⊗ Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

⊗ Qubit 1X dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q33230**

⊗ NextSeq PhiX Control Kit **Illumina, Inc. Catalog #FC-110-3002**

⊗ Qubit RNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32852**

⊗ NEBNext Ultra II RNA Library Prep Kit for Illumina - 24 rxns **New England Biolabs Catalog #E7770S**

⊗ NEBNext Multiplex Oligos for Illumina (Index Primers Set 1) - 24 rxns **New England Biolabs Catalog #E7335S**

⊗ NEBNext Multiplex Oligos for Illumina (Index Primers Set 2) - 24 rxns **New England Biolabs Catalog #E7500S**



Reduction of nucleic acids derived from non-virus (pretreatment)

- 1 Collect  400 μL or more virus culture medium in a 1.5 mL tube.

Note

If the viral particles are broken, the genome could be digested in this later process.

- 2 Centrifuge  00:03:00 at 17,000 x g and aspirate the supernatant with a 1 mL tuberculin syringe.

3m

3

Equipment

New Steradisc

NAME

0.45 μm filter 50pcs

TYPE

Kurabo

BRAND

S-1304

SKU

Filter the medium through a 0.45 μm filter into a 1.5 mL tube.

- 4  Micrococcal Nuclease - 320,000 gel units **New England Biolabs** **Catalog #M0247S**

 Benzonase® Nuclease 2.5ku **Catalog #70746-4CN**

Mix the following reagents in a new 1.5mL tube.

Component

Volume / sample

Micrococcal nuclease

 1 μL

Benzonase

 2 μL

Homemade buffer*

 7 μL



*see MATERIALS

5 Add  200 μL of filtrate into the tube, then mix by pipetting.

6 Incubate at  37 °C for  02:00:00 .

2h

7  High Pure Viral RNA Kit **Roche Catalog #11858882001**

Extract RNA from total volume ( 210 μL) and elute to  50 μL .

Note

If the concentration of RNA is measured here using the Qubit RNA HS Assay Kit or similar, it is not measurable due to low concentration. You can understand the large amount of viral genes are still included using real time PCR etc.

Fragmentaion of RNA and priming

8 0.2 mL PCR tubes are used to incubate mixtures.
For steps 9 to 47, refer to section 2 'Protocol for use with NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat)' in a manual of NEBNext Ultra II RNA Library Prep kit for Illumina (kit E7770). In our protocol, section 2.5 and beyond is referred.
Half volume of the reagent listed in the manual is used.

9 Mix the following components, keep  On ice

Component

Volume / sample

RNA

 2.5 μL

First Strand Synthesis Buffer (kit E7770)

 2.0 μL

Random Primers (kit E7770)

 0.5 μL

Total so far:  5 μL

10 Incubate in a thermal cycler set with the following program.
Keep the heat-lid at  105 °C .

10m

1.  94 °C for  00:10:00

2. Hold at  4 °C



Synthesis of 1st strand cDNA

- 11 Mix the following components, keep On ice

Component	Volume / sample
Product from step 10	5.0 µL
First Strand Synthesis Enzyme (kit E7770)	1.0 µL
Nuclease Free Water	4.0 µL

Total so far: 10 µL

- 12 Incubate in a thermal cycler set with the following program.

1h 15m

Keep the heat-lid at 80 °C .

1. 25 °C for 00:10:00

2. 42 °C for 00:50:00

3. 70 °C for 00:15:00

4. Hold at 4 °C

Synthesis of 2nd strand cDNA

- 13 Mix the following components, keep On ice

Component	Volume / sample
Product from step 13	10 µL
Second Strand Synthesis Buffer (kit E7770)	4.0 µL
Second Strand Synthesis Enzyme (kit E7770)	2.0 µL
Nuclease Free Water	24 µL

Total so far: 40 µL

- 14 Incubate in a thermal cycler set with the following program.

1h

Keep the heat-lid at 40 °C .

1. 16 °C for 01:00:00

2. Hold at 4 °C

Clean-up using magnetic beads



15 Clean-up products using

Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

Add 72 μL (1.8x) of AMPure XP per sample.

16 Incubate at Room temperature for 00:05:00

5m

17 Separate magnetic beads and remove supernatant.

18 To wash beads, add 200 μL of 80% ethanol, incubate for 00:00:30 , and remove supernatant (1/2)

30s

19 To wash beads, add 200 μL of 80% ethanol, incubate for 00:00:30 , and remove supernatant (2/2)

30s

20 Allow the beads to dry for 00:02:00 .

2m

21 Elute purified product in 26 μL of 0.1x TE (kit E7770).

22 Separate magnetic beads and transfer 25 μL of supernatant to a new 0.2 mL tube.

End Prep of cDNA Library

23 Mix the following components, keep On ice

Component	Volume / sample
Product from step 13	25 μL
End Prep Reaction Buffer (kit E7770)	3.5 μL
End Prep Reaction Enzyme (kit E7770)	1.5 μL

Total so far: 30 μL

24 Incubate in a thermal cycler set with the following program.

1h

Keep the heat-lid at 75 $^{\circ}\text{C}$.



1. 20 °C for 00:30:00
2. 65 °C for 00:30:00
3. Hold at 4 °C

Adaptor ligation

- 25 Mix the following components in a 1.5 mL low-binding tube, keep On ice

Component	Volume
NEBNext Adaptor for Illumina (E7335 or E7500)	1.0 µL
Adaptor Dilution Buffer(kit E7770)	199 µL

- 26 Mix the following components as master mix in a 1.5 mL tube, keep On ice

Component	Volume / sample
Ligation Enhancer (kit E7770)	0.5 µL
Ligation Master Mix (kit E7770)	15 µL

- 27 Mix the following components, in the order given, keep On ice

Component	Volume / sample
Product from step 24	30 µL
Diluted Adaptor (step 25)	1.25 µL
Master mix (step 26)	15.5 µL

Total so far: 46.75 µL

- 28 Incubate in a thermal cycler set with the following program.

Keep the heat-lid at 45 °C .

1. 20 °C for 00:15:00
2. 20 °C pose *
3. 37 °C for 00:15:00
4. Hold at 15 °C

* Add USER Enzyme (kit E7770) 1.5 µL / sample and mix

Total so far: 48.25 µL

30m

Clean-up using magnetic beads

29 Clean-up products using

 Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

Add  43.5 µL (0.9x) of AMPure XP per sample.

30 Incubate at  Room temperature for  00:10:00

10m

31 Separate magnetic beads and remove supernatant.

32 To wash beads, add  200 µL of 80% ethanol, incubate for  00:00:30 , and remove supernatant (1/2)

33 To wash beads, add  200 µL of 80% ethanol, incubate for  00:00:30 , and remove supernatant (2/2)

34 Allow the beads to dry for  00:02:00 .

35 Elute purified product in  8.0 µL of 0.1x TE (kit E7770).

36 Separate magnetic beads and transfer  7.5 µL of supernatant to a new 0.2 mL tube.

PCR Enrichment of Adaptor Ligated DNA

37 Mix the following components, keep  On ice

Component

Volume / sample

Adaptor Ligated DNA from step 36

 7.5 µL

Q5 Master Mix (kit E7770)

 12.5 µL

Index (X) Primer (E7335 or E7500)

 2.5 µL

Universal PCR Primer (E7335 or E7500)

 2.5 µL

For multiplex analysis of specimens fewer than seven, use 'Index oligo selector' to verify that the index combination is acceptable.

NEBNext Index oligo selector



38 Incubate in a thermal cycler set with the following program.

1m 55s

1. 98 °C for 00:00:30
2. 20 cycles x (98 °C for 00:00:10 , 65 °C for 00:01:15)
3. Hold at 4 °C

Clean-up of PCR product using magnetic beads and quantification of DNA

39 Clean-up products using

Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

Add 22.5 µL (0.9x) of AMPure XP per sample.

40 Incubate at Room temperature for 00:05:00

5m

41 Separate magnetic beads and remove supernatant.

42 To wash beads, add 200 µL of 80% ethanol, incubate for 00:00:30 , and remove supernatant (1/2)

43 To wash beads, add 200 µL of 80% ethanol, incubate for 00:00:30 , and remove supernatant (2/2)

44 Allow the beads to dry for 00:02:00 .

45 Elute purified product in 11.5 µL of 0.1x TE (kit E7770).

46 Separate magnetic beads and transfer 11 µL of supernatant to a new 0.2 mL tube.

47 Quantify the purified amplicon using fluorescent based method using

Qubit 1X dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q33230** .

Concentrations in the range of 10-100 ng/µL of purified amplicon are sufficient for the next section.



Library pooling

- 48 Take the purified PCR product from each tube and pool them into the 1.5 mL low-binding tube.

Adjust the volume to be pooled to average the amount of DNA in each sample.

Note

For example, if you want to obtain 100 ng from each sample,  2 μL of the 50 ng/ μL sample and  4 μL of the 25 ng/ μL sample should be aliquoted.

In the next section, purification with 0.6x AMPureXP to remove DNA of less than 300 bp resulted in the loss of much DNA (reduced to about one-eighth of the amount). Therefore, if the total amount of DNA obtained in this section is small (approximately 500 ng / 50 μL or less), it will be difficult to obtain a visible agarose gel electrophoresis in the section on library size estimation.

If the total amount of DNA from library pooling is small, removal of less than 300 bp of DNA should not be performed.

Briefly measure the volume of pooled mixture using pipette.

Add 0.1x TE (kit E7770) up to  50 μL of total volume.

Purification of the library for size selection

- 49 Clean-up products using

 Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

Add  30 μL (0.6x) of AMPure XP per sample.

- 50 Incubate at  Room temperature for  00:05:00

- 51 Separate magnetic beads and remove supernatant.

- 52 To wash beads, add  200 μL of 80% ethanol, incubate for  00:00:30 , and remove supernatant (1/2)



- 53 To wash beads, add  200 μ L of 80% ethanol, incubate for  00:00:30 , and remove supernatant (2/2)
- 54 Allow the beads to dry for  00:02:00 .
- 55 Elute purified product in  30 μ L of 0.1x TE (kit E7770).
- 56 Separate magnetic beads and transfer supernatant to a new 0.2 mL tube.
- 57 Quantify the purified amplicon using fluorescent based method using  Qubit 1X dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q33230** .
Concentrations 1.5 ng/ μ L or more of purified library is sufficient for the next section.

Estimation of library size

- 58 Quantify the purified library using  Qubit 1X dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q33230**

Note


More accurate molar concentrations can be determined using commercially available library quantification kits instead of the methods shown in this section.

- 59 Mix  5 μ L of the library with loading dye and electrophoresis on a 2% agarose gel alongside molecular markers.
- Obtain a smear image of the library.
- 60 Estimate approximate average library size (base pairs) on the smear image.
- The size of the most concentrated region can be read and used as an estimation.

Note

Image J is helpful to recognize distribution of the library size.
You can obtain a densitogram of the gel image.

<https://imagej.net/ij/>



Preparation of 50pM library for Illumina iSeq100

- 61 Calculate molar concentration of the library using the formula below.

$$Y \text{ (nM)} = X \text{ (ng/}\mu\text{L)} \div Z \text{ (base pairs)} \div 660 \text{ (g/mol)} \times 10^6$$

Y: molar concentration of the library

X: mass concentration of the library

Z: average library size

Note

See the Illumina website.

'Converting ng/ μ L to nM when calculating dsDNA library concentration'

https://knowledge.illumina.com/library-preparation/dna-library-prep/library-preparation-dna-library-prep-reference_material-list/000001240

Setting up the local run manager in iSeq100



- 62 To analyze libraries using NEBNext® multiplex oligos for Illumina, you need to load the index and other information into a local run manager.

Refer to the manuals on Illumina site

How to use a custom library prep kit in Local Run Manager v2

How to use custom library prep and index kits with Local Run Manager v3 and v4

Obtain a .tsv file for configuration from the NEBNext® multiplex oligo page on NEB.

Change DefaultReadLength1, 2 in the .tsv file from 251 to 151.

Start the Local Run Manager on iseq100 and open "Tools" on the dashboard.

In the drop-down menu, select "Index & Library Prep Kits", "Index Kit", and "Add Index Kit".

Select and load the modified .tsv file.

Note

If libraries are created using both E7335 and E7500, a .tsv file that has both index information can be used.

You can consolidate the information under [Indices] in the two .tsv files into one file, and change the "Name" and "Description" under [kit] as desired.



Protocol references

Conceição-Neto N, Zeller M, Lefrère H, De Bruyn P, Beller L, Deboutte W, Yinda CK, Lavigne R, Maes P, Van Ranst M, Heylen E, Matthijssens J. Modular approach to customise sample preparation procedures for viral metagenomics: a reproducible protocol for virome analysis. *Sci Rep*. 2015 Nov 12;5:16532. doi: 10.1038/srep16532.

Itokawa K, Sekizuka T, Hashino M, Tanaka R, Kuroda M.

nCoV-2019 sequencing protocol for illumina protocol V5.

<https://protocols.io/view/ncov-2019-sequencing-protocol-for-illumina-b2msgc6e.html>

Schneider CA, Rasband WS, Eliceiri KW.

NIH Image to ImageJ: 25 years of image analysis. *Nat Methods*. 2012 Jul;9(7):671-5. doi: 10.1038/nmeth.2089.

Illumina, Inc. Converting ng/μl to nM when calculating dsDNA library concentration

https://knowledge.illumina.com/library-preparation/dna-library-prep/library-preparation-dna-library-prep-reference_material-list/000001240

Accessed on May 3, 2024.

Illumina, Inc. How to use a custom library prep kit in Local Run Manager v2

https://knowledge.illumina.com/software/on-instrument-analysis-software/software-on-instrument-analysis-software-reference_material-list/000001317

Accessed on May 3, 2024.

Illumina, Inc. How to use custom library prep and index kits with Local Run Manager v3 and v4

https://knowledge.illumina.com/software/on-instrument-analysis-software/software-on-instrument-analysis-software-reference_material-list/000005897

Accessed on May 3, 2024.