


Jan 22, 2025

Version 3

Preparing Biolog Growth Plates V.3

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwj5ndlmk/v3>



Carlos Goller¹, Carly Sjogren¹

¹Biotechnology Teaching Program (BIT), North Carolina State University

BIT-Protocols

Tech. support phone: +91 95134-135 email: ccgoller@ncsu.edu



Nidhi Grover

North Carolina State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwj5ndlmk/v3>

Protocol Citation: Carlos Goller, Carly Sjogren 2025. Preparing Biolog Growth Plates. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwj5ndlmk/v3> Version created by **Carlos Carlos Goller**



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 22, 2025

Last Modified: January 22, 2025

Protocol Integer ID: 118848

Keywords: Biolog, Growth, multichannel, preparing biolog growth plates overview, biolog growth plates overview, proper use of biolog plate, biolog plate, bacterial densities analysis of growth data, special biolog plate, growth of bacteria, bacterial densities analysis, bacterial isolate, bacterium, following lab skill, lab safety, bacteria, lab skill, metabolic dye, lab, blood universal growth, metabolic capacity, growth data

Funders Acknowledgements:

Biotechnology Program

Abstract

Overview and Goals

Your bacterial isolate has been grown on agar plates. You practiced pipetting. Now, let's learn about how your organism grows! What nutrients help your bacterium grow? We will use plates with ninety-six wells with different conditions to examine the growth and metabolic capacity of your isolate. For this, we will use special Biolog plates and your Blood Universal Growth (BUG) media containing plates with your organism. **Biolog** plates allow researchers to test and visually (colorimetrically!) analyze the growth of bacteria in numerous conditions... up to 30 conditions in one plate! Biolog plates work by having powder chemicals in the wells along with an inoculating fluid and a dye that changes color when an organism is growing ("metabolic dye").

After completing this lab you will gain the following lab skills:

- Lab safety and proper personal protective equipment (PPE)
- Proper use of Biolog plates and multichannel pipettors
- Proper use of a turbidimeter to adjust bacterial densities
- Analysis of growth data of bacteria in 96-well plates

Image Attribution

Image created with BioRender.com



Guidelines

Work carefully and keep the tube capped. If necessary, use another swab if you think it has touched another surface.

Do not leave the liquid uncapped longer than necessary to prevent contamination. The tubes are glass: work carefully.

Materials

Activity 1:

- One 200 µl multichannel micropipette (p200)
- Micropipette tips for p200
- One 96-well plate
- Kimwipes
- One container with water with food coloring (yellow)
- One container with water with food coloring (green)
- Tip disposal container


Activities 2-3:

- One Inoculatorz™ swab
- IF-A inoculating fluid
- Tryptic Soy Agar (TSA) plate with bacteria
- Biolog GEN3 plate
- One 200 µl multichannel micropipette (p200)
- Micropipette tips for p200
- Kimwipes
- Tip disposal container

Before start


Review the figures below to learn about bacterial cell structure, shape, and growth.

Bacterial Structure



Bacterial structures:

- Pilus
- Flagellum
- Capsule
- Cell wall
- Plasma membrane
- Nucleoid
- Ribosome
- Plasmid




Bacterial shapes:

- Cocci
- Coccus
- Diplococci
- Diplococci (encapsulated)
- Staphylococci
- Streptococci
- Sarcina coccus
- Tetrad coccus
- Bacilli
- Coccobacillus
- Bacillus
- Diplobacilli
- Streptobacilli
- Palisades
- Appendaged bacteria
- Hypha
- Stalk
- Others
- Club rod
- Vibrio

- Spirillum
- Helical Form
- Spirochete
- Filamentous

Bacterial Growth



Bacterial growth phases over time:

- Lag phase (slow growth)
- Log phase (exponential growth)
- Stationary phase (growth plateau)
- Death phase (decline, depends on the microbe and growth conditions)

Next, read pages 2-5 of this document about the Biolog **GEN III MicroPlate™** (text, ~10 min). We will begin by practicing the use of a multi-channel pipette with containers ("reservoirs") with water and food coloring.



Activity 1: Multichannel pipetting

24m

- 1 Set your p200 multichannel pipet to  100 µL 1m



- 2 Load eight (8) tips onto each gasket of the multichannel pipet. 1m
- 3 Press the plunger TO the first stop. 
- 4 Submerge all 8 tips into the water with food coloring.
- 5 With thumb control, release the plunger to take up liquid. Check each tip to make sure each has the same volume of liquid. 

Note

Note the height of the liquid in the tips. Is it consistent?
- 6 Align each tip with one well across one column of the 96-well plate (refer to the image or the descriptions below) 1m
- 6.1 Yellow water should go in odd numbered columns 1, 3, 5, 7, 9, 11.
- 6.2 Green water should go in the even numbered columns 2, 4, 6, 8, 10, 12. 5m
- 6.3 If using colors other than yellow and green, make your own pattern.

Reference image:

Biolog GEN3 Bacterial Growth



Each researcher will use the multi-channel pipette to transfer 100 μ l of colored solutions to each numbered columns of the 96-well plate.

Check that you are delivering the appropriate volume to all wells.

Each researcher will use multi-channel pipette to transfer 100 μ L of colored solutions to each numbered columns of the 96-well plate. Check that you are delivering the appropriate volume to all wells.

- 7 Press the plunger **THROUGH** the first stop, raise your tips up out of the plate, and then you may release your plunger with thumb control.
- 8 Have each researcher do two to three columns of the plate to practice.

1m



15m

Activity 2: Turbidimeter

- 9 We will swab bacterial cells from your agar plates (BUG Agar: **BUG** (Biolog Universal Growth) Agar w/ 5% sheep blood), dilute in Biolog IF-A inoculation fluid, and check turbidity with a turbidimeter.

Calibrate the turbidimeter to a "blank" sample of inoculation fluid in a tube

- 9.1 Using the Biolog Turbidimeter, blank the turbidimeter with a clean tube containing uninoculated IF-A.

Note

Wipe the tube clean of dirt and fingerprints
Because the tubes used are not optically uniform, they should be blanked individually.



9.2 Set the 100% transmittance adjustment knob so that the meter reads **100%**.

9.3 Record reading.

10 Use a turbidimeter with bacterial sample in the inoculation fluid in a tube

10.1 Collect bacteria from your streak plate by touching the Biolog Inoculatorz™ swab to the bacteria

Safety information

Be gentle: do not push the swab through the agar media

10.2 Open your Biolog Inoculating Fluid tube (IF) tube and submerge the swab into the liquid. Swirl the swab around in the fluid for  00:00:15 to ensure inoculation.

15s

Safety information

Do not leave the liquid uncapped longer than necessary to prevent contamination. The tubes are glass: work carefully.

10.3 Cap the tube and gently invert the three (3) times to ensure thorough mixing.

10.4 The target cell density is **95%T** for our protocol.


1. Add more bacteria with the swab if necessary to get as close to this value as possible by continuing to add more bacteria with the swab to the tube.
2. Work carefully and keep the tube capped. If necessary, use another swab if you think it has touched another surface.

11 Once you've reached the target cell density at **95%T**, cap the tube and proceed to the next Activity.

Activity 3: Biolog plate preparation

12

Biolog GEN3 Bacterial Growth



Use the multi-channel pipette to transfer 100 μ L of bacterial solution to each well of the Biolog plate.

Check that you are delivering the appropriate volume to all wells.

Each researcher will use a multi-channel pipette to transfer 100 μ L of inoculating fluid and bacteria to each numbered column of the 96-well plate. Check that you are delivering the appropriate volume to all wells.

Pour the cell suspension into the multichannel pipette reservoir.

13 Load 8 tips onto each gasket of the multichannel pipet and fill the tips by drawing up

100 μ L of cell suspension from the reservoir.



14 Use the multichannel to add 100 μ L of diluted bacterial culture to each of the 96 wells of the **GEN III MicroPlate**TM.



Note

Use the technique practiced in Activity 1.

Safety information

If the tips touch the wells, replace your tips before delivering inoculation fluid to the next column. **Work carefully to avoid contamination!**



- 15 Eject the pipettor tips if any one tip touches the wells of the plate as this can cross contamination wells with chemicals. Reload tips as needed.
- 16 Check that all wells have liquid
- 17 Cover the GEN3 microplate with its lid and eject the pipette tips.
- 18 Incubate at 33 °C for 24:00:00 - 24:00:00 2d hours. For some plates, we will incubate in a heated Agilent/Biotek LogPhase 600 plate reader and take 590 nm/600 nm readings every 15 minutes for 2 days to learn about the growth patterns.




2d



Reference image:

Biolog GEN3 Bacterial Growth



Incubate the plate at 33°C for 1 to 2 days.

Illustration of Biolog GEN3 plate with several purple wells. Right: icon of calendar with text "Incubate the plate at 33°C for 1 to 2 days"

Clean Up Instructions

1. Biolog growth plates must be given to an instructor to place in our plate readers to measure viability and growth!



2. All liquids from Activity 1 (materials with no bacterial contamination) can be poured down the sink and plastic can go in the regular trash.
3. All liquids containing bacteria must be poured into a beaker for bleach decontamination (can be found on the blue bench paper by this sink labeled "liquid").
4. Glass tubes can be placed into a beaker (can be found on the blue bench paper by this sink labeled "glass").
5. All materials from Activity 2 that contact bacteria need to go into the biohazard trash.
6. All bacterial growth plates can be thrown away in the biohazard trash
7. Return all reusable materials to our storage drawers.
8. Wipe down all surfaces with Conflkt and then ethanol.
9. Remove PPE.
10. Wash hands.

Note

Critical Thinking Questions for Preparing Biolog Growth Plates

1. What are some advantages of using a multichannel pipette?
2. What are some disadvantages of using a multichannel pipette?
3. During the inoculation, describe how you ensure that you isolate a "pure culture." Why is a "pure culture" needed for this assay?
4. What is a turbidimeter, and why would measuring turbidity be important for this assay?

Protocol references

This protocol has been adapted from the Biolog GEN3 Protocol.

<https://www.biolog.com/products/microbial-identification-microplates/gen-iii-test-panel/>

Acknowledgements

We appreciate the support from the NC State Biotechnology Program and students who have helped improve this and other protocols.