



Oct 09, 2024

Protocols for stereotaxic injections into mouse brain and ex-vivo electrophysiology

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjn7qpgk5/v1>

Jyoti Gupta¹, Michael J. Higley¹

¹Yale University



Maria Matos

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjn7qpgk5/v1>

Protocol Citation: Jyoti Gupta, Michael J. Higley 2024. Protocols for stereotaxic injections into mouse brain and ex-vivo electrophysiology. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvjn7qpgk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: September 22, 2024

Last Modified: October 10, 2024

Protocol Integer ID: 109187

Keywords: ASAPCRN, protocols for stereotaxic injection, stereotaxic injection, synuclin pff, mouse brain, cell patch clamp recording, method for injection, injection, mice, mouse

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020616

Disclaimer

ETHICS DISCLAIMER

The [protocols.io](https://www.protocols.io) team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board

Abstract

This protocol describes the method for injection of α -Synuclin PFF and monomer into the mouse brain. The second part of the protocol describes preparation of acute slices from these mice and whole-cell patch clamp recordings.

Materials

Oxygenated (95% O₂/5% CO₂) ice-cold choline artificial cerebrospinal fluid (Choline ACSF) containing (in mM):

	A	B
	choline	110
	NaHCO ₃	25
	NaH ₂ PO ₄	1.25
	KCl	2.5
	MgCl ₂	7
	CaCl ₂	0.5
	glucose	20
	sodium ascorbate	11.6
	sodium pyruvate	3.1

ACSF ( 32 °C) containing (in mM):

	A	B
	NaCl	127
	NaHCO ₃	25
	NaH ₂ PO ₄	1.25
	KCl	2.5
	MgCl ₂	1
	CaCl ₂	2
	glucose bubbled with 95% O ₂ /5% CO ₂	20

Backfill glass electrodes (3-3.5 MΩ) with an internal solution containing (in mM):

	A	B
	cesium gluconate	126



	A	B
	HEPES	10
	sodium phosphocreatine	10
	magnesium chloride	4
	Na ₂ ATP	4
	Na ₂ GTP	0.4
	EGTA (pH 7.3 with cesium hydroxide)	1

Protocol for stereotaxic injections

- 1 Anesthetize the mouse with isoflurane by placing in the anesthesia chamber.
- 2 Place the anesthetized mouse in the stereotaxic frame and continually maintain isoflurane anesthesia using a vaporizer (Harvard Apparatus).
- 3 Locate the injection sites relative to bregma and make holes into the skull using a drill.
- 4 Lower the injection pipette into the brain at the predetermined coordinates. Perform the injection at a rate of 20 nL/s using a Nanoject III microinjector (Drummond Scientific). 
- 5 Following the injection, leave the pipette in place for 10 minutes to prevent backflow while withdrawing. 
- 6 Seal the incisions with sutures and allow the mice to recover on a heating pad until conscious.
- 7 Continue to monitor the mice and administer the post-operative analgesia.

Protocol for acute slice preparation and electrophysiology

30m

- 8 Anesthetize the mouse with isoflurane.
- 9 Transcardially perfuse the anesthetized mouse with oxygenated (95% O₂/5% CO₂) ice-cold choline artificial cerebrospinal fluid (Choline ACSF) containing (in mM):

	A	B
	choline	110
	NaHCO ₃	25
	NaH ₂ PO ₄	1.25
	KCl	2.5



A	B
MgCl ₂	7
CaCl ₂	0.5
glucose	20
sodium ascorbate	11.6
sodium pyruvate	3.1

10 Decapitated the mouse and quickly remove the brain.

11 Mount the brain in the vibratome and cut 300 µm thick slices in the coronal plane.

12 Transfer the slices to warm ACSF ( 32 °C) containing (in mM):



A	B
NaCl	127
NaHCO ₃	25
NaH ₂ PO ₄	1.25
KCl	2.5
MgCl ₂	1
CaCl ₂	2
glucose bubbled with 95% O ₂ /5% CO ₂	20

13 After incubating at  32 °C for  00:30:00 , transfer the slices to  Room temperature (RT).

30m



14 Transfer the slices to the slice chamber that is perfused with oxygenated ACSF.



15 Identify layer 2/3 pyramidal cells in the primary somatosensory cortex.

16 Backfill glass electrodes (3-3.5 M Ω) with an internal solution containing (in mM):

A	B
cesium gluconate	126
HEPES	10
sodium phosphocreatine	10
magnesium chloride	4
Na ₂ ATP	4
Na ₂ GTP	0.4
EGTA (pH 7.3 with cesium hydroxide)	1

17 Add **1M** 1 micromolar (μ M) tetrodotoxin to the bath to record miniature postsynaptic currents.



18 Voltage-clamp the cells at -70mV to record miniature excitatory postsynaptic currents (mEPSCs) and at 0mV to record miniature inhibitory postsynaptic currents (mIPSCs).

Note

All protocols were carried out in accordance with Yale IACUC guidelines.