

Jul 29, 2017

Purification of Nosema bombycis spores

PLOS One

DOI

dx.doi.org/10.17504/protocols.io.if3cbqn

Zhu Feng¹

¹Jiangsu University of Science and Technology



Chen Gong

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.if3cbqn

External link: https://doi.org/10.1371/journal.pone.0181703

Protocol Citation: Zhu Feng 2017. Purification of Nosema bombycis spores. protocols.io

https://dx.doi.org/10.17504/protocols.io.if3cbqn

Manuscript citation:

Chen G, Wang W, Chen H, Dai W, Peng X, Li X, Tang X, Xu L, Shen Z (2017) Functional characterization of an aquaporin from a microsporidium, Nosema bombycis. PLoS ONE 12(7): e0181703. doi: 10.1371/journal.pone.0181703



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 15, 2017

Last Modified: March 19, 2018

Protocol Integer ID: 6363

Keywords: purification of nosema bombyci, purified spore, spore, infected beetle, nosema bombyci, purification, sterile water, several times with sterile water, light microscope, microscope usb camera, ocular micrometer

Abstract

Infected beetles were collected and homogenized in sterile water. The homogenates were filtered through four layers of cheesecloth and centrifuged at 3000g for 15 min. The pellets were resuspended in sterile water, and the spores were purified by Percoll gradient centrifugation using 90% Percoll at 15,000g for 40 min. The spore band was collected and washed several times with sterile water. The purified spores (n = 50) were measured under a light microscope (IL/Leica Microsystems, Inc., Deerfield) with an ocular micrometer and photographed with the Microscope USB Camera.

Attachments



A new isolate of Nos...

538KB

Troubleshooting

