

Oct 18, 2020

Pyrolysis kinetics

DOI

dx.doi.org/10.17504/protocols.io.bnjbmcin

Jiaxin Li¹

¹South China University of Technology



Jiaxin Li

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bnjbmcin>

Protocol Citation: Jiaxin Li 2020. Pyrolysis kinetics. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bnjbmcin>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 18, 2020

Last Modified: October 18, 2020

Protocol Integer ID: 43331

Keywords: cells in planktonic culture, planktonic culture, phage, well microtiter plate, cell

Abstract

The sensitivity of cells in planktonic cultures to phages was tested in 96-well microtiter plates.

Troubleshooting



- 1 A 1:100 dilution of overnight bacterial culture in LB should be prepared and incubated at 37°C to reach an optical density of 0.5 at 600 nm (OD₆₀₀).
- 2 100 µl aliquots of diluted culture were distributed in microplate wells and 100 µl single or mixed phage suspensions (at MOI 10; 1; 0.1 and 0.01) were added.
- 3 Microplates were incubated for a further 24 h and optical density was measured at 600 nm in a plate reader every 870s.