


Aug 16, 2018

Version 1

Quick Fungal DNA extraction from colonies on plates V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.snzedf6>



Johannes Debler¹

¹Curtin University



Johannes Wolfram JWD Debler

Curtin University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.snzedf6>

Protocol Citation: Johannes Debler 2018. Quick Fungal DNA extraction from colonies on plates. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.snzedf6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 16, 2018

Last Modified: August 16, 2018

Protocol Integer ID: 14777

Keywords: fungal DNA extraction, quick fungal dna extraction from colony, colony pcr on fungal colony, quick fungal dna extraction, fungal colony, colony pcr, dna, pcr template, pcr

Abstract

Adapted from **Chi, M. H., Park, S. Y., & Lee, Y. H. (2009).**

Colony PCR on fungal colonies grown on plates does not work as well as it does for bacteria (it usually doesn't work at all). DNA therefore needs to be extracted first. As this DNA will only be used as a PCR template to check for presence / absence of individual genes I am not too concerned about high molecular weight or purity. This protocol is quick and can be used to process several samples at the same time.

Make stock solutions

1

1M TrisHCl pH 8.0	100 ml
12.14 g	Tris base or Trizma base
up to 100 ml	deionised H ₂ O
Adjust pH	concentrated HCl

2M KCl (Potassium chloride)	200 ml
15.11 g	KCl (Potassium chloride)
up to 200 ml	deionised H ₂ O

500 mM NaEDTA	100 ml
18.6 g	EDTA disodium salt dihydrate
up to 100 ml	deionised H ₂ O

Autoclave stock solutions at 121°C

Make Extraction Buffer

2

Extraction Buffer	100 ml
2M KCl	50 ml
1M Tris-HCl pH 8.0	10 ml
500 mM NaEDTA	2 ml
deionised H ₂ O	up to 100 ml

Preparation



- 3 Per sample prepare **2 × 1.5 ml Eppendorf tubes**. Label them with sample number and add:

Tube 1: **500 ul Extraction Buffer + 1 stainless steel bead**

Tube 2: **300 ul Isopropanol**

Harvest fungal material

- 4 With a **sterile toothpick or pipette tip** remove a fungal colony or a piece of mycelia of about 3-5 mm x 3-5 mm in size from a plate and put it into **Tube 1** (containing the Extraction Buffer and steel bead).

Homogenise fungal material

- 5 Homogenise samples for **5 min at 1500 rpm** in the MiniG 1600.

Equipment

new equipment

NAME

SPEX

BRAND

SP 1600

SKU

MiniG 1600 Automated Tissue Homogenizer and Cell Lyser

SPECIFICATIONS



Pellet Cell Debris

- 6 Pellet cell debris by centrifugation at **17.000 x g for 5 min**.

DNA Precipitation

- 7 ■ Transfer supernatant to **Tube 2** (Containing 300 ul Isopropanol)



- Vortex for 5 - 10 seconds
- Pellet DNA by centrifugation for **5 min at 17.000 x g**

Wash DNA

- 8
 - Discard supernatant
 - Add **800 ul 70% Ethanol**
 - Vortex mix sample, make sure pellet comes loose from bottom of the tube
 - Pellet DNA by centrifugation for **5 min at 17.000 x g**

Dissolve DNA

- 9
 - Discard supernatant
 - Remove as much of the ethanol as possible by drawing drops out with a pipette tip
 - Let leftover ethanol evaporate for 5-10 min
 - Dissolve pellet in **50 ul H₂O or 1x TE buffer**
 - Measure DNA concentration (e.g. Nanodrop)
 - store DNA at -20°C