

Dec 01, 2025

RNA extraction and purification from cultured diatom viruses

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr8k13lmk/v1>

Pauline Nogaret¹, Estelle Bigeard¹, Laure Arsenieff¹, Anne Claire Baudoux¹

¹Station biologique de Roscoff, FRANCE, UMR7144 CNRS Sorbonne Université

Estelle Bigeard: co-first author

Ecology of Marine Plankton (ECOMAP) team - Roscoff

Tech. support email: phytopk@sb-roscoff.fr



Pauline Nogaret

Station biologique de Roscoff, FRANCE, UMR7144 CNRS Sorbonne...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr8k13lmk/v1>

Protocol Citation: Pauline Nogaret, Estelle Bigeard, Laure Arsenieff, Anne Claire Baudoux 2025. RNA extraction and purification from cultured diatom viruses. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr8k13lmk/v1>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 27, 2024

Last Modified: December 01, 2025

Protocol Integer ID: 102510

Keywords: extraction, RT PCR, RNA virus, virus, chloroform purification, nucleic acid extraction, RNA viruses, marine ecology, environmental analysis, purification from cultured diatom virus, genome sequences of marine rna virus, marine rna virus, cultured diatom virus, rna extraction, infecting rna virus, rna virus, processing viral lysate, virus concentration by tangential flow filtration, virus concentration, sequencing genome, viral lysate, final nucleic acid quantification, genome sequence, based purification, purification

Funders Acknowledgements:

Agence Nationale de la Recherche

Grant ID: ANR-22-CE01-0018

Abstract

This protocol describes a workflow for processing viral lysates to obtain genome sequences of marine RNA viruses. The procedure includes lysate clarification, virus concentration by tangential flow filtration, chloroform-based purification, RNA extraction and conversion to cDNA, and final nucleic acid quantification. cDNA produced using this protocol is suitable for sequencing genomes of diatom-infecting RNA viruses. This protocol is adapted from the method described in Arsenieff et al. (2019).

Image Attribution

Transmission electron micrograph of *M.comicus* (RCC4660) infection with an Bacillarnavirus (RCC7291).

Guidelines

Guides and Manuals:

MasterPure Complete guide

SuperScript III Reverse Transcriptase guide

Materials

Equipments :

weighing machin & pH meter
Centrifuge Avanti J-25 XP & Rotor JA-10 (Bekman Coulter)
Centrifuge **5804 R** (Eppendorf)
Centrifuge **5427 R** (Eppendorf)
Chemical hood
Micropipets
Hand pump and pump
Eppendorf thermomixer & smartblock 1,5ml
Thermocycler
Microwave
Electrophoresis system & generator (**Sub Cell Model 96 cell**)
Image Quant LAZ4000 (Ge Healthcare)
Qubit 4 fluorimeter (**Invitrogen Q33238**)

Consumables:

Autoclaved 500mL Centrifuge bottles, PPCO/HDPE/PC/FEP, with screw cap, Nalgene® (VWR, cat no : 525-2313)
1L Ventilated Culture flask (Starlab, cat no : CC7682-4822)
0,45µm Filtration unit (VWR, cat no : 2566881-514-0343)
Vivaflow® 200 30Kda, MWCO, PES (Sartorius, cat no : VF20P2)
Vivaspin® 20 30Kda, MWCO, PES (Cytiva, cat no : 28932361)
Syringe filter 0,2µm 25mm, PES (Fisher scientific, cat no : 15392388)
Falcon Tubes 50mL, 15mL
Microtubes 1,5mL ,2mL, 0,2µL (Eppendorf-type)
Filter tips

Reagents:

KH₂PO₄ (Merck Sigma, cat no : 60218-100G)
Na₂HPO₄ (Merck Sigma, cat no : 71642)
Chloroform (Merck Sigma, cat no : C2432)
Isopropanol (Merck Sigma, cat no : I9516)
Absolute ethanol (Fisher scientific, cat no : 10342652)
Tris-Base (Merck Sigma, cat no : T6791)
Acetic acid (>99.7%) (Fisher Scientific, cat no : 11463473)
0,5 M EDTA pH 8,0 (Fisher Scientific, cat no : 15575-020)

Master Pure complete DNA and RNA purification kit (Epicentre, cat no : MC85200)
Qubit RNA BR Assay Kit (Invitrogen, cat no : Q10210)
Qubit dsDNA BR Assay Kit (Invitrogen, cat no : Q32853)

SuperScript III Reverse Transcriptase (Invitrogen, cat no : 18080044)



RNAse OUT Recombinant Rnase Inhibitor (Invitrogen, cat no : 10777019)

Random primers, 250ng (Invitrogen, cat no : 48190011)

dNTP mix, 10mM (Invitrogen, cat no : 18427013)

Primers *RdRp* Forward and Reverse (10μM)

DNA polymerase Taq platinum; MgCl₂ (50mM); Rxn PCR buffer (10X); dNTP mix (10mM) (Fisher scientific, cat no : 10966018)

Nuclease-Free Water (Invitrogen, cat no : AM9937)

Ice

Loading buffer

Agarose (Interchim, cat no : 31272L)

Ethidium bromide (Merck Sigma, cat no : E1510)

SmartLadder - 200 to 10000 bp (Eurogentec, cat no : MW-1700-10)

Before start

Solution preparation:

cold isopropanol at -20°C

cold 70% ethanol at -20°C

10mM phosphate solution (10 mM KH₂PO₄ + 10 mM Na₂HPO₄), pH 7.2 filtered through a 0,2μm filter

TAE Buffer (Tris Acetate EDTA) 50X :

Tris-Base (2M final) : 242.5 g (Merck Sigma, cat no : T6791)

Acetic acid (>99.7%, 2M final) : 7 ml (Fisher Scientific, cat no : 11463473)

0,5 M EDTA pH 8,0 (50 mM final) : 100 ml (Fisher Scientific, cat no : 15575-020)

dH₂O qsp 1L



RNA virus culture.

3w

- 1 One liter of diatom host culture is grown exponentially in K+Si medium (Keller et al., 1987) at 18 °C under a 12:12 h light:dark cycle with an irradiance of 100 $\mu\text{mol photons}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$, provided by white fluorescent light (Philips Master TL-D 18W/865). A freshly produced viral lysate, filtered through a 0.45 μm PES membrane, is added to the host culture at a final concentration of 10–20% (v/v). The infected culture is incubated under the same conditions as the host culture until complete lysis is observed.

1 L

RNA virus concentration.

3d 7h 38m

2 Lysate clarification :

30m

- Centrifuge the viral lysate (2×500mL) 5353 x g, 4°C, 00:30:00 (Avanti J-25XP, Beckman Coulter-JA-10 Rotor).
- Carefully transfer the supernatant into a 0.45 μm filtration unit (Stericup), avoiding disturbance of the pellet.
- Filter the supernatant completely using a hand pump.

3 Viral concentration :

1d 4h

- Concentrate the clarified viral lysate to 50 mL using a Vivaflow® 200 cartridge (30 kDa MWCO, PES, Sartorius).
- Further concentrate the sample to ~1 mL using a Vivaspin® 20 (30 kDa MWCO, PES, Cytiva) by centrifugation 6000 x g, 4°C, 00:10:00
- Keep the Vivaspin at 4°C overnight.
- Rinse the Vivaspin® 20 by adding 1 mL of 10 mM phosphate solution, and transfer the rinse to the Falcon tube containing the sample.
- Repeat the rinsing step until a final volume of 10 mL is reached.

4 Viral Purification :

2d 3h 8m

Under a chemical hood :

- Add 1 volume of chloroform (v/v, 10mL) into the Falcon tube containing the sample.
- Vortex thoroughly.
- Centrifuge. 2200 x g, 12°C

- Carefully transfer the aqueous phase (upper phase containing viruses) into a new Falcon tube.
- Add again 1 volume of chloroform
- Repeat previous steps a total of 7 times. ⌚ 03:00:00
- Leave the tube (containing 5mL of sample) open overnight under the chemical hood to evaporate residual chloroform. ⌚ Overnight
- Concentrate the sample in 1mL by using Vivaspin® 20 ⌚ 6000 x g, 12°C, 00:05:00
- Collect sample and add 1mL of nuclease-free water. Centrifuge in the same Vivaspin® 20 ⌚ 6000 x g, 12°C, 00:03:00 . Repeat this washing step 3 times.
- Incubate the Vivaspin® 20 overnight at 4°C and, then, collect the final sample (1mL).

Safety information

Chloroform is a CMR chemical reagent.



RNA extraction using MasterPure Complete DNA & RNA purification kit.

2h

- 5 Extraction of 500µL viral nucleic acids using MasterPure Complete DNA & RNA Purification Kit (Epicentre).

2h

Cellular lysis:

- Split in two tubes 250µL of viral sample : "tube 1" and "tube 2".
This two tubes will be pooled and the end of extraction for concentration.
- Thoroughly mix the several buffers to homogenize them before dispensing.
- Dilute 4µL of Proteinase K into 600µL of 2X T and C Lysis Solution.
- Vortex 10s.
- Add 250µL of 2X T and C Lysis Solution containing the Proteinase K to each tube of viral sample and mix thoroughly.

- Incubate at 65°C for 15 minutes in a Eppendorf thermomixer at 1000rpm with smartblock 1,5mL.
- Place the samples on ice for 5 minutes to inactivate the proteinase K.

Nucleic acids precipitation and extraction :

- Add 250 µL of MPC Protein Precipitation Reagent to 500 µL of lysed sample.
- Vortex vigorously for 10 seconds.
- Centrifuge.  11000 x g, 12°C, 00:10:00
- Transfer the supernatant to a clean microcentrifuge tube and discard the pellet.
- Centrifuge.  11000 x g, 12°C, 00:10:00
- Transfer the supernatant to a clean microcentrifuge tube and discard the pellet.
- Add 840 µL of cold isopropanol to the recovered supernatant. Invert tubes 40 times.
- Centrifuge tubes to pellet the total nucleic acid.  20800 x g, 12°C, 00:10:00
- Carefully pour off the isopropanol without dislodging the pellet.

Contaminating DNA digestion with DNase :

- Prepare DNase I solution (1/40) : 10µL Rnase Free Dnase I into 390µL Dnase Buffer.
- Completely resuspend the total nucleic acid pellet in 200 µL of DNase I solution (1/40).
- Incubate at 37°C for 10 minutes.
- Add 200 µL of 2X T and C Lysis Solution.
- Vortex for 5 seconds.
- Add 200 µL of MPC Protein Precipitation Reagent.
- Vortex 10 secondes.
- Place on ice for 4 minutes.
- Centrifuge  11000 x g, 12°C, 00:10:00 to pellet the debris.
- Transfer the supernatant containing the RNA into a clean microcentrifuge tube and discard the pellet.
- Centrifuge  11000 x g, 12°C, 00:03:00 keep and transfer the supernatant into a new tube.
- Add 840µL of cold isopropanol to the supernatant. Invert tubes 40 times.
- Centrifuge  20800 x g, 12°C, 00:10:00 and discard supernatant.

RNA precipitation:

- Rinse twice with cold 70% ethanol, being careful to not dislodge the pellet.
- Centrifuge.  11000 x g, 12°C, 00:00:30
- Remove all the ethanol by pipetting and keep at room temperature the tube open during 10 minutes for the ethanol evaporation .
- Resuspend the total nucleic acids (RNA) of "tube 1" in 35 µL of TE Buffer.



- Re use this solution for resuspend the total nucleic acids of the "tube 2" in the aim of concentrate sample.
- Keep 2µl for Qubit quantification and store RNA extract at -80°C.

RNA Quantification using Qubit 4 fluorimeter.

40m

- 6
- Quantify the RNA concentration using the Qubit 4 fluorimeter and the RNA BR Assay kit which provides a detection range of 10 and 1200ng of RNA.

Retrotranscription to cDNA.

2h

- 7
- By sample, add in 0,2 ml microtube :
- 1µL Random primers (250ng)
 - 11µL RNA extract (Total RNA should range between 10pg - 5µg - complete with RNA-free water if necessary)
 - 1µL 10mM dNTP
- qsp 13µL : Add Nuclease Free water to complete.

Thermocycler parameter : 65°C 5min.
2min in ice.

In the same microtube, add :

- 4µL 5X First strand buffer
- 1µL 0,1M DTT
- 2µL Superscript III RT (400U/µl)
- 1µL RNase OUT Recombinant Rnase Inhibitor.

Thermocycler parameter : 25°C 10min > 50°C 60min > 70°C 15min.

Store cDNA at -20°C.

cDNA Quantification using Qubit 4 fluorimeter.

- 8
- Quantify the DNA concentration using the Qubit 4 fluorimeter and the dsDNA BR Assay kit which provides a detection range of 4 and 2000ng of DNA.



Protocol references

Arsenieff L, Simon N, Rigaut-Jalabert F, Le Gall F, Chaffron S, Corre E, Com E, Bigeard E and Baudoux A-C (2019) First Viruses Infecting the Marine Diatom *Guinardia delicatula*. *Front. Microbiol.* 9:3235. doi: 10.3389/fmicb.2018.03235

Keller, M. D., Seluin, R. C., Claus, W., and Guillard, R. R. L. Media for the culture of oceanic ultraphytoplankton. *J. Phycol.* 23, 633–638. doi:10.1111/j.1529-8817.1987.tb04217