

May 11, 2024

RNA extraction from milk for HPAI surveillance

 Emerging Infectious Diseases

DOI

dx.doi.org/10.17504/protocols.io.n2bvjn6obgk5/v1

Andrew Lail¹, William C. Vuyk¹, David O'Connor¹, Isla Emmen²

¹University of Wisconsin-Madison; ²UW-Madison



Andrew Lail

University of Wisconsin-Madison

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjn6obgk5/v1>

External link: <https://doi.org/10.3201/eid3113.250087>

Protocol Citation: Andrew Lail, William C. Vuyk, David O'Connor, Isla Emmen 2024. RNA extraction from milk for HPAI surveillance. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvjn6obgk5/v1>

**Manuscript citation:**

Stachler E, Gnirke A, McMahon K, Gomez M, Stenson L, Guevara-Reyes C, Knoll H, Hill T, Hill S, Messer KS, Arizti-Sanz J, Albeez F, Curtis E, Samani P, Wewior N, O'Connor DH, Vuyk W, Khoury SE, Schnizlein MK, Rockey NC, Broemmel Z, Mina M, Madoff LC, Wohl S, O'Connor L, Brown CM, Ozonoff A, Park DJ, MacInnis BL, Sabeti PC Establishing Methods to Monitor Influenza (A)H5N1 Virus in Dairy Cattle Milk, Massachusetts, USA. Emerging Infectious Diseases 31(Suppl 1). doi: [10.3201/eid3113.250087](https://doi.org/10.3201/eid3113.250087)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 02, 2024

Last Modified: May 11, 2024

Protocol Integer ID: 99087

Keywords: Dairy, Milk, Influenza, Viral RNA, RNA Extraction, magmax rna extraction kit, rna extraction method, rna extraction, hpa1 rna segment, total rna yield, extraction, performing qpcr, rna, available milk sample, qpcr, hpa1 surveillance

Abstract

This is an RNA extraction method that isolates HPAI RNA segments from commercially-available milk samples. The protocol involves processing with a MagMax RNA extraction kit on a Kingfisher instrument. Then, the samples are concentrated using a Zymo Clean and Concentrate kit. The total RNA yields should be suitable for performing qPCR. After trying other methods involving centrifugation, filtration, and trizol, we determined this method to be the easiest and most effective.

Materials

MagMax Wastewater Ultra x96 kit

Zymo RNA Clean and Concentrator-5 kit

ThermoFisher TurboDNase kit

Heat block

Lab plasticware

Troubleshooting

Safety warnings



Performing this protocol with unpasteurized milk from an infected cow might expose the user to active flu virus. Please use pasteurized samples or take appropriate precautions.

RNA extraction

- 1 Use the MagMax Wastewater Ultra x96 kit¹ (or comparable magnetic bead isolation kit²⁻⁵) on the ThermoFisher Kingfisher Apex (or comparable instrument). For each sample, load 400 uL of pasteurized milk or heavy whipping cream into the binding plate.
- 2 Elute in 50 uL of elution solution or nuclease-free water.

Optional Sample Concentration and Cleaning

- 3 Pool the eluate from 6 extraction replicates of the same sample.
- 4 Take 200 uL forward into the DNase treatment. Using ThermoFisher TurboDNase, combine 200 uL of pooled eluate, 22 uL 10X TurboDNase buffer, and 2 uL TurboDNase.
- 5 Incubate at 37 C for 30 minutes.
- 6 Use a Zymo RNA Clean and Concentrator-5 kit⁷ to clean and concentrate the sample. Elute this into 30 uL final volume of nuclease-free water for downstream analysis.

Protocol references

1. MagMAXTM Wastewater Ultra Nucleic Acid Isolation Kit [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/A52606>
2. Guidelines and minimum requirements for diagnosis of H5N1 HPAIV infection in cattle [Internet]. Istituto Zooprofilattico Sperimentale delle Venezie; 2024 Apr [cited 2024 May 10]. Report No.: Version 2024.04.10. Available from: <https://www.izsvenezie.com/documents//reference-laboratories/avian-influenza/diagnostic-protocols/guidelines-diagnosis-h5n1-cattle.pdf>
3. MagMAXTM-96 Viral RNA Isolation Kit [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/AM1836>
4. MagMAXTM Viral/Pathogen Ultra Nucleic Acid Isolation Kit [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/A42356>
5. MagMAXTM Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/A42358>
6. TURBOTM DNase (2 U/ μ L) [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/AM2238>
7. ZYMO RESEARCH [Internet]. [cited 2024 May 10]. RNA Clean & Concentrator-5. Available from: <https://www.zymoresearch.com/products/rna-clean-concentrator-5>