



Feb 11, 2023

Version 1

## RNAi Library screen V.1

DOI

[dx.doi.org/10.17504/protocols.io.36wggje5ovk5/v1](https://dx.doi.org/10.17504/protocols.io.36wggje5ovk5/v1)

e.warren Warren<sup>1</sup>

<sup>1</sup>London Institute of Medical Sciences



e.warren Warren

LMS

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wggje5ovk5/v1>

**Protocol Citation:** e.warren Warren 2023. RNAi Library screen. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.36wggje5ovk5/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** February 08, 2023



**Last Modified:** February 11, 2023

**Protocol Integer ID:** 76666

**Keywords:** rnai knockdown library, rnai library screen protocol, phenotypic effect

## Abstract

Protocol for screening a RNAi knockdown library in for behavioural and phenotypic effects in *C. elegans*

## Materials

NGM agar plates

M9 Buffer

Tracking plates (UNIPLATE 650,96,PS,CLEAR:WHAT7701-1651)

Tracking plate lids (LID,POLYSTYRENE,CLEAR,UNIVERSAL PK100 : 512-1093)

96 well culture plates (Life technologies 268200)

Boekel pin replicator (140500 boekel)

Integra VIAFILL

Intrgra VIAFLO

Tips for VIAFLO (12.5 µl GRIPTIP, Sterile, LONG - 5 XYZ Racks of 384 Tips: 6404)

IPTG: 1M stock (dissolved in H2O)

Nystatin: 5000x

Carbenicillin: 100 mg/ml (1000x) (dissolved in H2O)

Tetracycline: 50mg/ml (5000x) (dissolved in 70% EtOH)

FUDR: 1g/10 ml (1000x) (dissolved in H2O)

## Troubleshooting



## Perparing *C. elegans*

### 1 **Prepare in advance**

Prepare ten 90 mm NGM agar plates for each replicate of the screen for worm maintenance.

To prepare NGM follow the steps in the protocol for "Making normal NGM" and pour 35 ml per plate.

#### Protocol

NAME

#### Making normal NGM

CREATED BY

Bonnie Evans

Preview

2 Seed each NGM agar plate with 1 ml of OP50 bacteria, and allow to dry.

### 3 **11 days before tracking**

Chunk worms onto ten seeded 90 mm NGM agar plates.

## Prepare tracking plates

### 4 **11 days before tracking**

Prepare 75 tracking plates (WHAT7701-1651, PK100 : 512-1093) containing NGM agar, prepared as above, with the addition of 1 mM IPTG, 100 µg/ml carbenicillin, 10 µg/ml tetracycline, 1x nystatin (60 Units/ml).

Dispense 200 µl per well using the protocol "Dispensing agar into multiwell plates".



### Protocol

NAME

## Dispensing agar into multiwell plates

CREATED BY

Ida Barlow

[Preview](#)

- 4.1 Dry the tracking plates in a drying oven (with lids removed) until they lose 2.5 to 3.5% of their total weight.

## Bleach synchronising *C. elegans*

### 5 8 days before tracking

Bleach synchronise worms from ten 90 mm plates according to the protocol "Bleach synchronisation of *C. elegans*".

Leave the eggs in M9 buffer in 15 ml falcon tubes on a rotator at 20 degrees for two days.

### Protocol

NAME

## Bleach synchronisation of *C. elegans*

CREATED BY

Ida Barlow

[Preview](#)

## Grow up RNAi library

- 6 We are using a version of the Ahringer RNAi library containing 6198 *C. elegans* genes with orthologues in humans.  
This library is stored as frozen glycerol stocks in seventy three 96 well plates.

### 8 days before tracking

Grow up the RNAi library.



#### Citation

Hernando-Rodríguez, B., Erinjeri, A.P., Rodríguez-Palero, M.J. et al. (2018)  
. Combined flow cytometry and high-throughput image analysis for the study of essential genes in *Caenorhabditis elegans*.  
BMC Biol .

<https://doi.org/10.1186/s12915-018-0496-5>

LINK

#### Citation

Kamath, R., Fraser, A., Dong, Y. et al. (2003)  
. Systematic functional analysis of the *Caenorhabditis elegans* genome using RNAi.  
Nature.

<https://doi.org/10.1038/nature01278>

LINK

- 7 Dispense LB broth (containing 100 µg/ml Carbenicillin, and 10 µg/ml Tetracycline), into seventy four 96 well plates (numbered 1-74) using the Integra VIAFILL.
- 8 Replicate all wells of the RNAi library into corresponding plates of LB (with carbenicillin and tetracycline).
- 8.1 Use a metal 96 pin replicator (Boekel) to transfer glycerol stock from stock plates to the LB culture plates.  
  
To do this remove plates a few at a time from the -80 on dry ice.  
Allow each plate around 1 minute to defrost slightly.  
Firmly press the replicator into the stock plate, making sure to pick up bacteria from each well.  
Transfer replicator to corresponding LB culture plate.
- 8.2 Between replicating each plate, the replicator is dipped in 100% ethanol and flame sterilised. This process is repeated twice to ensure full sterilisation.

One negative control plate is inoculated after double flame sterilisation of the replicator.

#### Safety information

Take caution when flame sterilising the replicator.  
Do not place flaming replicator back into the ethanol bath.  
Work in a clear area - remove any flammable items, papers, tissue, card board boxes etc  
Familiarise yourself with nearest safety blanket.

- 9 Place the LB plates containing the replicated library into large plastic boxes containing moistened paper towel, and incubate for 15 hours at 37 degrees.
- 10 Record OD600 values of 10 of the library plates using a plate reader (Tecan, Spark) to evaluate library growth.  
Visually check the whole library and manually record any wells which have not grown.

### Seeding tracking plates

#### 11 **7 days before tracking**

Seed the tracking plates with the bacterial cultures using the Integra VIAFLO.

Using the Integra VIAFLO with GRIPTIPs mix each well of the library by pipetting 12.5 µl up and down 6 times, before transferring 10 µl from the library to the tracking plates twice.

#### Note

This double transfer is required to pipette a total of 20 µl of bacterial culture to each well of the tracking plates since the max volume of this VIAFLO is 12.5 µl.

- 12 Dry the tracking plates in a biosafety cabinet with lids removed for 3 hours.
- 13 Stack the tracking plates (upside down) in a sealed plastic box containing damp paper towel and incubate at 20 degrees.



## Dispensing *C. elegans* to tracking plates

### 14 **6 days before tracking**

Dispense *C. elegans* onto the tracking plates using the protocol "Dispensing *C. elegans* to 96 well tracking plate using Integra VIAFILL".

We aimed to dispense an average of 4 worms per well, but will have some unavoidable variation between wells.

#### Note

This step should be carried out at the same time for each replicate, such that FUDR can be added 53 hours after dispensing worms.

For example: We dispensed worms at 4pm on Wednesday and add FUDR at 9 am on the following Friday.

#### Protocol

NAME

**Dispensing *C. elegans* to 96 well tracking plate using Integra VIAFILL**

CREATED BY

e.warren Warren

Preview

15 Leave plates to dry for 1 hour in a biosafety cabinet with lids removed.

16 Stack the tracking plates plates (upside down) in a sealed plastic box containing damp paper towel and incubate at 20 degrees.

## Addition of FUDR

17 **4 days before tracking,**



53 hours after dispensing worms, add FUDR to each well of the tracking plates to prevent hatching of progeny.

Dispense 10 µl of 20X FUDR to each well using the Integra VIAFILL with the small 8 channel cassette.

- 18 Leave plates to dry for 1 hour in a biosafety cabinet with lids removed.

#### Note

Check that plates are visibly dry.

- 19 Stack the tracking plates plates (upside down) in a sealed plastic box containing damp paper towel and incubate at 20 degrees.

## Tracking using Hydra rigs

### 20 On the day of tracking

Move tracking plates into the tracking room at 12 pm, one hour before tracking commences to acclimatise to temperature and humidity.

- 21 Record worm behaviour using the syngenta script (5 mins prestim; 6 mins bluelight with 60 sec OFF, [10sec ON, 90sec OFF] x 3 times; 5 mins postsim).

#### 21.1 Run the following script

```
~/scripts$ python run_syngenta_experiment_v2.py -f 291122_rnaifullscreen_run1 -r 01 02 03 04 05
```

-f : date\_file\_run number :

-r : specifies which rigs to start recording on

#### Note

Do not use blank spaces in file name  
Include the run number for each run of tracking

- 22 Move all files to the NAS and then transfer files from the NAS to the group drive on the network for processing.





## Citations

### Step 6

Hernando-Rodríguez, B., Erinjeri, A.P., Rodríguez-Palero, M.J. et al.. Combined flow cytometry and high-throughput image analysis for the study of essential genes in *Caenorhabditis elegans*.

<https://doi.org/10.1186/s12915-018-0496-5>

### Step 6

Kamath, R., Fraser, A., Dong, Y. et al.. Systematic functional analysis of the *Caenorhabditis elegans* genome using RNAi

<https://doi.org/10.1038/nature01278>