

Sep 04, 2025

Version 5

RSVAB WGS and GF protocols V.5

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3xbzqg25/v5>

Maria Iglesias¹

¹Laboratory of Reference and Research in Respiratory Viruses, National Centre for Microbiology, Instituto de Salud Carlos III, 28220 Majadahonda, Spain



Maria Iglesias

National Institute of Health Carlos III

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3xbzqg25/v5>

Protocol Citation: Maria Iglesias 2025. RSVAB WGS and GF protocols. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kqdg3xbzqg25/v5> Version created by [Maria Iglesias](#)

Manuscript citation:

Generic novel system for genomic characterization of Respiratory Syncytial Virus obtaining whole genome sequencing and a full-length G and F sequences.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 04, 2025

Last Modified: September 04, 2025

Protocol Integer ID: 226445

Keywords: novel rsv amplicon, viral genome, cdna from rsv, rsvab wg, viral nucleic acid extract, rsv, specific sequences of the main antigen, systems for genomic characterization, generating cdna, producing amplicon, genomic characterization, main antigen

Abstract

This SOP describes the procedure for generating cDNA from RSV viral nucleic acid extracts and subsequently producing amplicons tiling the viral genome using. We propose two systems for genomic characterization of RSV. First, a novel RSV amplicon-based system for WGS, and second, a method focused on obtaining the specific sequences of the main antigens, G and F.



Sample management

- 1 No prior treatment is necessary. However, in samples with high Cts, the use of DNaseI in the sample as a pretreatment can improve coverage.

dsCDNA generation:

- 2 During this step three master mixes will be prepared: MMI, MMII and MMIII.

Materials: **Kit Superscript III First Strand (Invitrogen)**
100% DMSO
RNAseH (Invitrogen)
Klenow fragment 3' → 5' exo (New England Biolabs)
Primer FR26RV-N : 5'GCC GGA GCT CTG CAG ATA TCNNNNNN 3'

Note

This step must be performed in a RNase free, pre-PCR environment in which post PCR RSV amplicons are not present, to minimise risk of sample contamination.

Citation

Díez-Fuertes F, Iglesias-Caballero M, García-Pérez J, Monzón S, Jiménez P, Varona S, Cuesta I, Zaballos Á, Jiménez M, Checa L, Pozo F, Pérez-Olmeda M, Thomson MM, Alcamí J, Casas I
(2021). A Founder Effect Led Early SARS-CoV-2 Transmission in Spain.

<https://doi.org/10.1128/JVI.01583-20>

LINK

- 3 **MMI Preparation:**

	A	B
	FR26RV-N (10uM)	2
	DMSO	0,5
	Total	2,5 ul

Mix thoroughly by vortexing.

4 **MMII Preparation:**

	A	B
	10x First Strand Buffer	2
	DTT 100 mM	2
	MgCl ₂ 25mM	4
	dNTPs	1
	RNaseOUT	0,5
	SSIII RT	0,5
	Total	10 ul

Kit Superscript III First Strand (Invitrogen)

5 **MMIII Preparation:**

	A	B
	Klenow 5'-3'	1
	RNaseH	0,5
	Total	1,5 ul

6 Defrost extracted RNA.



Maintain **on ice** the MMI,MMII and MMIII mixes.

7 **MMI Amplification:**

Add  5 μL  Sample in MMI mix

Place the tube on a thermocycler and run the following program:

A	B
65°C	5 min
4°C	2 min

Briefly tube centrifugation

8 **MMII Amplification:**

Addition of  10 μL from MMII in the tube with the MMI and the viral extraction.

Place the tube on a thermocycler and run the following program:

A	B
25°C	10 min
50°C	50 min
85°C	10 min
4°C	∞

Briefly tube centrifugation

9 **MMIII Amplification:**

Addition of  1.5 μL of MMIII into the tube with the previous mixes and the viral extraction

Place the tube on a thermocycler and run the following program:

A	B
37°C	60 min

A	B
75°C	15 min

Briefly tube centrifugation

10 **STOP POINT:** cDNA can be stored at 4°C (same day) or -20°C (up to a week).

RSVAB WGS protocol

11

Citation

Iglesias-Caballero M, Camarero-Serrano S, Varona S, Mas V, Calvo C, García ML, García-Costa J, Vázquez-Morón S, Monzón S, Campoy A, Cuesta I, Pozo F, Casas I (2023)

. Genomic characterisation of respiratory syncytial virus: a novel system for whole genome sequencing and full-length G and F gene sequences.

<https://doi.org/10.2807/1560-7917.ES.2023.28.49.2300637>

LINK

Materials:

2x MyTaqRed mix (Bioline)

Primers:

A	B
Primer ID	Sequence (5'-3')
Mix 1	
RSVCombinital	ACGCGAAAAAATGCGTACWACA
RSVWGS1R.2	GCKATTGCAGATCCWACACC

A	B
RSVWGS2F	CACTWACAATATGGGTGCC
RSVWGS2R.2	CRTTYCTTAARGTRGGCC
RSVWGS4R	CATGWTGWYTTATTTGCCCC
RSVWGS3.2F	ACATGGAAAGAYATYAGCC
RSVWGS3.2R	TTGCATCTGTAGCAGGAATGG
RSVWGS9F	GARCAACTCAAAGAAAATGG
RSVWGS9R	AYTGRAACATRGGCACCC
RSVCombending	ACGAGAAAAAAAGTGTCAAAAATAA
Mix 2	
RSVCombinitial	ACGCGAAAAAATGCGTACWACA
RSVWGS1R.2	AYTGRAACATRGGCACCC
RSVWGS5F	CATAATTAYTTTGAATGGC
RSVWGS5R	CAAACATTTAATCTRCTAAGGC
RSVWGS6F	TTATAYAGATATCAYATGGGTGG
RSVWGS6R	CCCTCTCCCAATCTTTTTC
RSVWGS13F	AAAAGATTGGGGAGAGGG
RSV OG1_21	GGGGCAAATGCAACCATGTCC
RSV GF_R	TTCGYGACATATTTGCCCC
RSVCombending	ACGAGAAAAAAAGTGTCAAAAATAA

	A	B	C	D	E	F	G
	1	EPI_ISL_1866820 1	-	13	RSVCombinitial	1	+
	2	EPI_ISL_1866820 1	2193	2212	RSVWGS9F	1	+
	3	EPI_ISL_1866820 1	2321	2340	RSVWGS4R	2	-
	4	EPI_ISL_1866820 1	3337	3355	RSVWGS_2F	1	+
	5	EPI_ISL_1866820 1	3349	3366	RSVWGS9R	2	-
	6	EPI_ISL_1866820 1	4656	4676	OG121	1	+
	7	EPI_ISL_1866820 1	6123	6142	RSVWGS_1R.2	2	-
	8	EPI_ISL_1866820 1	7647	7665	RSVGF-R	2	-
	9	EPI_ISL_1866820 1	7712	7730	RSVWGS_5F	1	+
	10	EPI_ISL_1866820 1	9374	9395	RSVWGS_5R	2	-
	11	EPI_ISL_1866820 1	9358	9376	RSVWGS_3.2F	1	+
	12	EPI_ISL_1866820 1	9986	1000 3	RSVWGS_2R.2	2	-
	13	EPI_ISL_1866820 1	10852	10874	RSVWGS_6F	1	+
	14	EPI_ISL_1866820 1	13076	1309 3	RSVWGS_13F	1	+
	15	EPI_ISL_1866820 1	13090	13109	RSVWGS_6R	2	-
	16	EPI_ISL_1866820 1	14267	1428 7	RSVWGS_3.2R	2	-
	17	EPI_ISL_1866820 1	15200	15225	RSVCombEnding	2	-
	1	EPI_ISL_1653999	1	22	RSVCombinitial	1	+
	2	EPI_ISL_1653999	2159	2178	RSVWGS9F	1	+
	3	EPI_ISL_1653999	2288	2307	RSVWGS4R	2	-

	A	B	C	D	E	F	G
	4	EPI_ISL_1653999	3311	3329	RSVWGS_2F	1	+
	5	EPI_ISL_1653999	3323	3340	RSVWGS9R	2	-
	6	EPI_ISL_1653999	4631	4651	OG121	1	+
	7	EPI_ISL_1653999	6102	6121	RSVWGS_1R.2	2	-
	8	EPI_ISL_1653999	7618	7636	RSVGF-R	2	-
	9	EPI_ISL_1653999	7699	7717	RSVWGS_5F	1	+
	10	EPI_ISL_1653999	9344	9365	RSVWGS_5R	2	-
	11	EPI_ISL_1653999	9328	9346	RSVWGS_3.2F	1	+
	12	EPI_ISL_1653999	9956	9973	RSVWGS_2R.2	2	-
	13	EPI_ISL_1653999	10822	10844	RSVWGS_6F	1	+
	14	EPI_ISL_1653999	13062	13079	RSVWGS_13F	1	+
	15	EPI_ISL_1653999	13060	13059	RSVWGS_6R	2	-
	16	EPI_ISL_1653999	14237	14257	RSVWGS_3.2R	2	-
	17	EPI_ISL_1653999	15209	15222	RSVCombEnding	2	-

Primer scheme with RSA and RSB RefSeq

Note

The protocol is based in the RSV genome amplification in two separate mixes with an unique amplification program. The mixes that will be mixed at the end of cycling.

12 ..2Preparation of RSV Amplification Mix 1:

	A	B
	MyTaq Red 2x	15
	H2O	8,4



A	B
RSV Combinital (5 uM)	0,2
RSVWGS1R2.2 (5uM)	0,2
RSVWGS2F (5 uM)	0,2
RSVWGSW2R.2 (5 uM)	0,2
RSVWGS4R (5 uM)	0,2
RSVWGS3.2F (5 uM)	0,2
RSVWGS3.2R (5 uM)	0,2
RSVWGS9F (5uM)	0,2
RSVWGS9R (5 uM)	0,2
RSV Combending (5uM)	0,2
Total	25

13 Preparation of RSV Amplification Mix 2:

A	B
2x My Taq Red	15
H2O	7,6
RSV Combinital (5uM)	0,2
RSVWGS1R.2 (5 uM)	0,2
RSVWGS5F (5 uM)	0,2
RSVWGS5R (5 uM)	0,2
RSVWGS6F (5 uM)	0,2
RSVWGS6R (5 uM)	0,2
RSVWGS9F (5 uM)	0,2
RSVWGS9R (5 uM)	0,2
RSVWGS13F (5 uM)	0,2
RSV OG1_21 (5 uM)	0,2
RSV GF_R (5 uM)	0,2


A	B
RSV Combending (5 uM)	0,2
Total	25 ul

14 Addition of  5 μ L of the previous prepared double stranded cDNA on each mix.

15 Amplification protocol:

A	B	C
95°C	1 min	
95°C	30 seg	x45
55°C	8 min	
72°C	2 min	
72°C	5 min	
12°C	∞	

16 To assess PCR performance, the amplicons can be loaded onto a 1% agarose gel for electrophoresis.



1-4: <22
5-8: 22-28
9-12: >30



- 17 Finally, mix in one single tube both mixes and proceed to purification and library preparation.

RSVAB GF protocol starting from ds cDNA

- 18 Due to the significance of achieving accurate RSV genomic characterization, it was developed the RSVAB-GF PCR to complement the genomic coverage of both antigenic major proteins in cases where WGS encounters difficulties, and to provide a simpler and more cost-effective method of obtaining the sequences of both antigens.

- 19 Materials:

2x MyTaqRed mix (Bioline)

Primers:

A	B
RSVWGS_2F	CACTWACAATATGGGTGCC
RSVGF-R	TTCGYGACATATTTGCCCC

- 20 **Preparation of cDNA GF amplification mix:**

A	B
H2O	5,5
2X MyTaqRed	12,5
RSVWGS_2F	1
RSVGF-R (10 uM)	1
Total	20 ul

- 21 Addition of  5 µL of the previous prepared double stranded cDNA on the mix.

- 22 **Amplification protocol cDNA GF:**



	A	B	C
	95°C	1 min	
	95°C	30 seg	x35
	60°C	3 min	
	72°C	2 min	
	72°C	5 min	
	12 °C	∞	

RSVAB GF protocol starting from viral extraction

23 Materials:
 Qiagen OneStep RT-PCR kit.

24 **Preparation of GF amplification mix:**

A	B
H2O	10
5xQ PCR MM	6
dNTPs	1
RSVWGS_2F (10uM)	1
RSVGF-R (10 uM)	1
RT-PCR mix	1
Total	20 ul

25 Addition of  10 µL of the viral extraction

26 **Amplification protocol GF:**



	A	B	C
	48°C	60 min	
	95°C	15 min	
	95°C	30 seg	x 45
	55°C	2 min	
	72°C	1 min	
	72°C	5 min	
	12°C	∞	

Citations

Step 11

Iglesias-Caballero M, Camarero-Serrano S, Varona S, Mas V, Calvo C, García ML, García-Costa J, Vázquez-Morón S, Monzón S, Campoy A, Cuesta I, Pozo F, Casas I. Genomic characterisation of respiratory syncytial virus: a novel system for whole genome sequencing and full-length G and F gene sequences.

<https://doi.org/10.2807/1560-7917.ES.2023.28.49.2300637>

Step 2

Díez-Fuertes F, Iglesias-Caballero M, García-Pérez J, Monzón S, Jiménez P, Varona S, Cuesta I, Zaballos Á, Jiménez M, Checa L, Pozo F, Pérez-Olmeda M, Thomson MM, Alcamí J, Casas I. A Founder Effect Led Early SARS-CoV-2 Transmission in Spain.

<https://doi.org/10.1128/JVI.01583-20>