

May 07, 2020

Sampling and viral concentration for SARS-CoV-2 detection in wastewater

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.bf5mjg46>

Veronica Antelo¹, Gregorio Iraola¹

¹Microbial Genomics Laboratory, Institut Pasteur Montevideo, Uruguay

Coronavirus Method De...



Iraola Lab

Institut Pasteur Montevideo

IRAOLA
LAB

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bf5mjg46>

Protocol Citation: Veronica Antelo, Gregorio Iraola 2020. Sampling and viral concentration for SARS-CoV-2 detection in wastewater . protocols.io <https://dx.doi.org/10.17504/protocols.io.bf5mjg46>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol.

Created: May 07, 2020

Last Modified: May 07, 2020

Protocol Integer ID: 36749

Keywords: viral concentration for sar, viral concentration, sar

Guidelines

This protocol has been effective at recovering enveloped viruses genes for qPCR detection, but not appropriate when infective viruses are desired.

Materials

STEP MATERIALS

✖ 200 mL to 1 L plastic bottles with screw cap

✖ 0.25 N glycine buffer (pH 9.5)

✖ 2× phosphate-buffered saline (PBS, pH 7.2)

Protocol materials

✖ 200 mL to 1 L plastic bottles with screw cap

✖ 0.25 N glycine buffer (pH 9.5)

✖ 2× phosphate-buffered saline (PBS, pH 7.2)

Sampling wastewater

1

Autoclave plastic containers and expose to UV for 10 min.

200 mL to 1 L plastic bottles with screw cap

2

Fill bottles with wastewater at sampling points using sterile gloves.

Note

Sampling points may be diverse and accessibility to water sources will be sometimes difficult depending on site infrastructure. Consider these aspects before starting sampling.

Note

If site is inaccessible by hand, use a rope tied to the bottle to reach narrow or deep places. When recovering the bottle, immediately tap it and sanitize the outside using EtOH 70%.

3

Place bottles in a refrigerated container (i.e. with ice flakes or gel packs) until returning to the lab. Our sampling days usually take 4-6 hours since first sample is collected.

Concentration by ultracentrifugation

1h

4

Prepare your reagents and materials:

Equipment

Ultracentrifuge

NAME

Beckman Coulter

BRAND

XPN-90

SKU

Your ultracentrifuge must reach 100,000 x g and cool down to 4° C SPECIFICATIONS



Equipment

Centrifuge

NAME

Bench centrifuge

TYPE

Eppendorf

BRAND

5424

SKU

⊗ 0.25 N glycine buffer (pH 9.5)

⊗ 2× phosphate-buffered saline (PBS, pH 7.2)

- 5 Place 42 mL of sewage water in an appropriate ultracentrifugation tube and ultracentrifugate as follows:

1h

⦿ 100000 x g

🌡 4 °C

⌚ 01:00:00

- 6 Discard supernatant resuspend pelleted viral particles as follows:

- 6.1 Discard supernatant.

- 6.2 Resuspend viral particles in  3.5 mL glycine buffer .

5m

- 6.3 Incubate on ice for  00:30:00

30m



7 Add  3.5 mL PBS and gently mix.

5m

8 Clarify supernatant by centrifugation as follows:

15m

 12000 x g

 00:15:00

9 Finally recover viruses by ultracentrifugation as follows:

1h

 100000 x g

 4 °C

 01:00:00

10 Resuspend viruses in  200 µL PBS .

11 Process immediately for nucleic acid extraction or store at  -80 °C .