


Mar 10, 2023

🌐 SARS-CoV-2 RNA extraction with Ceres Nanotrap and Zymo Environ Water



DOI

<https://dx.doi.org/10.17504/protocols.io.14egn26qqg5d/v1>

Amanda Windsor¹, kathryn.judy¹, Tamara Walsky¹, Padmini Ramachandran¹, Chris Grim¹, Maria Hoffmann¹

¹FDA CFSAN

GenomeTrakr

Tech. support email: genomeTrakr@fda.hhs.gov



Amanda Windsor

FDA

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn26qqg5d/v1>

Protocol Citation: Amanda Windsor, kathryn.judy , Tamara Walsky, Padmini Ramachandran, Chris Grim, Maria Hoffmann 2023. SARS-CoV-2 RNA extraction with Ceres Nanotrap and Zymo Environ Water . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn26qqg5d/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: December 06, 2022

Last Modified: March 10, 2023

Protocol Integer ID: 73627

Keywords: SARS-CoV-2, Ceres Nanotrap, Zymo Environ Water, RNA, wastewater, rna extraction, ceres nanotrap, based virus capture, virus capture, zymo environ water, zymo dnase step, wastewater, extraction, dnase, rna, virus

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **[protocols.io](#)** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **[protocols.io](#)**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol uses the Ceres Nanotrap® particle-based virus capture and concentration method for 10mL of wastewater followed by extraction with the Zymo Environ™ Water RNA extraction kit with a Zymo DNase step.

Guidelines

When developing this protocol, we extracted at least 2 replicates of each wastewater sample to ensure we had ample RNA for downstream processes

RNA extraction is performed at room temperature and centrifugation at 10,000-16,000 x g for 30s



Materials

1. Wastewater sample

Reagents and Kits

1. Ceres Nanosciences Nanotrap[®] Magnetic Viral Particles (Ceres Nanosciences: SKU 44202)
2. Ceres Nanosciences Nanotrap[®] Enhancement Reagent 2 (Ceres Nanosciences: SKU 10112)
3. Magnetic separator for 15 mL conical tubes, such as Invitrogen[™] DynaMag[™]-2 Magnet (ThermoFisher Cat# 12-301-D)
4. Magnetic separator for 2mL micro centrifuge tubes, such as Invitrogen[™] DynaMag[™]-2 Magnet (ThermoFisher Cat# 12-321-D)
5. Zymo Environ[™] Water RNA Kit (Zymo Research: R2042)
6. Zymo DNA/RNA Shield[™] (Zymo Research: R1100-50 or R1100-250)
7. Zymo DNase Set 1 (Zymo Research: E1010)

Equipment

1. Programable Heat Block
2. Mini vortex mixer
3. Mini Centrifuge (Max capable of 16,000 x g & fits 1.5/2mL tubes)
4. tube rotator (e.g. Fisherbrand Mini Tube Rotator Cat 88-861-05 or similar)

Consumables

1. 100% absolute ethanol
2. DNase/RNase Free Water
3. 15mL conical tubes
4. 1.5 or 2mL microcentrifuge tubes
5. 100-1000uL pipette
6. 20-200uL pipette
7. 100-1000uL filtered pipette tips
8. 20-200uL filtered pipette tips
9. serological pipetting aid
10. 10mL serological pipettes
11. 5mL serological pipettes
12. 1mL serological pipettes

Safety warnings

 Please refer to Safety Data Sheets (SDS) for health and environmental hazards.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.



Before start

- Store the reagents separately from RNA/TNA (total nucleic acid) samples.
- Use a clean designated work area and separate pipettes for pre- and post-extraction steps to minimize the potential for cross-contamination
- Wear a lab coat and protective eyewear.
- Wear gloves and change them often.
- Prevent contamination by using aerosol-resistant pipette tips.



Before you start

- 1 Turn on heat block  95 °C
- 2 Ensure appropriate volume of DNase I is available (5µL per sample), or make new aliquots
- 2.1 Add  275 µL DNase/RNase-Free water to reconstitute lyophilized DNase I (1U/µL)

Note

Aliquot reconstituted DNase I in volumes appropriate for your lab's throughput (e.g., 15-20µL aliquots in 0.5mL microcentrifuge tubes) to avoid multiple freeze/thaw cycles

Viral Capture with Nanotrap® Particles

30m

- 3 Shake wastewater bottle to mix then let sit  00:00:45 45s
- 4 Using a 10mL serological pipette, carefully pipette  10 mL of wastewater into a 15mL conical tube

Note

Input volumes of 20, 30, and 40mL have also been tested. See note in step 7 for volume of nanotrap particles to add based on starting volume

- 5 Add  100 µL of Nanotrap® Enhancement Reagent 2 (ER2) and invert 15mL tube 2-3 times to mix
- 6 Re-suspend Nanotrap® particles by inverting the bottle 5 times
- 7 Add  150 µL Nanotrap® particles to the sample

**Note**

the volume of nanoparticles for different starting volumes of wastewater are as follows:

	Input Wastewater (mL)	Nanotrap Particles (μL)
	20	300
	30	450
	40	600

volume of nanoparticles for different starting volumes of wastewater

- 8 Incubate samples  Room temperature  00:10:00 with constant rotation 10m
- 9 Place samples on magnetic rack to separate Magnetic Nanotrap® particles from the sample - at least  00:02:00 2m
- 10 After beads have settled, use a 5mL serological pipette to remove all of the supernatant without disturbing the pelleted beads
- 11 Add  1 mL of DNase/RNase Free water to the tube
- 11.1 Remove tube from magnet and re-suspend the pelleted beads using a 100-1000uL pipette
- 11.2 Transfer suspended beads to a 1.5mL microcentrifuge tube
- 12 Place microcentrifuge tube on the 2mL tube-compatible magnetic rack
- 12.1 Incubate until the beads have settled - at least  00:02:00 2m



13 Remove supernatant with a 100-1000uL pipette without disturbing the pellet. Remove any small amount of remaining supernatant with a smaller pipette tip (e.g. 2-20uL pipette)

14 Remove the tubes from the magnet and re-suspend the pellet with  375 μ L Zymo DNA/RNA Shield and  125 μ L Zymo DNase/RNase-Free water from the Zymo Environ Water RNA Kit

15 Incubate the samples at  95 °C for  00:05:00

5m

15.1 While samples are incubating, add  400 μ L Zymo RNA Binding Buffer from the Zymo Environ Water RNA Kit to new 1.5mL tubes (one per sample)

16 Remove tubes from heat block and place on a magnetic rack and allow beads to settle until supernatant is clear - at least  00:02:00

2m

Note

Collect any liquid from caps by brief centrifugation prior to placing the tubes on the magnetic rack

16.1 Reset heat block temperature to  27 °C

17 Transfer  400 μ L of supernatant to the corresponding tube prepared in step 15.1 and mix by gentle pipetting

Note

sample tubes volume = 800 μ L

Zymo Environ™ Water RNA Kit

5m



- 18 Transfer entire sample to a Zymo-Spin™ IIICG Column in a clean collection tube and centrifuge  10000 x g, Room temperature, 00:00:30 and **keep the flow-through**

30s

Note

Label both the spin column and collection tube

- 19 Add 800uL of ethanol (95-100%) to the flow-through in the collection tube from step 18 and mix well by gentle pipetting

Note

sample volume = 1600µL

- 20 Transfer  800 µL into a **new** Zymo-Spin™ IIICG Column in a clean collection tube and centrifuge  10000 x g, Room temperature, 00:00:30 and **discard the flow through**

30s

- 20.1 Repeat step 19 with the remaining  800 µL of sample using the same collection tube

- 21 Add  400 µL of RNA Prep Buffer to the column and centrifuge  10000 x g, Room temperature, 00:00:30 and **discard the flow-through**

30s

- 21.1 Transfer column to an RNase-Free 1.5mL microcentrifuge tube

- 22 Add  100 µL Zymo DNase/RNase-Free water directly to the column matrix and centrifuge  10000 x g, Room temperature, 00:00:30 and **keep the flow-through for step 24**

30s

- 23 Place a Zymo-Spin III-HRC Filter into a new collection tube and add  600 µL Prep Solution Centrifuge  8000 x g, Room temperature, 00:03:00 **discard the flow-through**

3m

- 23.1 transfer the column to an RNase-Free 1.5mL microcentrifuge tube



- 24 Transfer the eluted RNA from step 22 into the Zymo-Spin III-HRC filter prepared in step 23 and  16000 x g, Room temperature, 00:03:00 **keep the flow-through** 3m
- 25 Add  200 μL RNA Binding Buffer to the **filtrate** and mix well by gently pipetting up and down.
- 25.1 Add  300 μL of ethanol (95-100%) to the filtrate + RNA Binding Buffer and mix well by gently pipetting up and down.
- 26 Transfer the mixture into a Zymo-Spin IC column in a collection tube and  10000 x g, Room temperature, 00:00:30 **discard the flow-through** 30s
- 27 Add  400 μL RNA Wash Buffer to column and  10000 x g, Room temperature, 00:00:30 **Discard flow-through** 30s
- 28 Add  5 μL DNase I and  35 μL DNA Digestion Buffer to the column matrix
- 28.1 Incubate  27 °C  00:20:00 20m
- 29 Add  400 μL of RNA Prep Buffer to the column and  10000 x g, Room temperature, 00:00:30 **discard the flow-through.** 30s
- 30 Add  700 μL of RNA Wash Buffer to the column and  10000 x g, Room temperature, 00:00:30 **discard the flow-through.** 30s
- 31 Add  400 μL of RNA Wash Buffer to the column and  10000 x g, Room temperature, 00:02:00 to ensure complete removal of the wash buffer. Transfer the column carefully into an RNase-free tube. 2m
- 32 Add  50 μL of DNase/RNase-Free Water directly to the column matrix and  10000 x g, Room temperature, 00:00:30 The eluted RNA can be used immediately or stored at  -70 °C . 30s

