

Jul 24, 2025

Screening protocol for colistin resistance using CHROMagar™ COL-APSE medium

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlky4o1g5r/v1>



Soe Yu Naing¹

¹Utrecht University

Soe Yu Naing: Project Manager | PhD Candidate



Soe Yu Naing

Utrecht University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlky4o1g5r/v1>

External link: <https://youtu.be/KC1x9u4d02g>

Protocol Citation: Soe Yu Naing 2025. Screening protocol for colistin resistance using CHROMagar™ COL-APSE medium. protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlky4o1g5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 23, 2025

Last Modified: July 24, 2025

Protocol Integer ID: 223073

Keywords: LMIC, Antimicrobial resistance, Colistin, Screening, COINCIDE, Microbiology, Affordable, Easy, CHROMagar, CHROMagar™ COL-APSE, protocol for colistin resistance, colistin resistance, apse medium colistin resistance, antimicrobial resistance monitoring, chromagar screening, scalable tool for antimicrobial resistance monitoring, colistin, resistant isolate, differentiation between colistin, resistant gram, screening protocol, selective chromagar col, large numbers of isolate, practical for routine surveillance effort, negative bacteria, routine surveillance effort, isolate

Funders Acknowledgements:

JPIAMR

Grant ID: COINCIDE

ZonMw

Grant ID: COINCIDE

International Development Research Centre (IDRC)

Grant ID: COINCIDE

the Research Foundation – Flanders (FWO)

Grant ID: COINCIDE

Canadian Institutes of Health Research (CIHR)

Grant ID: COINCIDE

Abstract

Colistin resistance presents a major challenge in combating multidrug-resistant Gram-negative bacteria. To address this concern, a screening protocol employing the selective CHROMagar COL-APSE medium has been developed. This medium enables differentiation between colistin-susceptible and colistin-resistant *E. coli* strains based on growth characteristics and distinct colony coloration. The method is well-suited for screening large numbers of isolates, making it practical for routine surveillance efforts. Its applicability across both low- and middle-income countries (LMICs) and high-income countries (HICs) enhances its global relevance, offering a cost-effective and scalable tool for antimicrobial resistance monitoring. Additionally, a tailored spot-test protocol was developed to complement the CHROMagar screening, providing a rapid and adaptable approach for preliminary identification of resistant isolates.



Materials

- CHROMagarTM COL-APSE medium plates prepared according to the manufacturer's protocol (Appendix: Manufacture leaflet) Plates can be stored for up to 1 month under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration
- fresh 24-hour cultures of the E.coli isolates grown on blood agar
- MCR-1 control strain for quality control purposes.
- Tubes with sterile saline NaCl (0.85%) to prepare 0.5 McFarland suspension
- Vortex
- Densitometer for measuring bacterial density
- Sterile wooden-cotton swab



1. Background

- 1 Colistin resistance has become a significant concern in the fight against multidrug-resistant Gram-negative bacteria. To address this issue, a screening protocol utilizing CHROMagar™ COL-APSE medium has been developed. This selective agar medium allows for the differentiation of colistin-susceptible and colistin-resistant strains based on growth and distinctive colony colors.

2. Aim

- 2 To efficiently screen a large batch of *E. coli* isolates obtained from fecal samples for colistin resistance using CHROMagar™ COL-APSE medium.

3. Materials needed

- 3
 - CHROMagar™ COL-APSE medium plates prepared according to the manufacturer's protocol
Plates can be stored for up to 1 month under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration
 - fresh 24-hour cultures of the *E. coli* isolates grown on blood agar
 - MCR-1 control strain for quality control purposes.
 - Tubes with sterile saline NaCl (0.85%) to prepare 0.5 McFarland suspension
 - Vortex
 - Densitometer for measuring bacterial density
 - Sterile wooden-cotton swab

4. Procedures

- 4 Note: If the agar plate has been refrigerated, allow to warm to room temperature before inoculation
 - 4.1 Work with batches of strains as a total of 8 strains can be tested per plate
 - 4.2 Include 1 positive control and 1 negative control strain on a plate on the day of testing regardless of the number of test strains
 - 4.3 Make a pie of 8 on CHROMagar plates and label them carefully



- 4.4 Prepare 0.5 McFarland bacterial suspension from fresh 24-hour cultures of the *E.coli* isolates, grown on blood agar
- 4.5 Use a wood-cotton swab to immerse into a 0.5 McFarland bacterial suspension, and make sure to remove excess liquid by pressing and turning the swab against the inside of the tube
- 4.6 Use the wood-cotton swab to spot inoculate the bacterial suspension on CHROMagar plates (8 grid as a pie)
- 4.7 Incubate in aerobic conditions at 35-37 °C for 18-24 hours
- 4.8 On the next day, check the results and take a photograph of the plates
- 4.9 Record positive and negatives strains and include the results of the control strains
- 4.10 If the results of the control strains are not correct, all strains of that day has to be repeated
- 4.11 In case of any doubt or possible contamination, retest the strain(s)

The scheme for colistin resistance screening using selective CHROMagar COL-APSE medium to detect *E. coli* isolates resistant to colistin

5





Figure 1. Screening protocol using a selective CHROMagar COL-APSE medium to detect *E. coli* isolates resistant to colistin.

5. Interpretation of colistin resistance

- 6 Growth of *E. coli* with dark pink to reddish color (as in figure below) is considered positive. Absence of growth is considered negative.



Example of colony appearance of colistin-resistant *E. coli*

6. Quality Control

- 7 Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the following ATCC strain: *MCR-1 E. coli* NCTC 13846 displaying typical colony appearance of dark pink to reddish.

Please perform the calibration of densitometer regularly.

7. Video Protocol

- 8 [COINCIDE VideoProtocol ColistinScreening Chromagar_Naing](#)

Protocol references

Abdul Momin MHF, Bean DC, Hendriksen RS, Haenni M, Phee LM, Wareham DW. CHROMagar COL-APSE: a selective bacterial culture medium for the isolation and differentiation of colistin-resistant Gram-negative pathogens. *J Med Microbiol.* 2017;66(11):1554-1561. doi:10.1099/jmm.0.000602



Acknowledgements

The COINCIDE study is funded by the Joint Programming Initiative on Antimicrobial Resistance (JPIAMR)'s 2021 call on "One Health interventions to prevent or reduce the development and transmission of AMR. The participating countries are supported by the following agencies: the International Development Research Centre (IDRC) for Indonesia, the Netherlands Organisation for Health Research and Development (ZonMw) for the Netherlands, the Research Foundation – Flanders (FWO) for Belgium, and the Canadian Institutes of Health Research (CIHR) for Canada.