


Jul 14, 2023

Sequencing of Canine Parvovirus (CPV) from Rectal Swab Samples V.1

DOI

dx.doi.org/10.17504/protocols.io.36wgq3563lk5/v1



Sara França de Araújo dos Santos¹, Ueric José Borges de Souza¹, Martha Trindade Oliveira², Jairo Jaime³, Fernando Rosado Spilki⁴, Ana Cláudia Franco², Paulo Michel Roehe², Fabrício Souza Campos²

¹- Bioinformatics and Biotechnology Laboratory, Campus of Gurupi, Federal University of Tocantins, Gurupi 77410-570, Brazil;

²Virology Laboratory, Department of Microbiology, Immunology, and Parasitology, Institute of Basic Health Sciences, Federal University of Rio Grande do Sul, Porto Alegre 90035-003, Brazil;

³Universidad Nacional de Colombia, Sede Bogotá. Facultad de Medicina Veterinaria y de Zootecnia, Departamento de Salud Animal. Centro de Investigación en Infectología e Inmunología Veterinaria (CI3V). Carrera 30 # 45-03, Bogotá D.C. CP 11132. Colombia;

⁴Molecular Microbiology Laboratory, Feevale University, Novo Hamburgo 93525-075, Brazil



Martha Oliveira

UFRGS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgq3563lk5/v1>

Protocol Citation: Sara França de Araújo dos Santos, Ueric José Borges de Souza, Martha Trindade Oliveira, Jairo Jaime, Fernando Rosado Spilki, Ana Cláudia Franco, Paulo Michel Roehe, Fabrício Souza Campos 2023. Sequencing of Canine Parvovirus (CPV) from Rectal Swab Samples V.1. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgq3563lk5/v1>

Manuscript citation:

Santos SFA, Souza UJB, Oliveira MT, Jaime J, Spilki RS, Franco AC, Roehe PM, Campos FS. Recovery of complete genomes of canine parvovirus from clinical samples. bioRxiv preprint, 25p. 2023. <https://www.biorxiv.org/content/10.1101/2023.07.12.548703v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 30, 2023

Last Modified: July 14, 2023

Protocol Integer ID: 84318

Keywords: CPV, Canine Parvovirus, sequencing, viral DNA, rectal swab, DNA extraction, canine parvovirus, viral genome, full length of the viral genome, pathogenicity for puppy, sequencing, genome assembly, genome, reliable coverage of the genome, viral enteritis, puppies with clinical signal, contagious viral disease, multiplex pcr, cpv, pcr, pathogenicity, larger numbers of read, resulting read

Funders Acknowledgements:

This work was partially supported by grants from Rede Corona-ômica BR MCTI/FINEP (<http://www.corona-omica.br-mcti.Incc.br>, accessed on 25 June 2023) affiliated by RedeVírus/MCTI

Grant ID: FINEP 01.20.0029.000462/20; CNPq 404096/2020-4

This work was partially supported by grant from the National Council for Scientific and Technological Development–CNPq

Grant ID: Process number: 443215/2019-7

Abstract

Canine parvovirus (CPV) is a highly contagious viral disease that affects dogs, especially puppies. CPV-2 is recognized for its resilience in contaminated environments, ease of transmission among dogs, and pathogenicity for puppies. In this protocol, we have adapted the methodology to allow the recovery of complete CPV-2 genomes directly from clinical samples (dry swabs) from puppies with clinical signals of viral enteritis. A multiplex PCR was designed with primers targeting fragments of 400 to 1,000 base pairs (bp) along the full length of the viral genome. The resulting reads were compared after sequencing with the Nanopore technology. Genome assembly revealed that the smaller fragments generated larger numbers of reads, allowing a more reliable coverage of the genome than those attained with primers targeting larger amplicons. Both new methodologies were efficient in amplification and sequencing.

Materials

Consumables:

Isopropanol

Beta-mercaptoethanol

Ethanol

Nuclease-free water

(primers from Table 1 and 2)

R9.4 Oxford MinION flowcell (FLO-MIN106)(Oxford Nanopore)

Commercial Kits:

- * Quick-DNA/RNA™ Viral MagBead (Zymo Research)
- * Q5® High-Fidelity 2X Master Mix (New England Biolabs)
- * AMPure XP beads (Beckman Coulter)
- * Ligation Sequencing kit SQK-LSK-109 (Oxford Nanopore)
- * Native Barcoding kits EXP-NBD104 and EXP-NBD114 (Oxford Nanopore)

Equipment:

Thermocycler

Magnetic stand (or automated extractor with magnetic block)

MinION Mk1B device (Oxford Nanopore)

Protocol materials

 Quick-DNA/RNA Viral MagBead **Zymo Research Catalog #R2140 / R2141**

 Q5 High-Fidelity PCR Kit - 200 rxns **New England Biolabs Catalog #E0555L**

 Ampure XP beads **Beckman Coulter Catalog #A63881**

 Qubit dsDNA HS Assay kit **Thermo Fisher Scientific Catalog #Q32854**

Troubleshooting

Nucleic Acid Extraction using Quick-DNA/RNA Viral MagBead (Zymo Research)

1 This protocol uses

 Quick-DNA/RNA Viral MagBead Zymo Research Catalog #R2140 / R2141 for extraction.

All kit's solutions need to be prepared beforehand, as described by the kit. If you've already prepared the solution, go to step 2. **This is a modified version of the kit's protocol.** The original kit's protocol can be found [here](#).

1.1 DNA/RNA Buffer:

Add  500 μ L of **beta-mercaptoethanol** (user supplied) per 100 ml Viral DNA/RNA **Buffer**, (final 0.5% (v/v))

1.2 DNA/RNA Wash 1:

Kit R2140 → Add  20 mL of **isopropanol** (2-propanol PA, user supplied) to MagBead DNA/RNA **Wash 1** concentrate

Kit R2141 → Add  80 mL of **isopropanol** (2-propanol PA, user supplied) to MagBead DNA/RNA **Wash 1** concentrate

1.3 DNA/RNA Wash 2:

Kit R2140 → Add  30 mL of **isopropanol** (2-propanol PA, user supplied) to MagBead DNA/RNA **Wash 2** concentrate

Kit R2141 → Add  120 mL of **isopropanol** (2-propanol PA, user supplied) to MagBead DNA/RNA **Wash 2** concentrate

1.4 DNA/RNA Shield™:

Kit R2140 → Add  25 mL of **nuclease-free water** (user supplied) to 2X DNA/RNA Shield™ concentrate

Kit R2141 → Add  125 mL of **nuclease-free water** (user supplied) to 2X DNA/RNA Shield™ concentrate

1.5 Proteinase K (20 mg/ml):

Kit R2140 → Add  1.04 mL of **Proteinase K Storage Buffer** to lyophilized Proteinase K ( 20 mg)

Kit R2141 → Add  3.12 mL of **Proteinase K Storage Buffer** to lyophilized Proteinase K ( 60 mg)

Mix by vortexing. Use immediately or store frozen aliquots.

2 Elute each dry  Sample into  400 µL of 1X viral DNA/RNA **buffer**.

2.1 Our samples were dry rectal swabs of animals confirmed positive by a rapid test. Samples were collected and stored in a centrifuge tube (without transporting media) at -20 °C until processing.

3 For each sample, add  10 µL of beads and  4 µL of proteinase K.

Obs: when working with many samples, one can always prepare a mixed solution with both components. Just make to keep the beads suspended while pipetting.

3.1 Transfer the plate (or your tubes) to **magnetic stand** (user supplied). After the beads have pelleted, aspirate and discard the supernatant.

This protocol can be **fast tracked** by the use of an **automated extractor**. We use EXTRACTA 96 (**Loccus**).

4 Perform one wash of the pellet (beads) with  250 µL of MagBead DNA/RNA **Wash 1** solution.



- 5 Perform one wash of the pellet (beads) with  250 μL of MagBead DNA/RNA **Wash 2** solution.
- 6 Add  250 μL of **80% ethanol** and mix well.
Transfer the sample (liquid and beads) to a new place/tube.
- 7 Pellet the beads and discard the the supernatant.
- 8 Add  50 μL of **nuclease-free water** and mix well (this is the elution step).
- 9 Pellet the beads (magnetically) and transfer the supernatant (containing DNA/RNA) to a new plate/tube.

Multiplex PCR to obtain CPV genome fragments

- 10 This protocol uses 2 sets of oligos: one set amplifies ~400 bp sequences that overlap about 100 nucleotides with each other (Table 1); the other amplifies ~1000 bp sequences that overlap about 100 nucleotides with each other (Table 2). Primers are listed in the tables below and were design based on Canine parvovirus reference sequence **NC_001539.1** (GenBank).

Primers should be mix into 2 pools for each set. To prepare a pool, add  5 μL of each primer (at  100 micromolar (μM)) as indicated. This will be your 10X (stock) solution. Working solution should be dilute to 1X concentration ( 10 micromolar (μM)).

Name	Pool	Sequence (5' – 3')	Size (nt)	%GC	Tm
CPV-400_0.7_LEFT*	1	ATGTCTGGCAACCAGTATACTG	22	45.50	60.30
CPV-400_1_LEFT	1	AACCAACTGACCAAGTTCACG T	22	45.45	60.80



Name	Pool	Sequence (5' – 3')	Size (nt)	%GC	Tm
CPV-400_1_RIGHT	1	GTTCCAGCGAACATCCTTTCCA	22	50.0 0	61.3 1
CPV-400_1.5_LEFT*	1	AGGTGGCGGGCTAATTGTG	19	57.9 0	59.5 0
CPV-400_2_LEFT	2	AAGAAACATGCAGAAAATGAAG CATT	26	30.7 7	59.7 3
CPV-400_2_RIGHT	2	CGTAGCCATTTACCAGTTGCTT G	23	47.8 3	60.6 7
CPV-400_3_LEFT	1	ATGGGGAAAAGATCAAGGCTG G	22	50.0 0	60.8 1
CPV-400_3_RIGHT	1	AGTGTGCTGACAATTTGTCTGT C	23	43.4 8	59.9 4
CPV-400_4_LEFT	2	GGGTGACTATATTAACATACAG ACATAAGC	30	36.6 7	60.5 9
CPV-400_4_RIGHT	2	TCCTGGTTGTGCCATCATTTC	22	45.4 5	60.6 8
CPV-400_5_LEFT	1	ACTTTGCGGGACTTGGTTAGTA	22	45.4 5	59.8 1
CPV-400_5_RIGHT	1	ACAACCAACATTACCCACAGCT	22	45.4 5	60.6 1
CPV-400_6_LEFT	2	CAGTTCTTTTTTCATGGACCAGC A	23	43.4 8	59.9 3
CPV-400_6_RIGHT	2	AAACCAAAGTCTCCTGGAAGC T	22	45.4 5	60.0 1
CPV-400_7_LEFT	1	TGGATGTGAAGAAAGACCTGAA CA	24	41.67	60.4 7
CPV-400_7_RIGHT	1	AACGCCAAGTTGGTTTGATTGT	22	40.9 1	60.0 1



Name	Pool	Sequence (5' – 3')	Size (nt)	%GC	Tm
CPV-400_8_LEFT	2	AGTGGACCTTGCACTGGAAC	20	55.00	60.20
CPV-400_8_RIGHT	2	GCTTCGTCGTGTTCTTTTGCAG	22	50.00	61.33
CPV-400_9_LEFT	1	AAATATCTTGGGCCTGGGAACA	22	45.45	59.87
CPV-400_9_RIGHT	1	ACTGCTCCATCACTCATTGGTG	22	50.00	60.80
CPV-400_10_LEFT	2	ACCACCTCATATTTTCATCAATCTTGC	27	37.04	60.91
CPV-400_10_RIGHT	2	TCAACCAATGACCAAGGTGTTACA	24	41.67	60.95
CPV-400_11_LEFT	1	GTGGTTGTAAATAATATGGATAA AACTGCA	30	30.00	60.00
CPV-400_11_RIGHT	1	TGTTCTATCCCATTGAAAATAA TATCTCCA	30	30.00	59.80
CPV-400_12_LEFT	2	TGCCATTTACTCCAGCAGCTAT	22	45.45	60.01
CPV-400_12_RIGHT	2	TCCCATTTGAGTTACACCACGT	22	45.45	60.08
CPV-400_13_LEFT	1	TTTGCCTCAATCTGAAGGAGCT	22	45.45	60.14
CPV-400_13_RIGHT	1	ATCATTCGTTACAGGAAGGTAA AAGTT	27	33.33	59.83
CPV-400_14_LEFT	2	ACACCTGAGAGATTACATATAT AGCACA	29	34.48	60.63
CPV-400_14_RIGHT	2	ACCTTTCCACCAAAAATCTGAG TAAG	26	38.46	60.18

Name	Pool	Sequence (5' – 3')	Size (nt)	%GC	Tm
CPV-400_15_LEFT	1	CAAATGGTCAAATTTGGGATAA AGAATTTG	30	30.0 0	60.1 5
CPV-400_15_RIGHT	1	TTCTAGGTGCTAGTTGATATGTA ATAAACA	30	30.0 0	59.5 6
CPV-400_16_LEFT	2	TGTTTATTACATATCAACTAGCA CCTAGAA	30	30.0 0	59.5 6
CPV-400_16_RIGHT	2	TCTAAGGGCAAACCAACCAAC C	22	50.0 0	61.2 0
CPV-400_17_LEFT	1	AGGTTTGTAGATGGTATACAAT AACTGT	29	31.0 3	59.6 1
CPV-400_17_RIGHT	1	AGCTTTAAATACTAATTTACCTT TCCACCA	30	30.0 0	60.5 0
CPV-400_17.5_RIGHT *	1	AAGTATCAATCTGTCTTTAAGG GG	24	37.5 0	60.1 0
CPV-400_18_LEFT*	2	TATAAGGTGAACTAACCTTACC ATA	25	32.0 0	59.2 0
CPV-400_18_RIGHT*	2	TTAATATAATTTTCTAGGTGCTA GTTG	27	25.9 0	59.2 0

Table 1. Primers targeting 400 bp amplicons

Name	Pool	Sequence (5'-3')	Size (nt)	%GC	Tm
CPV-1000_1_LEFT	1	CTGACCAAGTTCACGTACGTAT GA	24	45.8 3	60. 93
CPV-1000_1_RIGHT	1	TGTTCAAGTGTAAGTGTGCTGA CA	24	41.6 7	61.1 8
CPV-1000_2_LEFT	2	GTGAATGGGTGACTATATTAACA TACAGAC	30	36.6 7	60. 73

Name	Pool	Sequence (5'-3')	Size (nt)	%GC	Tm
CPV-1000_2_RIGHT	2	ACCAAACCAAAGTCTCCTGGAA G	23	47.8 3	60. 95
CPV-1000_3_LEFT	1	AAGCAAATTGAACCAACTCCAG T	23	39.13	59. 55
CPV-1000_3_RIGHT	1	GGTGGTGGTTTACTTCTTTTAGT TGG	26	42.3 1	60. 90
CPV-1000_4_LEFT	2	CTAAGGACGCTAAAGATTGGGG G	23	52.17	61.0 0
CPV-1000_4_RIGHT	2	G TTCCTGTAGCAAATTCATCAC CTG	25	44.0 0	60. 77
CPV-1000_5_LEFT	1	CCATCTCATACTGGAAGTAGTG GC	24	50.0 0	60. 82
CPV-1000_5_RIGHT*	1	TGGATTCCAAGTATGAGAGGCT CT	24	45.8 3	61.2 1
CPV-1000_6_LEFT	2	AACCAAGACTTCATGTAAATGC ACC	25	40.0 0	60. 66
CPV-1000_6_RIGHT*	2	TGGATTCCAAGTATGAGAGGCT CT	24	45.8 3	61.2 1

Table 2. Primers targeting 1000 bp amplicons

11 We used

 Q5 High-Fidelity PCR Kit - 200 rxns **New England Biolabs Catalog #E0555L** to perform PCR.

Reactions should be set up independent **for each pool**, with a final volume of  25 µL, as bellow:

Reactions for 400 bp products:

* Master Mix 2X -  12.5 µL



* Primer pool 1 **or** 2 - 1.6 μL

* Nuclease-free H₂O - 5.9 μL

* DNA - 5.0 μL

Reactions for 1000 bp products:

* Master Mix 2X - 12.5 μL

* Primer pool 1 **or** 2 - 0.4 μL

* Nuclease-free H₂O - 7.1 μL

* DNA - 5.0 μL

- 11.1 The final concentration in a reaction should be 15 nanomolar (nM) per primer. Therefore, the volume of primers used in a reaction will vary according with the number of primers in the pool. More about this subject can be found in [Quick et al. \(2017\)](#).

- 12 PCR's run method should be set as:

8m 15s

* 98 °C for 00:03:00

* 35 cycles of (98 °C for 00:00:15 ; 63 °C for 00:05:00)

* Hold at 4 °C

Purifying PCR Products with AMPure XP beads (Beckman Coulter)

- 13 Mix amplification products of Pool 1 and Pool 2 (final volume 50 μL).

- 14 We purify the fragments with Ampure XP beads **Beckman Coulter Catalog #A63881** , using a modified protocol (original protocol ca be found [here](#)).

- 14.1 This is step can also be fast-tracked by the use of an automated system.

- 15 Add an equal volume of AMPure XP per sample (50 μL).

- 16 After biding, on a magnetic stand, wash beads once with 80 % (v/v) ethanol.

17 Elute DNA with  20 µL of (sequencing appropriate) buffer.

18 Transfer DNA to a new plate/tube.

19 Quantify your DNA.

19.1 We used  Qubit dsDNA HS Assay kit Thermo Fisher Scientific Catalog #Q32854

Equipment

Qubit™ 3 Fluorometer

NAME

Fluorometer for nucleic acid quantitation

TYPE

Invitrogen

BRAND

Q33216

SKU

<https://www.thermofisher.com/order/catalog/product/br/en/Q33216>^{LINK}

Fluorometer for nucleic acid quantitation

SPECIFICATIONS

Library preparation & sequencing

20 To prepare a library, use the Ligation and Barcoding kits (following manufacturer's instructions).

We used **Ligation Sequencing kit SQK-LSK-109** and **Native Barcoding kits EXP-NBD104 and EXP-NBD114** (Oxford Nanopore)

21 Load library on your flow cell and sequence your sample.



- 21.1 We used **R9.4 Oxford MinION flow cell** (FLO-MIN106) and **MinION Mk1B** device for sequencing.

Protocol references

Quick, J., Grubaugh, N. D., Pullan, S. T., Claro, I. M., Smith, A. D., Gangavarapu, K., Oliveira, G., Robles-Sikisaka, R., Rogers, T. F., Beutler, N. A., Burton, D. R., Lewis-Ximenez, L. L., de Jesus, J. G., Giovanetti, M., Hill, S. C., Black, A., Bedford, T., Carroll, M. W., Nunes, M., Alcantara, L. C., Jr, ... Loman, N. J. (2017). Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples. *Nature protocols*, 12(6), 1261–1276. <https://doi.org/10.1038/nprot.2017.066>