



Mar 19, 2025

🌐 SeV cbVG high (HD) stock preparation from recombinant fluorescent reporter virus

📖 [Journal of Virology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk9xqvv5r/v1>

Carolina Lopez¹

¹Washington University



Carolina Lopez

Washington University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk9xqvv5r/v1>

External link: <https://doi.org/10.1128/jvi.00894-25>

Protocol Citation: Carolina Lopez 2025. SeV cbVG high (HD) stock preparation from recombinant fluorescent reporter virus . protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlk9xqvv5r/v1>

**Manuscript citation:**

Pye SE, Achouri E, Yang Y, Musa A, López CB Validation of diverse and previously untraceable Sendai virus copyback viral genomes by direct RNA sequencing. Journal of Virology 99(8). doi: [10.1128/jvi.00894-25](https://doi.org/10.1128/jvi.00894-25)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 20, 2025

Last Modified: March 19, 2025

Protocol Integer ID: 123173

Keywords: generating cbvg high virus stock, cbvg high virus stock, stock preparation from recombinant fluorescent reporter virus, recombinant fluorescent reporter virus, rescued recombinant virus, recombinant virus, cbvg, stock preparation, virus

Abstract

Protocol for generating cbVG high virus stocks from newly rescued recombinant viruses

Materials

1. Cultured LLCMK2 cells (ATCC, #CCL-7)
2. Tissue culture media and Infection media (See media prepare protocol)
3. 0.05% Trypsin (Thermo, #25300054)
4. PBS (Phosphate-buffered saline)
5. Cell counter and Trypan blue
6. T75 flask, 96 well plate
7. Pipettes
8. 2.0ml Cryogenic vials (Corning, #430488)
9. TPCK-Trypsin (for Sendai virus infection, see SeV infection protocol)

Seed LLCMK2 cells to T75 flask

1. Check LLCMK2 cells to ensure they are in good condition, approximately 90% confluent.
2. Wash cells twice with PBS.
3. Add  3 mL 0.05% trypsin-EDTA to the T75 flask, incubate at room temperature for 3-4 minutes.
4. Add  7 mL TCM to the flask to stop trypsinization and collect cells.
5. Combine all cells together.
6. Mix  10 μ L cells and  10 μ L Trypan blue, then add  10 μ L mixture to a cell counter slide and count the cells.
7. Plate 2,000,000 cells in  10 mL of TCM media to T75 flask. One more flask will be needed for cell counting.
8. Incubate cells at 37 °C overnight.

Note: The size of the flask depends on the final volume you need. You can scale down to use a T25 flask or a 6-well plate. However, since high cbVG (HD) virus stocks usually have a lower titer, and the process to create a HD stock is time-consuming, a T75 flask is recommended.

Infect LLCMK2 cells with new rescued virus with a MOI of 10

2. 1. Prepare infection media containing TPCK-trypsin: Add TPCK-trypsin to the infection media at a final concentration of 2 μ g/mL. Prepare  15 mL of the media for each flask.
2. Collect the cells from T75 flask seeded for counting cells and collect cells number as above.
3. Using the number of cells in T75 flask and the MOI of your virus stock to calculate the virus volume needed for infection: $\text{Virus Volume } (\mu\text{L}) = \text{MOI} \times \text{Number of Cells} / \text{Virus Titer (TCID}_{50}/\text{mL})$.
4. Wash cells twice with PBS.
5. Add the virus and prepared infection media containing TPCK-trypsin at final volume of  5 mL to cells.
6. Incubate the flask at 37 °C for 1 hour and gently shake it every 15 minutes.
7. Remove the infection media from flask.
8. Wash the cells twice with PBS.
9. Incubate the cells at 37 °C for 5 days.



Note: In this protocol, we are using rSeV^{eGFP} as an example.

Virus stock titration

- 3 Follow the protocol "SeV fluorescent reporter virus titration by TCID₅₀" to do the titration

Protocol

NAME

SeV fluorescent reporter virus titration by TCID₅₀

CREATED BY

Carolina Lopez

Preview

Harvest the passage 1 HD virus stock

- 4
 1. Check the infected cells, you should be able see eGFP signal in all cells at the Passage 1. During later Passage HD stocks infection, eGFP signal may only be seen in some of the cells not all cells.
 2. Collect the supernatant, and spin it at  2000 rpm, 4°C for 5 minutes to remove cells debris.
 3. Aliquot the virus, the volume of each tube is depended on how you will use it. Generally,  500 µL per tube is recommended.
 4. Quickly freeze the virus using a dry ice-ethanol mixture.
 5. Label them and store them at  -80 °C

Make P2, P3... HD stocks

- 5 Repeat the whole process to make P2, P3, ... HD stocks